


**Vitalograph<sup>®</sup>**

micro

Model 6300

# User Training Manual



**Vitalograph Ltd. UK**

Maids Moreton, Buckingham, MK18 1SW, England

**Phone:** (01280) 827110

**Fax:** (01280) 823302

**e-mail:** sales@vitalograph.co.uk

**www.vitalograph.co.uk**

**Vitalograph Ltd. Export**

Maids Moreton, Buckingham, MK18 1SW, England

**Phone:** +44 1280 827120

**Fax:** +44 1280 823302

**e-mail:** sales@vitalograph.co.uk

**www.vitalograph.eu**

**Vitalograph Ltd. Hong Kong**

Unit 230, Corporation Park, 11 On Lai Street, Shatin, Hong Kong

**Phone:** +852 2117 2678

**Fax:** +852 2117 2679

**Email:** sales@vitalograph.hk

**www.vitalograph.hk**

**Vitalograph GmbH**

Rellinger Straße 64a, D-20257 Hamburg, Germany

**Phone:** (040) 547391-0

**Fax:** (040) 547391-40

**e-mail:** info@vitalograph.de

**www.vitalograph.de**

**Vitalograph Inc.**

13310 West 99th Street, Lenexa, Kansas, 66215, USA.

**Toll Free:** 800 255 6626

**Phone:** (913) 730 3200

**Fax:** (913) 730 3232

**e-mail:** vitcs@vitalograph.com

**www.vitalograph.com**

**Vitalograph (Ireland) Ltd.**

Gort Road Business Park, Ennis, Co. Clare, Ireland

**Phone:** +353 65 6864100

**Fax:** +353 65 6829289

**e-mail:** sales@vitalograph.ie

**www.vitalograph.ie**

© Copyright Vitalograph 2015

Current Edition (Issue 3)

Cat. No. 07687

 **vitalograph** is a registered trademark

## Table of Contents

### **1. Description of the Vitalograph micro**

1. Indications for Use

### **2. Contraindications and Precautions for use**

### **3. Main Components of the Vitalograph micro**

### **4. Features of the Vitalograph micro**

### **5. Getting the Vitalograph micro ready for use**

### **6. Power Management in the Vitalograph micro**

1. Battery Pack
2. Battery Low Detect
3. Power Save Mode

### **7. Operating the Vitalograph micro**

1. Main Menu
2. New Subject Information
3. Performing a Test Session
4. Performing a Post Test Session
5. Configuration Menu
6. Reports and Printing
7. Storing Results
8. Accuracy/Calibration
9. Setting up a New Flowhead

### **8. Cleaning Instructions**

1. Cleaning and Low Level Disinfection of the Vitalograph micro
2. Disassembling and Cleaning of the Fleisch Flowhead
3. Reassembling and Low Level Disinfection of the Fleisch Flowhead

### **9. Fault Finding Guide**

### **10. Customer Service**

### **11. Consumables and Accessories**

### **12. Explanation of Symbols**

### **13. Technical Specifications**

### **14. CE Notice**

### **15. FDA Notice**

### **16. Declaration of Conformity**

### **17. Guarantee**

## 1. DESCRIPTION OF THE VITALOGRAPH MICRO

The Vitalograph micro is a handheld spirometer designed for use by trained professionals in the doctor's office, clinic, hospital department, etc. for measuring and archiving tests on human subjects. Demographic data is entered via a keypad and stored, together with spirometry test data. Current test data can be viewed on the screen and downloaded/printed via a PC. There are a variety of configuration options available.

Information about the device can be obtained from the About box. This information can be used if any queries are made to Vitalograph or a service agent.

To access the About box:

1. Press the **Configuration** button from the Main Menu .



2. Press the **About** button.



### 1.1 Indications for use

The device is a spirometer which measures patient respiratory parameters including FVC, FEV1, FEV6, PEF, MVV and VC. The Vitalograph micro is a handheld spirometer designed for lung function testing for use on adults and pediatrics, 5 years and older, in a variety of environments such as hospital wards, health centers and private homes under the supervision of a healthcare provider.


## 2. CONTRAINDICATIONS AND PRECAUTIONS FOR USE

1. No modification of this equipment is allowed. Any unauthorised changes to the Vitalograph micro device may compromise product safety and/or data and as such Vitalograph cannot be held responsible and the device will no longer be supported.
2. The Vitalograph micro is not designed as a sterile device. Always follow the safety guidelines given by the manufacturer of cleaning and disinfectant chemicals.
3. For spirometry testing a new Bacterial Vital Filter (BVF) should be used for each subject to protect both device and subject. A BVF is for single subject use only. A delay of at least 5 minutes should be allowed between subjects to allow settling of previously aerosolized particles in the measuring device.
4. Spirometry may support or exclude diagnosis, but it cannot make one (ATS/ERS 2005).
5. Take care not to block the mouthpiece with the tongue or teeth during testing. A 'spitting' action or cough will give false readings.
6. All values displayed by the Vitalograph micro are expressed as BTPS values.
7. Time zero is determined using the back-extrapolated method, from the steepest part of the curve.
8. Subject fatigue may occur during spirometry testing depending on the subjects characteristics e.g age, health status. To avoid this the subject may sit during the testing or alternatively take a break between tests. A subject fatigue warning will appear after 8 manoeuvres and the maximum number of allowed manoeuvres in one session is 20.
9. All spirometry standards recommend checking the accuracy of lung function measuring devices at least daily with a 3-L syringe to validate that the instrument is measuring accurately. The Vitalograph micro should never be outside accuracy limits. In normal use, calibration traceability certification is recommended as a part of the routine annual service. Accuracy should also be checked after cleaning or disassembling the spirometer for any reason, after adjusting calibration or if the flowhead or device has been dropped.
10. After fitting a New Flowhead to the Vitalograph micro it is necessary set-up the flowhead with the device. The procedure is outlined below.
11. Service and repairs should be carried out only by the manufacturer or by Service Agents specifically approved by Vitalograph.
12. Per section 14 below Medical Devices may be affected by cellular telephones and other personal or household devices not intended for medical facilities. It is recommended that all equipment used near the Vitalograph product comply with the medical electromagnetic compatibility standard. If interference is suspected or possible,

switching off the offending device is the normal solution, as is required in aircraft and medical facilities. NON MEDICAL EQUIPMENT MUST BE KEPT OUTSIDE THE PATIENT ENVIRONMENT. Portable and mobile RF communications equipment can affect medical electrical equipment.

### 3. MAIN COMPONENTS OF THE VITALOGRAPH MICRO

The main components for the Vitalograph micro are shown below.



A	<i>PC Software</i>
B	<i>Flowhead</i>
C	<i>Flowhead Release Button</i>
D	<i>LCD/Touch Panel Display</i>
E	<i>Mini USB Port</i>
F	<i>On/Off Button</i>
G	<i>LED</i>
H	<i>Battery compartment (4 x 1.5V AAA Batteries)</i>

### 4. FEATURES OF THE VITALOGRAPH MICRO

The Vitalograph micro features include:

- Very high accuracy, linearity and stability
- Printing to the Vitalograph Report Utility, software included with the device
- Storage of tests and demographic information
- Fleisch pneumotachograph


- Removable flowhead
- Touch screen color display
- Clear sounds for audio feedback
- Choice of predicted values
- Real-time test quality prompts

## 5. GETTING THE VITALOGRAPH MICRO READY FOR USE

1. Insert 4 x 1.5V AAA batteries into the battery compartment on rear of device. Alternatively the Vitalograph micro device may be powered from the USB cable supplied with the device. Connect one end of the USB cable into an available USB connection on a PC and the other end into the USB connection on the device.
2. Connect a BVF to the flowhead.
3. Press the On/Off switch on the front face of the instrument to turn the device on and the Vitalograph micro is ready for use.

## 6. POWER MANAGEMENT IN THE VITALOGRAPH MICRO

The Vitalograph micro can be powered from the PC via the USB cable or from the internal batteries. The LED on the front face of the device and the battery power icon show the power status of the device.

When powered from USB power this icon  will be displayed on the status bar at the top of the screen.

### 6.1 Batteries

The Vitalograph micro uses 4 x 1.5V non-rechargeable AAA alkaline batteries. This allows the device to be used without the USB cable connected to the device.

### 6.2 Battery Power Indications

The Vitalograph micro has a number of battery power indications:



When the batteries are full the Green Battery full Icon is displayed on the Main Menu screen device.  
When the USB power is connected a Plug Icon is displayed.  
The LED will stay Green for both full batteries and USB power.



When the batteries start to run low the Battery Low Icon is displayed and both it and the LED will turn orange.  
You will be allowed to continue to use the device. It is advised that you prepare to change the batteries or plug in the USB cable connected to a computer and continue testing.



When the batteries are approaching fully discharged the Battery Discharged icon will be displayed on the full screen on power up and on the main screen status bar and both it and the LED will turn red.  
It is advised that you change the batteries or attach to a PC using a USB cable to continue testing

### 6.3 Power Save Mode


In order to improve battery life, when powered by batteries only and left unused the screen on the Vitalograph micro device will dim after 30 seconds, go blank after 60 seconds and the device will auto power down if left unused for 2 minutes.

When running off of the USB the screen on the Vitalograph micro device will go blank if left unused for 5 minutes. The device will not auto power down if powered by USB.

Pressing the On/Off button will bring the device out of power save mode.

## 7. OPERATING THE VITALOGRAPH MICRO

After turning on the device for the first time you are presented with the Setup screen for the Date and Time.



Enter the current time and date using the up/down buttons.



Press the enter button to save these settings and continue onto the Main Menu screen.

## 7.1 Main Menu

The Main Menu screen includes the following options - New Subject, VC Test, FVC Test and Post Test.



When turning the device on for the first time the test screen icons will appear greyed out and can't be selected until a subject is created. The Post Test icon will stay greyed out until an FVC pre-test is performed.

In addition to displaying the time and the battery icon the status bar at the top on the screen will show various icons to indicate the following:

1. V - indicates a VC test has been performed
2. F - indicates an FVC has been performed
3. P - indicates post mode

The icons will only appear after the test has been completed.



## 7.2 New Subject Information



1. Select the New Subject button on the Main Menu, to bring you into the New Subject screen .

2. The Subject information fields available are as follows -




3. Age, Height and Gender are on the first screen and are enabled by default. Weight and Population Group are on the second screen and are not enabled by default. They have enable options in the configuration

menu  - Subject Options .

4. To enter information for Age, Height and Weight just touch on the box and type in the information using the touch panel keypad. Units will automatically switch between cm/kg and in/lbs. The Gender is selected by pressing the appropriate option male  or female .

The Population Group is selected from a list by pressing the appropriate option on screen. To access additional Population Groups select the arrow on the right of the screen.

5. Press the Enter button to save the subject to the database and return to the Home Menu.

6. If a value is not entered for Age, Height or Gender then an Error Icon  will appear next to the

empty field when the enter button is pressed, this is to indicate that the predicted values will not appear in the results of any testing done. To exit the new subject screen press the enter button again.

**Note:** With the exception of the first use of the device or after clearing the memory, when the user enters the New Subject Screen the details will appear greyed out. If the user hits the Age, Height or Weight box it will jump



into the Data Entry screen and will clear the data to allow the user to enter new data. When the user returns to the New Subject screen all the other fields will be cleared. The user will need to enter all subject details or an error icon will appear.

## 7.3 Performing a Test Session



### 7.3.1 Checks to Make before Performing a Test Session


Before starting a test session, there are a number of checks which should be made:

1. Ensure that the accuracy of the Vitalograph micro unit was checked recently. (Refer to the section on Checking Accuracy)
2. Ensure a subject is created and the required demographic information is entered for the subject. A test session can be performed without filling in any details for the subject but is advised to do so.
3. Fit a disposable BVF mouthpiece to the flowhead. The use of a disposable noseclip is also recommended.

### 7.3.2 Performing a VC Test

Perform the VC test as follows:

1. Select the 'VC Test' option  from the **Main Menu** .
2. Wait for the 'Exhale to Begin' icon to appear . This indicates that the Vitalograph micro unit is ready to accept a blow.

3. The results may be viewed as either a Volume/time (V/t) or Volume Bar graph by pressing the graph button on the side of the test screen. These are not enabled during test.  Volume/time

(V/t)  Volume

4. The graph may be changed to a full screen graph by using the zoom button on the side of the test screen. To return to normal mode select the zoom in button. These are not enabled during test.



5. Perform VC test. Additional guidance on performing Spirometry Manoeuvres can be found on the Vitalograph website.
6. Repeat the blow three times or more to obtain good test quality.
7. The results summary on the bottom of the screen shows the VC of the last blow. The number of blows is shown in a separate box next to the last test VC.
8. If you wish to delete the current blow you can do this as follows:

- a. Following completion of the blow you wish to clear, select the Delete option




from the menu on the side of the test screen.


- b. Two Delete icons will appear in the same area on the side of the test screen, one Delete with a green tick and the other Delete icon with a red 'X'. To confirm the deletion of the blow press the Delete icon with the green tick. To cancel the deletion select the Delete icon with the red 'X'.

- c. If you wish to clear/delete all the sessions on the device you may connect to Vitalograph Reports as outlined in section 6.5 and move all data to Vitalograph Reports and clear it from the device.
9. After performing the VC tests press the enter button to exit the **VC Test** screen. This brings you back to the **Main Menu**.

### 7.3.3 Performing an FVC Test

1. Select the 'FVC Test' option  from the **Main Menu**.


2. Wait for the 'Exhale to Begin' icon to appear . This indicates that the Vitalograph micro unit is ready to accept a blow.

3. The results may be viewed as either a Volume/time (V/t) or Flow/Volume (F/V) graph by pressing the graph button on the side of the test screen. These are not enabled during test.  Volume/time

(V/t)  Flow/Volume (F/V)

4. The graph may be changed to a full screen graph by using the zoom button on the side of the test screen. To return to normal mode select the zoom in button. These are not enabled during test.



5. Perform FVC test. Additional guidance on performing Spirometry Manoeuvres can be found on the Vitalograph website.
6. Repeat the blow three times or more to obtain good test quality.
7. The results summary on the bottom of the screen shows the FVC and FEV1 of the last blow.
8. The number of usable blows and bad blow indicator (!) are shown in a separate box next to the last test FVC and FEV1. There will also be an audible sound at the end of test to indicate a bad blow.
9. The best three tests are shown on the graph in order of rank (best 1, 2, 3...). A legend is shown on the top of the graph to explain the order of the tests.
10. To view results select the results button  on the side of the test screen.

- You can select the test results you want to view by using the left/right arrows.
- You can scroll through the results for each test by using the up/down arrows. The number of parameters displayed will depend on the configured parameters.
- The tests are shown in order of rank (best is ranked number 1 then 2, 3...).
- The results screen has several different columns, arranged in a similar manner to the printout. The first column will display the parameter name, the second displays the units, the third the test value and the fourth column shows the %Pred or Z-value depending on the configuration.

11. If you wish to delete the current blow you can do this as follows:

- Following completion of the blow you wish to clear, select the Delete option



from the menu on the side of the test screen.

- Two Delete icons will appear in the same area on the side of the test screen, one

Delete with a green tick and the other Delete icon with a red 'X'. To confirm the deletion of the blow press the Delete icon with the green tick. To cancel the deletion select the Delete icon with the red 'X'.

- c. If you wish to clear/delete all the sessions on the device (with the exception of the latest FVC Pre-test), you can connect to Vitalograph Reports as outlined in section 6.5 and move all data to Vitalograph Reports and clear it from the device.

12. After performing the FVC tests press the enter button to exit the **FVC Test** screen. This brings you back to the **Main Menu**.

**Note:** Different tests conducted during the same session i.e VC, FVC, Post will be treated as a single session and will be printed as one report. If more than one test is required for the same subject the device should be switched off and on again between tests so that they are registered as separate sessions and separate reports can be generated.

**Note:** A session ends and is saved when one of the following occur - turn off the device, create a new subject, connect to Vitalograph Reports.

## 7.4 Performing a Post Test Session

A Post test session can be performed on the last FVC pre-test session performed. The device will retain the last pre-test even when it is turned off and on again and/or the data has been transmitted to Vitalograph Reports.

To perform a Post test:

1. Select 'Post Mode' from the **Main Menu**




2. A Post test session can be performed on the last FVC pre-test session performed.
3. Perform the Post FVC test as outlined in section 6.2.3 Performing an FVC Test.

**Note:** A Post Test may only be selected if a FVC Pre-test has been completed. When you leave the Post FVC test screen and return to the main menu you will not be able to select either the VC or FVC test as you are still in Post mode. These options will be greyed out.

## 7.5 Configuration Menu

To access the Configuration menu press the  icon on the top right corner of the main screen.


In the Configuration menu there are four options as shown. These are for the 1. Subject options (top left icon), 2. Device settings (top right icon), 3. Accuracy and Calibration (bottom left icon) and the 4. About box (bottom right icon).



1. The Subject option allows you to configure the following

- a. Posture - sets the Posture recorded for the session, sitting or standing.
- b. Weight - turn on to enable the option to enter Weight in the subject screen.

- c. Population Group - turn on to enable the option to enter Population Group in the subject screen.





By default Posture will be set to sitting, Weight and Population Group are set to Off.


2. The Device settings allow you to configure the following -


Device Option  . Select this to configure -

- Select % Predicted or Z-score - the parameter selected will then be displayed in the results screen.
- Audio - you can put the device in silent mode by turning this option Off.
- User Passcode - you can turn this on to lock the device and will be prompted for a passcode.
- Set Temperature (up to 2 decimal places) - you can enter the temperature here. The default setting is 23°C.



- Parameters  . Choose the Parameters option to select what parameters you want to display on the results screen, use the left/right arrows to navigate between the screens. A maximum of 8 parameters can be selected.

- Date/Time.  . Select this option to set or change the Date and/or Time. Date format is YYYY/MM/DD, the Time is fixed at 24Hr format. Use the up/down arrows to edit these fields.

- Service mode/Technician  . This option is for servicing and technicians, an 8 digit passcode is required to enter this mode.

- The following is a list of all available parameters and definitions of the parameters:

•

<b>Parameter</b>	<b>Definition</b>
VC	Vital capacity (L)
FVC	Forced vital capacity (L)
FEV1	Forced expiratory volume after 1 second (L)
FEV1R	FEV1 divided by the largest VC from the VC or FVC manoeuvre.
PEF L/s	Peak expiratory flow (L/sec)
PEF L/min	Peak expiratory flow (L/min)
FEF25-75	Maximal mid expiratory flow: the mean FEF in the time interval between 25% and 75% of the FVC (L/sec)
FEF75-85	Forced late expiratory flow: the mean FEF in the time interval between 75% and 85% of the FVC (L/sec)
EVC	Expiratory vital capacity (L)
IVC	Inspiratory vital capacity (L)
FIVC	Forced inspiratory vital capacity (L)
FIVC/FVC	Ratio FIVC of FVC
FEV.5	Forced expiratory volume after 0.5 seconds (L)
PIF L/s	Peak inspiratory flow (L/sec)
FMFT	Forced mid-expiratory flow time (sec)
FET	Forced expiratory time (sec)
FEV.5/FVC	Ratio FEV 0.5 of FVC
FEV.75	Forced expiratory volume after 0.75 seconds (L)
FEV.75/FVC	Ratio FEV 0.75 of FVC
FEV1/VC	Ratio FEV1 of VC
FEV1/IVC	Ratio FEV1 of IVC
FEV1/FVC	Ratio FEV1 of FVC
FEV1/FIVC	Ratio FEV1 of FIVC
FEV1/FEV6	Ratio FEV1 of FEV6

FEV1/PEF	Ratio FEV1 of PEF
FEV3	Forced expiratory volume after 3 seconds (L)
FEV3/FVC	Ratio FEV3 of FVC
FEV6	Forced expiratory volume after 6 seconds (L)
FEF25	Forced expiratory flow at 25% of the FVC (L/sec)
FEF50	Forced expiratory flow at 50% of the FVC (L/sec)
FEF75	Forced expiratory flow at 75% of the FVC (L/sec)
FEF0.2-1.2	Mean forced expiratory flow in the volume interval between 0.2 and 1.2 L of the test (L/sec)
FEF 25-75/FVC	Ratio FEF25-75 of FVC
FIV1	Forced inspiratory volume after 1 second (L)
FIV1/FVC	Ratio FIV1 of FVC
PIF L/min	Peak inspiratory flow (L/min)
FIF25	Forced inspiratory flow at 25% of the FVC (L/sec)
FIF50	Forced inspiratory flow at 50% of the FVC (L/sec)
FIF75	Forced inspiratory flow at 75% of the FVC (L/sec)
FIF50-FEF50	Ratio FIF 50% of FEF 50%
FIF50-FEF50	Ratio FEF 50% of FIF 50%
MVVind	Maximum voluntary ventilation indirectly calculated from the FEV1 (L/min)
Rind	Airways Resistance Indirect measurement.
Vext	Extrapolated volume (L)
Vext/FVC	Ratio Vext to FVC
FEV1/EVC	Ratio FEV1 to EVC

3. Accuracy/Calibration. Spirometry standards recommend checking the accuracy of lung function measuring devices at least daily with a 3-L syringe to validate that the instrument is measuring accurately. See section

7.8 for full details on performing an Accuracy Check/Calibration.

4. About Box. Information about the software can be obtained from the About box. This information can be used if any queries are made to Vitalograph or a service agent. This information includes the Model number (6300), Serial number of the device, the Software reference number, the date of the last Accuracy check and the Service Due date.



## 7.6 Reports and Printing

The Vitalograph micro can print pdf reports by connecting it to a PC running the Vitalograph Reports application. The micro has to be in the home screen to connect.

1. To produce PDF reports from the micro, ensure that the software is installed on the PC you wish to report to by running the Vitalograph Reports CD supplied with the micro and following the on-screen instructions.
2. The micro can be connected to the PC either using the USB cable supplied with the device, or for the Bluetooth version a Bluetooth dongle may be used.
3. Ensure that Vitalograph reports is open and the micro is switched on and in the home screen.
4. Additional guidance on using Vitalograph Reports can be found in the Vitalograph Reports help menu or in the 07339 User Training Manual on the supplied CD.

**Note:** When the micro is connected to Vitalograph Reports it will move, not copy, the stored sessions with the exception of the latest FVC Pre-session.

**Note:** Different tests conducted during the same session i.e VC, FVC, Post will be treated as a single session and will be printed as one report. If more than one test is required for the same subject the device should be switched off and on again between tests so that they are registered as separate sessions and separate reports can be generated.

## 7.7 Storing Results

The Vitalograph micro has the capacity to store 750 subject entries with the corresponding session data. It is possible to perform up to 20 blows per session, however only a maximum of the best 3 blows will be stored with each session. The session information will also include the subject details and the best pre-test if it is a Post test session.



The Vitalograph micro is intended to be used as a temporary storage device. The Vitalograph micro can be connected to Vitalograph Reports to produce pdf reports of the session data, this will move all subject/sessions to Vitalograph Reports and flag them as deleted on the device with the exception of the last FVC Pre-test done. This is outlined in more detail in section 6.5.


If more than 750 subject/session entries are stored on the device the existing subject/sessions entries will be deleted on a First In First Out (FIFO) basis i.e the first session entered will be the first to be deleted. Any subjects/sessions previously flagged as deleted by Vitalograph Reports will be counted in the total storage available.

## 7.8 Accuracy/Calibration

The Vitalograph micro should never be outside accuracy limits unless damaged or in a fault condition. In this event, see the fault-finding guide. In normal use, calibration traceability certification is recommended as a part of the routine annual service. It is recommended that you perform an accuracy check on the device and check the temperature daily. ISO26782 recommendations require that the difference between the volume measured by the spirometer and the volume pumped into the spirometer from a syringe is within 3%.

Follow these steps to check the accuracy of the unit.

1. Select the Configuration button on the top right corner of the main menu screen. 
2. Select the Accuracy/Calibration icon. 
3. Enter the Syringe volume and reference using the touch panel keypad.
4. Enter the ambient temperature using the touch panel keypad.
5. Pump air through the flowhead to bring it to ambient temperature. If the flowhead has very recently been used for testing or has come from a cold environment, its temperature should be equilibrated with ambient by pumping air through it from the syringe several times.
6. Press the 'Enter' key to bring you into the Accuracy Check screen and follow the on-screen instructions.



7. The Accuracy Check result is shown in % in the top right corner of the screen. If it is reproducible and within 3% a green tick pass icon will be shown and pressing the Enter key will return you to the main Configuration menu and the Accuracy pass is recorded.
8. In the unlikely event the Accuracy Check result is outside 3% the error icon will be shown and pressing the enter key proceeds to the Calibration Update routine to update the Calibration.
9. The Calibration Update screen will show the volume (L) on the top left corner of the screen, next to the number of strokes.
10. The procedure is the same as for the Accuracy Check. If two of the strokes are within 3% of the reference volume pressing the Enter key will return you to the Configuration menu and the Calibration factor is not updated and a pass is recorded. If outside 3% the error icon is shown and pressing the Enter key will return you to the Configuration menu, the Calibration factor is updated and the Calibration update is recorded.

**Note:** To exit the Accuracy Check screen without performing a check, press the Enter key again to return to the Configuration Menu screen. The accuracy check will not be logged to the Vitalograph micro memory in this case.

When to check accuracy -

1. In accordance with your own established procedures
2. After annual maintenance checks
3. After cleaning or disassembling spirometer for any reason
4. After adjusting calibration
5. If the flowhead or device has been dropped
6. If a new flowhead has been fitted

## 7.9 Setting up a New Flowhead

After fitting a New Flowhead to the Vitalograph micro it is necessary set-up the flowhead with the device.

Follow these steps when setting up a new flowhead.

1. Perform an accuracy check / calibration as per section 7.8 above.



2. Select the Accuracy/Calibration icon. 
3. Enter the Syringe volume and reference using the touch panel keypad.
4. Enter the ambient temperature using the touch panel keypad.
5. Pump air through the flowhead to bring it to ambient temperature. If the flowhead has very recently been used for testing or has come from a cold environment, its temperature should be equilibrated with ambient by pumping air through it from the syringe several times.
6. Press the 'Enter' key to bring you into the Accuracy Check screen.
7. Withdraw the syringe fully.
8. Select the New Flowhead icon 
9. Slowly push in the syringe fully and then withdraw fully, keeping the flowrate below the 0.75 L/sec limit lines on the graph. If this is performed correctly new limit lines of 2.50 L/sec will appear.
10. Using a medium speed push in the syringe fully and then withdraw fully, keeping the flowrate between the 0.75 L/sec and 2.50 L/sec limit lines on the graph. If this is performed correctly new limit lines of 10.00 L/sec will appear.
11. Using a fast stroke push in the syringe fully and then withdraw fully, keeping the flowrate between the 2.50 L/sec and 10.00 L/sec limit lines on the graph.
12. The result is shown in % on the bottom of the screen. If it is within 3% a green tick icon will be shown and pressing the Enter key will return you to the main Configuration menu. If this happens then no adjustments are required for the new flowhead.
13. If the results is outside 3% the error icon will be shown and pressing the enter key proceeds to the new flowhead set-up screen.
14. The procedure and screen is the same. Using a slow, medium and fast stroke push in and withdraw the syringe fully. At the end of the procedure a green tick icon will appear and pressing the enter key will return you to the Configuration menu, the new flowhead has been successfully set-up.

**Note:** : To exit the New Flowhead screen, press the Enter key again to return to the Configuration Menu screen. The result will not be logged to the Vitalograph micro memory in this case.

## 8. CLEANING INSTRUCTIONS

### 8.1. Cleaning and Low Level Disinfecting the Vitalograph micro

The parts of the Vitalograph micro that make up the flowhead, which comes into contact with subjects being tested, require **low level disinfection**. The body of the device may be **cleaned** with an alcohol wipe.

A spirometer is not designed as a 'sterile' device.

A new BVF should be used with each subject to prevent cross contamination, the BVF protects both the subject and the device.

One cleaning cycle should be performed weekly.

It is recommended that the flowhead be replaced annually.

#### Table of Cleaning/Disinfection Methods.

Part	Clean/Low Level Disinfection	Recommended Cleaning/Low Level Disinfection
Case Exterior	Clean	Wiping with a 70% isopropyl alcohol impregnated cloth.  For the screen, lightly wipe the surface with cotton pad or other soft material. <b>NOTE: DO NOT use any chemicals. DO NOT wipe in a circular motion. Strokes should be either up/down or over/back.</b>
Remote Flowhead Attachments (not supplied as standard)	Clean	
White Flowhead Tube (not supplied as standard)	Clean	
Screen	Clean	
Flowhead body	Clean & Low Level Disinfection	See section 8.2 for disassembling and cleaning of flowhead. See section 8.3 for reassembling and low level disinfection of the flowhead
Flowhead cone	Clean & Low Level Disinfection	

## 8.2. Disassembling and Cleaning of the Fleisch Flowhead

1. Hold the device body firmly in your left hand.
2. Hold the flowhead with your right hand, at the same time press and hold the button firmly on the front of the fleisch flowhead.
3. Slide the flowhead away from the device from left to right.
4. Remove the flowhead cone from the flowhead, by twisting and pulling it away from the flowhead.
5. To clean the flowhead body swill vigorously in warm soapy water. Do not attempt to “rub” or “scrub” at capillaries.
6. To clean the flowhead cone wash in warm soapy water. Rub surface to remove any visible soiling.
7. Examine all parts to ensure they are visibly clean. If not visibly clean repeat the cleaning process.
8. Rinse all parts in distilled water.
9. Leave to dry completely before reassembling. Drying the fleisch element may require placing it in a warm place overnight. A drying cabinet is ideal.

*Always follow the safety guidelines given by the manufacturer of cleaning and disinfectant chemicals.*

## 8.3. Reassembling and Low Level Disinfection of the Fleisch Flowhead

1. Examine the fleisch element to ensure that no liquid or particles remain in the holes, grooves or pressure tappings.
2. Examine the rubber grommets at the top of the device to ensure no liquids or particles remain in the holes. Also ensure the grommets are not damaged.
3. Fit the flowhead cone to the flowhead.
4. Slide the flowhead into the grooves in the top cover. The Vitalograph logo and button on the flowhead should be on the same face as the LCD when assembled.
5. Wipe all external surfaces of the flowhead with a 70% isopropyl alcohol impregnated cloth.

It is recommended that an accuracy check is carried out following reassembly to verify

correct operation and accuracy.

## **9. FAULT FINDING GUIDE**

Problem Fault Symptoms:	<ul style="list-style-type: none"> <li>• Accuracy check variations &gt; +/-3%</li> <li>• False readings suspected</li> </ul>
Possible Causes: (In probable order)	<ul style="list-style-type: none"> <li>• Recheck Accuracy/Calibration</li> <li>• Was the correct syringe volume selected?</li> <li>• An accuracy check is required after cleaning/disinfecting the flowhead assembly.</li> <li>• Flowhead body pressure port holes/grommets blocked.</li> <li>• Flowhead fleisch element not dried thoroughly.</li> <li>• Flowhead fleisch element assembly blocked.</li> <li>• Main PCB failure – contact support.</li> </ul>
Problem Fault Symptoms:	<ul style="list-style-type: none"> <li>• Test begins automatically</li> <li>• Volume accumulates automatically without the subject blowing.</li> <li>• Very small VC or FVC test displayed</li> </ul>
Possible Causes: (In probable order)	<ul style="list-style-type: none"> <li>• Flowhead and/or tubing when using remote flowhead not stationary at the start of test. Hold them steady until the 'Blow Icon' appears.</li> <li>• Return to Main Menu and re-enter the test routine.</li> </ul>
Problem Fault Symptoms:	<ul style="list-style-type: none"> <li>• No volume measurements.</li> </ul>
Possible Causes: (In probable order)	<ul style="list-style-type: none"> <li>• Ensure that the grommets on flowhead are not pinched or trapped.</li> </ul>
Problem Fault Symptoms:	<ul style="list-style-type: none"> <li>• Cannot print to PC (Vitalograph Reports Application).</li> <li>• Corrupt or missing data on printout.</li> </ul>
Possible Causes: (In probable order)	<ul style="list-style-type: none"> <li>• Check USB cable is connected between Vitalograph micro and the PC.</li> <li>• Check to ensure the Vitalograph Reports Application is correctly installed.</li> <li>• Check to ensure the required software drivers are installed on the PC.</li> <li>• Main PCB failure – contact support.</li> </ul>
Problem Fault Symptoms:	<ul style="list-style-type: none"> <li>• Cannot read screen.</li> </ul>
Possible Causes: (In probable order)	<ul style="list-style-type: none"> <li>• The batteries may be low. Plug in the USB cable and switch on the device.</li> <li>• LCD failure – contact support.</li> <li>• Main PCB failure – contact support.</li> </ul>

## 10. CUSTOMER SERVICE








Service and repairs should be carried out only by the manufacturer or by Service Agents specifically approved by Vitalograph.

For the names and addresses of approved Vitalograph Service Agents or to arrange spirometry workshops, please refer to the contact information at the start of this manual.



























## 11. CONSUMABLES AND ACCESSORIES




























Cat. no	Description
20303	Nose Clips (200)
28350	BVF - Bacterial/Viral Filters (50)
36020	3-L Precision Syringe
79158	Flow Cone (10)
41422	1M USB Cable
79191	Flowhead Complete
79192	Flowhead Connection Tube
79163	Remote Flowhead Adapter Kit
65030SPR	Vitalograph Reports Application
79166	Stylus (2)




## 12. EXPLANATION OF SYMBOLS/ICONS

	Type BF equipment
	Class II
VA	Power rating (in Watts)
	Voltage DC
	Attention (reference relevant section in manual)
	Manufacturer
	Year of Manufacture
	Attention (reference relevant section in manual)

## micro Icons

	Subject		VC Test
	FVC Test		Post Test
	Settings		Enter
	Subject Options		Device Settings
	Accuracy / Calibration		About
	Age		Height
	Gender- Male		Gender- Female
	Posture - Sitting		Posture - Standing
	Weight On		Weight Off
	Population Group - On		Population Group - Off
	Results Options		Parameters
	Time/Date		Service Mode
	VC Volume-Time Graph		VC Volume Graph
	Zoom Out		Zoom In

			
	Test Results		Blow Now
	FVC Volume-Time Graph		FVC Flow - Volume Graph
	Delete		Error
	Serial Number		Software Number
	Syringe Volume		micro Device
	USB Power		Battery Full
	Battery Low		Battery Empty
	Z Score		% Predicted
	Sound On		Sound Off
	User Passcode - On/Locked		User Passcode - Off/Unlocked
	Temperature		Bluetooth
	Accuracy / Calibration Fail		Accuracy / Calibration Pass
	Strokes not repeatable		Syringe in stroke, Syringe out stroke

			
	<p>Error at last Accuracy Check (shown on start up)</p>		<p>New Flowhead Setup</p>

**Other Labels**



Bluetooth



Power input connector



USB connector



The device must be taken to separate collection at the product end-of-life. Do not dispose of these products as unsorted municipal waste.


**13. TECHNICAL SPECIFICATIONS**



Product	Vitalograph micro
Model	6300
Flow detection principle	Fleisch type pneumotachograph
Back pressure	Less than 0.1kPa/L/second @ 14L/s, complies with ISO26782
Volume detection	Flow integration sampling @ 100Hz
Maximum test duration	90 seconds
Maximum displayed volume	10 L
Volume accuracy	Better than $\pm 3\%$
Linearity	Better than $\pm 3\%$
Voltage/Frequency	5V USB Power/6V AAA Alkaline Batteries
Power	5 Watts
Accuracy when operated in operating temperature range conditions	Flow $\pm 10\%$ Max. flow rate $\pm 16$ L/s Min. flow rate $\pm 0.02$ L/s
Operating temperature range	ISO26782 limits: 17–37°C Design limits: 10–40°C
Performance standards the Vitalograph micro meets or exceeds	EN ISO 23747:2009, ISO 26782 2009 & ATS/ERS 2005
Safety standards	EN 60601 (IEC 60601)
QA/GMP standards	ISO 13485:2003, FDA 21 CFR 820, CMDR SOR 98/282 and JPAL
Size	83 mm x 91 mm x 32 mm net
Weight	0.25 kg net
Storage Temperature	0–50°C
Storage Relative Humidity	10%–95%


Printer	No direct print. Use Vitalograph reports to print to pdf.
Communications	USB as standard. Variant with additional Bluetooth comms available
<p><i>Notes:</i></p> <ul style="list-style-type: none"> <li>• All values displayed by the Vitalograph micro are expressed as BTPS values.</li> <li>• Take care not to block the mouthpiece with the tongue or teeth. A 'spitting' action or coughing will give false readings.</li> <li>• Time zero is determined using the back-extrapolated method, from the steepest part of the curve.</li> </ul>	

#### 14. CE NOTICE

Marking by the symbol  indicates compliance of the Model 6300 micro to the Medical Devices Directive of the European Community. Such marking is indicative that the Model 6300 micro meets or exceeds the following technical standards:

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The Model 6300 micro is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 6300 micro should assure that it is used in such an environment		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The Model 6300 micro uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	N/A as Battery or USB Powered	The Model 6300 micro is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A as Battery or USB Powered	
Voltage Fluctuations/Flicker emissions IEC61000-3-3	N/A as Battery or USB Powered	

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The Model 6300 micro is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 6300 micro should assure that it is used in such an environment			
<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1 kV for input/output lines	N/A as not connected to mains, Battery or USB Powered	
Surge IEC 61000-4-5	±1kV differential mode ±2 kV common mode	N/A as not connected to mains, Battery or USB Powered .	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % 100V (>95% dip in 100V) for 0.5 cycle  40 % 100V (60% dip in 100V) for 5 cycles  70 % 100V (30 % dip in 100V) for 25 cycles  <5 % 100V (>95 % dip in 100V) for 5 sec	N/A as not connected to mains, Battery or USB Powered.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The Model 6300 micro is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 6300 micro should assure that it is used in such an environment			
<b>Immunity test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
<p>Conducted RF. IEC 61000-4-6</p> <p>Radiated RF. IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz in ISM bands</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3Vrms from 150kHz top 80kHz</p> <p>3V/m from 80MHz top 2.5GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p><math>d = 1.2 \sqrt{P}</math> 80MHz to 800 MHz  <math>d = 2.3 \sqrt{P}</math> 800 MHz to 2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

<b>Recommended separation distances between portable and mobile RF communication equipment and the Model 6300 micro</b>			
The Model 6300 micro is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model 6300 micro can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 6300 micro as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.1m	0.1m	0.2m
0.1	0.4m	0.4m	0.7m
1	1.2m	1.2m	2.3m
10	3.8m	3.8m	7.3m
100	12m	12m	23m

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

**Note 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Medical Devices may be affected by cellular telephones and other personal or household devices not intended for medical facilities. It is recommended that all equipment used near the Vitalograph product comply with the medical electromagnetic compatibility standard and to check before use that no interference is evident or possible. If interference is suspected or possible, switching off the offending device is the normal solution, as is required in aircraft and medical facilities.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided,

Portable and mobile RF communications equipment can affect medical electrical equipment.

## 15. FDA NOTICE

Caution: Federal Law restricts this device to sale by, or on the order of a physician.

## 16. DECLARATION OF CONFORMITY

Product: 6300 micro™

Vitalograph hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

- European Medical Devices Directive {MDD} 93/42/EEC, as amended.

This device is classified as 2a per Annex IX of the MDD also meets the provisions of the Essential Requirements, Annex I, via compliance with Annex II of the Medical Devices Directive as per Article 11, section 3a, excluding point 4 of Annex II.



- Canadian Medical Device Regulation {CMDR SOR/98-282}
- FDA Quality System Regulation {QSR} 21 CFR 820.
- EN ISO 13485: 2003. Medical devices. Quality management systems. Requirements for regulatory purposes.

Certifying Body: British Standards Institute {BSI}.

{For 93/42/EEC and CMDR}.

BSI Notified Body #: 0086

Certificate Nos. CE 00772, CE 85553, MD 82182, FM 83550

Signed on behalf of Vitalograph (Ireland) Ltd.

A handwritten signature in black ink, appearing to read 'B. R. Garbe'.

B. R. Garbe.

Group Managing Director

## 17. GUARANTEE

Subject to the conditions listed below, Vitalograph Ltd. and its associated companies, (hereinafter called the Company) guarantee to repair or at its option replace any component thereof, which, in the opinion of the Company is faulty or below standard as a result of inferior workmanship or materials.

The conditions of this Guarantee are:

This Guarantee shall only apply to hardware defects which are notified to the Company or to its accredited distributor within 1 year of the date of purchase of the equipment, unless otherwise agreed in writing by the Company.

Software (meaning computer software, or user installable modules) is guaranteed for 90 days from the date of purchase.

The Company warrants that the software when correctly used in conjunction with the hardware will perform in the manner described in the Company's literature and user manuals. The Company undertakes to rectify at no expense to the customer any software failure notified within the period stated above, provided that the failure can be recreated and the software has been installed and used in accordance with the user manual. Notwithstanding this clause, the software is not warranted to be free of errors.

This Guarantee does not cover any faults caused by accident, misuse, neglect, tampering with the equipment, use of consumable items or parts not approved by the Company, or any attempt at adjustment or repair other than by personnel accredited by the Company, nor does it cover reinstatement of any configuration changes caused by the installation of any software.

If a defect occurs please contact the supplier from it was purchased for advice. The Company does not authorize any person to create for it any other obligation or liability in connection with Vitalograph® equipment.

This Guarantee is not transferable and no person, firm or company has any authority to vary the terms or conditions of this guarantee.

To the maximum extent permitted by law, the Company does not accept liability for any consequential damages arising out of the use of, or inability to use any Vitalograph® equipment.

This Guarantee is offered as an additional benefit to the Consumer's statutory rights and does not affect these rights in any way.