

VacSacc®

BactiClear®

ANTIMICROBIAL Disposable Suction Liner Range



Our **commitment** is to **quality**

PART OF **DCC VITAL**
PHARMA DEVICES LOGISTICS

BactiClear® ANTIMICROBIAL

Disposable Suction Liner Range

The only system where all external surfaces are totally antimicrobial, designed to minimise the risk of cross infection.

Liners are available in 1 Litre, 2 Litre and 3 Litre versions. With canisters which fit any bracket within the hospital, suitable for ALL Clinical Areas.

Use as a completely closed system with the standard or PreGelled suction liners. Sample port versions for the addition of solidifying gel or for using with a sample sock are also available. All our liners with the exception of the cascade versions have a hydrophobic bacterial shut off filter for complete overflow protection.

Our unique patented designed liner is folded for ease of storage and begins to inflate immediately the vacuum is attached. With the integral carry handles, the liner allows for the safe removal from the canister, standing upright, whilst full of fluid, this makes

BactiClear® Antimicrobial Disposable Suction Range from **VacSax®** the product of choice



Cascade Stand detail



8V 2 Litre Cascade Stand set up



Cascade Stand detail


BactiClear® Range
(see reverse for specifications)



1, 2 & 3 Litre Liners and Canisters.



Manufactured in the UK to exacting standards.



VacSax BactiClear® Antimicrobial Suction Liner with Hydrophobic Bacterial Filters

Part No.	Description	Quantity
9910-340	1 Litre Antimicrobial Liner	box of 25
VAL-102	1 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-201	2 Litre Antimicrobial Liner	box of 30
VAL-301	3 Litre Antimicrobial Liner	box of 20
VAL-202	2 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-104	1 Litre Antimicrobial Liner with sample port	box of 30
VAL-204	2 Litre Antimicrobial Liner with sample port	box of 30
VAL-304	3 Litre Antimicrobial Liner with sample port	box of 20



VacSax BactiClear® Antimicrobial Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
VAL-203	2 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 30
VAL-303	3 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 20
VAL-205	2 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 30
VAL-306	3 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 20



BactiClear®

ANTIMICROBIAL

VacSax BactiClear® Antimicrobial PreGelled Liners with Hydrophobic Bacterial Filters

Part No.	Description	Quantity
9810-101	1 Litre PreGelled Antimicrobial Liner	box of 25
GAL-211	2 Litre PreGelled Antimicrobial Liner	box of 30
GAL-311	3 Litre PreGelled Antimicrobial Liner	box of 20
GAL-212	2 Litre Wide Bore PreGelled Antimicrobial Liner	box of 30



VacSax BactiClear® Antimicrobial PreGelled Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
GAL-213	2 Litre PreGelled Antimicrobial Cascade Liner including Cascade Tubing	box of 30
GAL-313	3 Litre PreGelled Antimicrobial Cascade Liner including Cascade Tubing	box of 20



BactiClear Liners - Patent No. GB2497360

ECOMED

Serving The HealthCare Industry

Australia T: 1800 67 5432 | F: 1800 67 5433
 E: sales@ecomед.com.au | www.ecomed.com.au
 New Zealand T: 0800 300 421 | F: 0800 300422
 E: sales@ecomед.co.nz | www.ecomed.co.nz