



## Disposable Pulse Oximeter Sensors

- Affordable alternative
- Direct Connection
- Gentle for all patients
- Compatible with all popular oximeters







# Disposable Pulse Oximeter Sensors

Ecomed disposable pulse oximeter sensors are designed for comfort and efficiency, featuring a unique pre-aligned design for accurate placement.

Available to suit most oximeter products, keeping up to date with the latest technology and features.

- Affordable Oxygen Saturation Measurement
- Direct connection to your existing monitor
- Effective infection control, single patient use.
- Compatible with all popular oximeters
- Hypoallergenic adhesion, gentle for all ages
- High quality manufacture with accurate results +/- two digits (80% - 100%)

Compatibility	Monitors	Foam Tape Sensors			
		Adult 45cm Patients > 30kg	Paediatric 45cm Patients 10- 50kg	Infant 90cm Patients 1-20kg	Neonatal 90cm Patients < 3kg
					
BCI	3100, 3101, 3300, 3301, 3302, 3303, 3304, 3401, 6100, 9110, Advisor, Autocorr, Minitorr, Minitorr Plus	PMA-POD-0055	PMA-POD-0056	PMA-POD-0057	PMA-POD-0058
Criticare	503 Spot, 504Plus, 504Dx, 506DXN, 501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	PMA-POD-0059	PMA-POD-0060	PMA-POD-0061	PMA-POD-0062
Datex	AS/3, CS/3, CH/S, CH/RS, CG/2S, Cardiocap II, Ultima, Light Monitor, M-ESTPR, Oscar II, Satelite Trans Plus, ULTSVI, ULTSVO, ULTSIO, ULTS00, OSP-200	PMA-POD-0063	PMA-POD-0064	PMA-POD-0065	PMA-POD-0066
Nellcor	All monitors without OxiMax technology	PMA-POD-0047	PMA-POD-0048	PMA-POD-0049	PMA-POD-0050
Nellcor	All Nellcor Models	PMA-POD-0047OX	PMA-POD-0048OX	PMA-POD-0049OX	PMA-POD-0050OX
Nonin	8500, 8600, 8700, 8800	PMA-POD-0071	PMA-POD-0072	PMA-POD-0073	PMA-POD-0074
Novamatrix	505, 510, 511, 515 B/C, 520, 7000, 7100, Oxypleth	PMA-POD-0067	PMA-POD-0068	PMA-POD-0069	PMA-POD-0070
Ohmeda	3700, 3710, 3770, 3775, Biox3740, 3760, 3800, 3900, 4700, Tuffsat, Oxicap, 5250 RGM, Rascal II	PMA-POD-0051	PMA-POD-0052	PMA-POD-0053	PMA-POD-0054