

| | élance 5 | élance 5i | élance 5c | élance 7 | élance 7i | élance 7c |
|---------------------------------|----------|-----------|-----------|----------|-----------|-----------|
| Patient Type | | | | | | |
| Pediatric | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| Adult | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| Display | | | | | | |
| 10.2" Widescreen LCD | ▼ | ▼ | ▼ | | | |
| 12.1" Widescreen LCD | | | | ▼ | ▼ | ▼ |
| Touchscreen Controls | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| 800 x 480 resolution | ▼ | ▼ | ▼ | | | |
| 1280 x 800 resolution | | | | ▼ | ▼ | ▼ |
| Waveforms | up to 4 | up to 5 | up to 5 | up to 6 | up to 6 | up to 6 |
| Parameters | | | | | | |
| ECG - 3 & 5 Lead | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| Resp | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| NIBP | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| SpO ₂ | | | | | | |
| Spacelabs Std | Opt | Opt | Opt | Opt | Opt | Opt |
| Masimo | Opt | Opt | Opt | Opt | Opt | Opt |
| Nellcor | Opt | Opt | Opt | Opt | Opt | Opt |
| Temp - 2 ch | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| IBP - 2 ch | | ▼ | ▼ | | ▼ | ▼ |
| EtCO ₂ - Microstream | | | ▼ | | | ▼ |
| Features | | | | | | |
| ST Segment Analysis | Opt | Opt | Opt | Opt | Opt | Opt |
| Enhanced Arrhythmia Analysis | Opt | Opt | Opt | Opt | Opt | Opt |
| 120 Hrs. Trends | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| Ext. 2 Ch Recorder | Opt | Opt | Opt | Opt | Opt | Opt |
| Networking | Opt | Opt | Opt | Opt | Opt | Opt |
| Power | | | | | | |
| Internal Power Supply | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| Li-ion, Internal Battery | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |
| 5-Hour Battery Time | ▼ | ▼ | ▼ | ▼ | ▼ | ▼ |

Distributed by



☎ 1800 67 5432

📄 1800 67 5433

✉ sales@ecomед.com.au

🌐 www.ecomed.com.au