

seca 201 & seca 203

Ergonomic Measuring Tape

product sheet [seca 201 & seca 203

new



- Retractable measuring tape for precise body measurements.
- WHR calculator for quick determination of body fat distribution (only 203).
- WHR scale of metal for long life (only 203).

seca 201 & seca 203:

Two high-quality medical products.

In the seca 201 and 203 you've found exactly the right products if you're in search of a measuring tape that is clearly a cut above the rest. Or if you're looking for a measuring tape with millimeter markings that can accurately measure to five millimeters.



seca 201:

The seca 201 is a true medical product. It's very easy to pull out the tape, which just as easily clicks back into place. Handling is extraordinarily convenient. Made of durable materials, the scale features a tape with a top measurement of 205 centimeters and its own low weight of only 55 grams.

seca 203:

The seca 203 also has an extra instrument to calculate the Waist-to-Hip Ratio (WHR). This ratio is an indirect method of determining the amount of fatty tissue in the abdomen and consequently, the

distribution of fat reserves. The high-quality scale is made of metal so that you'll be able to depend on this valuable additional function for a long time.



seca 201:
High-quality materials for precise functionality.



Millimeter graduation for accurate measurements.



seca 203:
Display to determine body fat distribution.

seca 201 & seca 203

Technical Data

seca 201:

- Measurement range: 205 cm
- Division: 1 mm
- Diameter: 70 mm
- Height: 22 mm
- Weight of device: about 55 g

seca 203:

- Measurement range: 205 cm
- Division: 1 mm
- Diameter: 70 mm
- Height: 28 mm
- Weight of device: about 75 g
- Function: WHR calculator