


seca 211

Disposable measuring tape
for head circumference in
practical wall dispenser

new



- Hygienic disposable measuring tape made of sturdy material to prevent cross infection.
- Always at your fingertips thanks to the exclusive wall dispenser that can be installed simply anywhere.
- Exact results can be read easily.
- Scale from 12 to 59 cm and from 5 to 23 inches.

seca 211:

The measuring tape for daily use.

Our experts are constantly striving to make your routine work as easy as possible. That's why seca now offers the environmentally friendly and high-quality disposable measuring tape seca 211 in a practical wall dispenser for improved hygienic conditions. Now you have the desired measuring tape at hand when you need it. The small, space-saving wall dispenser is easy to install and to refill.



Always at hand in a wall dispenser.

The seca 211 is available when needed. One hundred measuring tapes fit into this small wall dispenser that releases them one at a time. The dispenser can be mounted easily anywhere at all. Refilling is also a simple task.



Simple refilling of up to 100 measuring tapes.

Optimum handling.

The seca 211 has a reading aid with highly legible measurement points for precise results. Thanks to the carefully chosen weight of the paper, the measuring tape is strong enough to do the job and light enough to tear easily from the wall dispenser. It's an environmentally friendly solution that can be recycled simply.



Highly legible markings allow easy and precise reading of results.

Hygienic one-time use.

The disposable measuring tape seca 211 fulfills demanding medical hygiene standards. One-time use of the stable measuring tapes to prevent cross infection.

seca 211

Technical Data

- Measuring range: 12–59 cm / 5–23"
- Graduation: 1 mm / 1/16"
- Weight of device: 1.6 g / 0.06 oz
- Dimensions of wall dispenser (WxHxD): 50 x 73 x 36 mm / 1.97 x 2.87 x 1.42"
- Quantity: 500 / 1.000