

seca 777 series

seca
Precision for health

Time to say good buy!



ECOMED
Serving The HealthCare Industry

seca 797 Premium



The all-inclusive EMR validated scale

-
-
-
-
-
-
-
-
-
-
-
-
-
-
-



Integrated handrail seca 477

- + Enhanced patient safety without affecting the accuracy of each measurement
- + Extra support and stability during the weighing process
- + Enables a firm stance for patients in need of assistance



seca integration solutions

A clear improvement to your workflow

The seca 787 and seca 797 can be integrated into your information system via the existing network. At the push of a button, measured values can be sent directly to the electronic medical record (EMR), where they are immediately available to all users.

Connected with an optional USB barcode scanner, both patient and user are identified and the measurement results are always correctly assigned. This saves time, increases process reliability and eliminates transcription errors.

Get further information on seca integration solutions at: integration.seca.com



The seca integration workflow is supported by several system providers:



Technical data

	seca 777	seca 787	seca 797
Functions	TARE, HOLD, lbs/kg switch-over lock, BMI	TARE, HOLD, Auto-HOLD, Damping (motion filter), lbs/kg switch-over lock, automatic BMI, auto-off	TARE, HOLD, Auto-HOLD, Damping (motion filter), kg/lbs switch-over lock, automatic BMI, USB-port and holder for scanner
Data transmission	USB interface for customer specific solution	LAN or Wi-Fi via interface module seca 452 (optional)	Integrated Wi-Fi module
Power supply	Batteries, power adapter (optional)	Batteries, power adapter (optional)	Power adapter, batteries (optional)
Accessories	Handrail seca 477, Panda seca 459, Power adapter 401	Handrail seca 477, Panda seca 459, Power adapter 401	Handrail seca 477, Panda seca 459
Weight	28.0 lbs / 12.7 kg	28.4 lbs / 12.9 kg	28.9 lbs / 13.1 kg
Measuring range (total)	4 – 90 inch / 10 – 230 cm	4 – 90 inch / 10 – 230 cm	4 – 90 inch / 10 – 230 cm
Measuring range (digital)		48 – 90 inch / 122 – 230 cm	48 – 90 inch / 122 – 230 cm
Capacity		550 lbs / 250 kg	
Increments		0.2 lbs / 100 g	
Dimensions (WxHxD)		16.1 x 53.4 x 22.7 inch / 41.0 x 135.6 x 57.6 cm	
Platform (WxHxD)		13.8 x 2.3 x 16.5 inch / 35.0 x 5.9 x 41.8 cm	
Electromagnetic radiation resistance		60601-1-2 compliant	

We reserve the right to make modifications.