

Digital flat scale with Bluetooth interface and high capacity

seca

Precision for health

ECOMED
Serving The HealthCare Industry

seca 813 bt

- + Solid steel frame, up to 440 lbs capacity
- + Non-slip rubber base for safe standing during measuring
- + Large LCD digits for easy reading
- + Transmits weight via Bluetooth
- + Compatible with Google Fit and Apple Health for simple weight trend tracking

 **Bluetooth®**

 Works with
Apple Health

 **Google Fit**



seca 813 bt

Digital flat scale with Bluetooth interface and high capacity



Pleasant and comfortable weighing

All platforms are not created equal. For this weighing platform seca has chosen a high quality, two-component raised profile material which is soft like rubber and pleasant to the touch. But don't be fooled – this material is also tough. It fulfills two requirements at the same time: comfortable weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

A robust, sturdy scale – with a touch of elegance

seca has designed the seca 813 bt with a Bluetooth interface to transfer the Auto-HOLD value automatically to an external device. Thanks to its solid steel frame, the seca 813 bt can withstand high weight and features large LCD digits for easy reading. Furthermore, its wide platform with a low-profile facilitates trouble-free weighing.

Capacity	440 lbs / 200 kg
Graduation	0.2 lb / 100 g
Dimensions (W x H x D)	17 x 1.9 x 14.7" / 433 x 47 x 373 mm
Net weight	6.4 lbs / 2.9 kg
Power supply	Batteries
Interfaces	Bluetooth
Functions	Tap-on automatic switch-on, Auto-HOLD, automatic switch-off, kg/lbs switch-over, Bluetooth connectivity
Scope of delivery	4 x AAA batteries



Paying attention to details

Weight and girth are no problem for the seca 813 bt. An issue with many scales is the display. Thanks to the 1.1 inch high LCD digits, it is extremely easy to read the measurements. You see, a scale from seca makes life more simple by paying attention to even the slightest details.

We reserve the right to make modifications.

201181304021