

# seca 876

## Flat scale for mobile use



- Mother/child-function for trouble-free weighing of small children.
- Large leveling base for high stability.
- Easy to transport at only 9.3 lbs.
- Battery powered.

## seca 876:

### A lightweight for heavyweight use.

In the seca 876 you've found the right solution if you're looking for a lightweight scale (only 9.3 lbs) that's stable, compact, easy to clean and equipped with a mother/child-function.



#### Great stability.

The seca 876 is very light at only 9.3 lbs, but also very tough. Equipped with a steel frame, the scale also boasts a top-quality ribbed surface, extra-strength two-component ABS plastic and Plexiglas®. Small children and adults up to 550 lbs can be weighed accurately. The large leveling bases keep the scale stable on any type of flooring.



Flat construction and slip-proof surface for a comfortable foothold.



Large display for easy reading.



Tip-on automatic.

#### Light but impact resistant.

The combination of steel, ABS plastics and Plexiglas is not only very stable but also especially light. Just as important, the seca 876 is mobile and its large leveling base allows the scale to be positioned without any trouble regardless of the surface underneath it. It's easy to step onto the scale, thanks to the flat construction and slip-proof standing surface.



Large leveling base on the bottom of the scale for safe placement.

#### Mother/child-function.

The scale is easy to operate even with the display and operating elements integrated in the base. Just a light tap of the toe on the scale surface turns on the scale. Infants and small children can be easily weighed while being held by the mother. All she has to do is weigh herself separately first. Then the mother/child-function ascertains the tare of the weight.

## seca 876

### Technical data

- Capacity: 250 kg / 550 lbs
- Graduation: 100g < 150 kg > 200 g / 0.2 lbs < 330 lbs > 0.5 lbs
- Dimensions (WxHxD): 321 x 60 x 356 mm / 12.6 x 2.4 x 14.2"
- Weight of scale: 4.2 kg / 9.3 lbs
- Power: batteries
- Functions: mother/child-function, tip-on automatic switch-on, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, carrying case seca 414 (for combination with stadiometer and adapter element), carrying case seca 421

new



seca 415

Carrying case for the flat scales seca 874, seca 876 and seca 803.

Page 61

---

330 x 380 x 70 mm  
13 x 15 x 2.8"

---

800 g / 1.8 lbs