

Your requirements for cuffs vary by care area and patient need. Statcorp Medical cuffs may be ordered with connectors that make them compatible with most brands of automated blood pressure monitors and wall or portable manometers.

Durable, leak-proof, comfortable, accurate, and economical

Features

Three styles are available:

UltraCheck Reusable

- PVC / DEHP free
- Anti-microbial protection built into the fabric
- Fluid barrier

SoftCheck

- Soft fabric prevents perspiration under cuff

SoftCheck Vinyl

- Soft vinyl is comfortable yet can be wiped off
- Now available in Yellow Vinyl for Isolation areas

- Meets AHA/ AAMI Standards
- Color coded
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Latex free
- Fabric fold at elbow is more comfortable than a weld
- Statcorp adapters allow for compatibility with monitors worldwide





UltraCheck® & SoftCheck®

Fittings



DM



HP



ML



FL



MQ



MQML



MQMF (MQ + MQML)

Manual bulbs



UltraCheck



UltraCheck with gauge



SoftCheck

Yellow Vinyl



For Isolation Patients

Description	Single tube	Dual tube	Limb circumference
UltraCheck Infant	US0814--	UD0814--	08 – 14 cm
UltraCheck Child	US1320--	UD1320--	13 – 20 cm
UltraCheck Small Adult	US1826--	UD1826--	18 – 26 cm
UltraCheck Adult	US2635--	UD2635--	26 – 35 cm
UltraCheck Adult Long	US2938--	UD2938--	29 – 38 cm
UltraCheck Large Adult	US3242--	UD3242--	32 – 42 cm
UltraCheck Large Adult Long	US3544--	UD3544--	35 – 44 cm
UltraCheck Thigh	US4250--	UD4250--	42 – 50 cm
SoftCheck Infant	ST0814--	DT0814--	08 – 14 cm
SoftCheck Child	ST1320--	DT1320--	13 – 20 cm
SoftCheck Small Adult	ST1826--	DT1826--	18 – 26 cm
SoftCheck Adult	ST2635--	DT2635--	26 – 35 cm
SoftCheck Adult Long	ST2938--	DT2938--	29 – 38 cm
SoftCheck Large Adult	ST3242--	DT3242--	32 – 42 cm
SoftCheck Large Adult Long	ST3544--	DT3544--	35 – 44 cm
SoftCheck Thigh	ST4250--	DT4250--	42 – 50 cm
SoftCheck Vinyl Infant	VST0814--	VDT0814--	08 – 14 cm
SoftCheck Vinyl Child	VST1320--	VDT1320--	13 – 20 cm
SoftCheck Vinyl Small Adult	VST1826--	VDT1826--	18 – 26 cm
SoftCheck Vinyl Adult	VST2635--	VDT2635--	26 – 35 cm
SoftCheck Vinyl Adult Long	VST2938--	VDT2938--	29 – 38 cm
SoftCheck Vinyl Large Adult	VST3242--	VDT3242--	32 – 42 cm
SoftCheck Vinyl Lrg Adult Long	VST3544--	VDT3544--	35 – 44 cm
SoftCheck Vinyl Thigh	VST4250--	VDT4250--	42 – 50 cm
SoftCheck Yellow Vinyl Small	YVST1826--	VDT1826--	18 – 26 cm
SoftCheck Yellow Vinyl Adult	YVST2635--	VDT2635--	26 – 35 cm
SoftCheck Yellow Vinyl Large	YVST3242--	VDT3242--	32 – 42 cm

--To order cuffs for a specific monitor, add the appropriate suffix to the product number.

For Critikon Dinamap®, add DM, for Datascope®, Spacelabs®, and Philips®, add HP, for GE Healthcare®, add MQ, and for Datex-Ohmeda®, add MQMF.

Statcorp, Inc.
an OSI Optoelectronics company
14476 Duval Place West, Suite 303
Jacksonville, FL 32218
800-992-0014, 904-786-5113
www.statcorpmedical.com

Product is  Marked
BSI
ISO 13485
Certificate Number FM 40458

ECOMED
Serving The HealthCare Industry

Australia:

☎ 1800 67 5432
☎ 1800 67 5433
✉ sales@ecomед.com.au
🌐 www.ecomed.com.au

New Zealand:

☎ 0800 300 421
☎ 0800 300 422
✉ sales@ecomед.co.nz
🌐 www.ecomed.co.nz

Statcorp Medical reserves the right to make changes to this sheet and the product at any time without notice.

Velcro® is a registered trademark of VELCRO USA, Inc. All other product names mentioned herein may be trademarks and/or registered trademarks of their respective owners.

BPCSHT Rev 02