

Petite Elite

Portable Suction Pump - High Suction.

- Robust, quiet, high performance pump.
 - Suitable for use anywhere in or out of the hospital.
 - High performance >25Lpm at 500+ mmHg.
 - Lightweight and portable with optional carry bag.
- Single x 1000ml Standard Jar Assembly, Autoclaveable
 - Robust Moulded Case Construction
 - Encapsulated Bacterial/Hydrophobic Filter Disk
 - Luminous Touch Control Panel
 - Maximum Vacuum Override Button
 - Battery Level/Charge Indicator
 - High Suction – High Flow 80 to 500mmHg. Flow Rate >25Lpm
 - 115v – 240v AC
 - Battery duration up to 1hr 20m
 - Battery 12v DC, 4.3 mAh, Rechargeable
 - 4.59Kg



PSP001: Petite Elite portable Suction Pump with reusable autoclavable cansiter.

PSP004: Petite Elite Portable Suction Pump with Disposable VacSax Canister.

PSP050: Petite Elite Carry Bag

PSP500: 1000ml Standard Jar Assembly

PSP600: Filter Disks