

SUPPLIES & ACCESSORIES PRODUCT CATALOG



SPACELABS
HEALTHCARE

Spacelabs Healthcare reserves the right to make changes to the product at any time without notice.
All brands and product names and trademarks of their respective owners.

Product may not be available for sale in all locations. Check with your Spacelabs representative.

Order Code	Description	UOM	Page(s)
1103	Unifusor Plus, 500mL w/trigger gauge inflation system	ea	40
1104	Unifusor Plus, 1000mL w/trigger gauge inflation system	ea	40
1105	Unifusor Plus, 3000mL w/trigger trigger gauge inflation system	ea	40
1199	Unifusor Plus, trigger gauge, assembly only	ea	40
2100	Jackson-Rees pediatric circuit	ea	107
2150	Rendell Baker pediatric mask set, 4 sizes	ea	107
27311	Aria download cable	ea	102
27380	Aria belt replacement clip	ea	101
42100	Malleable face mask set, 7 sizes	ea	107
58100	Child Limb Clamp Electrodes	pk	104
10110090	Adult Patient flow sensor	ea	106
10110102	3M Patient flow sensor tubing assembly	ea	106
10110103	6M Patient flow sensor tubing assembly	ea	106
10110104	1.8M Patient flow sensor tubing assembly kit, single ended	ea	106
10110106	1.8M Patient flow sensor & tubing assembly kit, single ended	ea	106
10110122	2.5M Patient flow sensor tubing assembly	ea	106
10110218	2M Patient flow sensor tubing assembly	ea	106
10110219	1.5M Patient flow sensor tubing assembly	ea	106
10400001	Basic Y piece breathing circuit w/limb and bag 2010	ea	107
54100026	22 mm connector - 30f scavenging 1971	ea	107
54100027	Filter-Guard 1944 standard bacterial/viral breathing filter, 22mm OD x 22mm ID	ea	107
54100031	1.8M Magill-Mapleson, open circuit	ea	107
54100032	Individual deluxe Bain coaxial circuit 2115	ea	107
54100033	1.6M Basic Y piece breathing circuit 200	ea	107
54100035	Pediatric Y piece breathing circuit 2010-022	ea	107
54100036	Infant T piece mapleson f breathing circuit 2120	ea	107
54100037	Large disposable face mask	ea	107
54100038	Medium disposable face mask	ea	107
54100039	small disposable face mask	ea	107
54100040	Infant disposable face mask	ea	107
54100041	0.5 L Anesthetic Reservoir Bag, 15mm female neck	ea	107
54100042	1.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	107
54100043	2.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	107
54100044	3.0 L Anesthetic Reservoir Bag, 22mm female neck	ea	107
54100063	Flow sensor caps	50/pk	106
70300001	O2 probe with tee	ea	106
000-0027-02	ABP English Patient Diary, 24 hour Monitoring	ea	98
006-0196-00	Z-fold thermal paper, 2- channel printing	10/bx	89
006-0210-02	ABP Thermal paper, 57 mm x 32 ft	6/bx	98
010-1649-00	Masimo RD SET SpO2 Cable, 12 ft.	ea	81, 94
010-1980-00	CO2 mainstream sensor, 92517	ea	71
010-2146-00	Masimo RD Set SpO2 Adapter Cable, RD SET to LNOP, 1.5 ft.	ea	81, 94
012-0283-00	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 46 cm/18 in	1 Set	52
012-0283-02	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 102 cm/40 in	1 Set	52
012-0283-11	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, AAMI, 61 cm/24 in	1 Set	52
012-0284-00	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 46 cm/18 in	1 Set	52
012-0284-02	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 102 cm/40 in	1 Set	52
012-0284-11	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, AAMI, 61 cm/24 in	1 Set	52
012-0285-00	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 46 cm/18 in	1 Set	52
012-0285-02	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 102 cm/40 in	1 Set	52
012-0285-11	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, AAMI, 61 cm/24 in	1 Set	52, 55
012-0286-00	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 46 cm/18 in	1 Set	52
012-0286-02	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 102 cm/40 in	1 Set	52
012-0286-11	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, AAMI, 61 cm/24 in	1 Set	52, 55
012-0295-00	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, red	3/pk	53

Order Code	Description	UOM	Page(s)
012-0295-02	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, red	3/pk	53
012-0295-05	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 102 cm/40 in, red	3/pk	65
012-0295-06	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, red	3/pk	53, 55
012-0295-07	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 61 cm/24 in, red	3/pk	65, 66
012-0295-10	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, white	3/pk	53
012-0295-12	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, white	3/pk	53
012-0295-15	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 102 cm/40 in, white	3/pk	65
012-0295-16	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, white	3/pk	53, 55
012-0295-17	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 61 cm/24 in, white	3/pk	65, 66
012-0295-20	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, black	3/pk	53
012-0295-22	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, black	3/pk	53
012-0295-25	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 102 cm/40 in, black	3/pk	65
012-0295-26	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, black	3/pk	53, 55
012-0295-27	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 61 cm/24 in, black	3/pk	65, 66
012-0295-30	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, green	3/pk	53
012-0295-32	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, green	3/pk	53
012-0295-35	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 102 cm/40 in, green	3/pk	65
012-0295-36	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, green	3/pk	53, 55
012-0295-38	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 61 cm/24 in, green	3/pk	65, 66
012-0295-40	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 46 cm/18 in, brown	3/pk	53
012-0295-42	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 102 cm/40 in, brown	3/pk	53
012-0295-46	ECG Lead Wires, TruLink, DIN unshielded, snap, AAMI, 61 cm/24 in, brown	3/pk	53, 55
012-0295-48	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 61 cm/24 in, yellow	3/pk	65, 66
012-0295-52	ECG Lead Wires, TruLink, DIN unshielded, snap, IEC, 102 cm/40 in, yellow	3/pk	65
012-0296-00	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, red	3/pk	54
012-0296-02	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, red	3/pk	54
012-0296-04	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, red	3/pk	54, 56
012-0296-10	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, white	3/pk	54
012-0296-12	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, white	3/pk	54
012-0296-14	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, white	3/pk	54, 56
012-0296-20	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 46 cm/18 in, black	3/pk	54
012-0296-22	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, black	3/pk	54
012-0296-24	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, black	3/pk	54, 56
012-0296-32	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, green	3/pk	54
012-0296-34	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, green	3/pk	54, 56
012-0296-42	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 102 cm/40 in, brown	3/pk	54
012-0296-44	ECG Lead Wires, TruLink, DIN unshielded, pinch, AAMI, 61 cm/24 in, brown	3/pk	54, 56
012-0325-01	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, IEC, 102 cm/40 in	ea	65
012-0325-10	ECG Lead Wire Set, TruLink, DIN unshielded, 3L pinch, IEC, 61 cm/24 in	ea	65
012-0326-01	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, IEC, 102 cm/40 in	ea	65
012-0326-10	ECG Lead Wire Set, TruLink, DIN unshielded, 3L, snap, IEC, 61 cm/24 in	ea	65
012-0327-01	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, IEC, 102 cm/40 in	ea	65
012-0327-10	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, pinch, IEC, 61 cm/24 in	ea	65, 66
012-0328-01	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, IEC, 102 cm/40 in	ea	65
012-0328-10	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, IEC, 61 cm/24 in	ea	65, 66
012-0498-00	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, red	3/pk	49
012-0498-01	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, green	3/pk	49
012-0498-02	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, yellow	3/pk	62
012-0498-03	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, white	3/pk	49
012-0498-04	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, black	3/pk	49
012-0498-05	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, brown	3/pk	49
012-0498-06	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, red	3/pk	49
012-0498-07	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, green	3/pk	49
012-0498-08	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, white	3/pk	49
012-0498-09	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, black	3/pk	49

Order Code	Description	UOM	Page(s)
012-0498-11	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, yellow/brown	3/pk	49
012-0498-12	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, green/brown	3/pk	49
012-0498-13	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, blue/brown	3/pk	49
012-0498-14	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, orange/brown	3/pk	49
012-0498-15	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, violet/brown	3/pk	49
012-0498-21	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, yellow/white	3/pk	62
012-0498-22	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, green/white	3/pk	62
012-0498-24	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, black/white	3/pk	62
012-0498-25	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, violet/white	3/pk	62
012-0498-26	ECG Lead Wires, TruLink, shielded, snap, IEC, 102 cm/40 in, yellow	3/pk	62
012-0498-27	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, brown	3/pk	49
012-0498-34	ECG Lead Wires, TruLink, shielded, snap, AAMI, 61 cm/24 in, black	5/pk	49
012-0498-39	ECG Lead Wires, TruLink, shielded, snap, AAMI, 102 cm/40 in, black	5/pk	49
012-0605-10	ECG Lead Wire Set, TruLink, DIN unshielded, 5L, snap, 64 cm/25.2 in	1 Set	55
012-0698-00	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, red	3/pk	62
012-0698-01	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, green	3/pk	62
012-0698-03	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, white	3/pk	62
012-0698-04	ECG Lead Wires, TruLink, shielded, snap, IEC, 61 cm/24 in, black	3/pk	62
012-0698-06	ECG Lead Wires, TruLink, shielded, snap, IEC, 102 cm/40 in, red	3/pk	62
012-0698-07	ECG Lead Wires, TruLink, shielded, snap, IEC, 102 cm/40 in, green	3/pk	62
012-0698-08	ECG Lead Wires, TruLink, shielded, snap, IEC, 102 cm/40 in, white	3/pk	62
012-0698-09	ECG Lead Wires, TruLink, shielded, snap, IEC, 102 cm/40 in, black	3/pk	62
015-0067-04Q	ABP Cuff, Spacelabs Cloth, self-applied, sm adult, 17-26 cm	ea	96
015-0068-05Q	ABP Cuff, Spacelabs Cloth, self-applied, adult, 24-32 cm	ea	96
015-0070-00	Cuff support harness	ea	98
015-0118-04Q	ABP Cuff, Spacelabs Cloth, self-applied, child, 12-20 cm	ea	96
015-0312-01	Anesthetic agent sampling line + gas sampling tee, 90515/91518	50/bx	75
015-0313-01	Anesthetic agent sampling line + gas sampling tee, 90515/91518	50/bx	75
015-0321-01	Dixtal SPO2 universal sensor ear clip, reusable	10/pk	78
015-0330-00	SPO2 universal foam wrap	25/pk	78
015-0500-02	Telemetry pouch, single use	100/bx	57
015-0660-00	SPO2 finger sensor, adult, reusable	ea	58, 68, 77, 92, 111
015-0661-00	SPO2 Y multi-site sensor, universal, reusable	ea	77, 92
015-0662-00	SPO2 foam sensor, adult, single use	24/bx	58, 68, 78, 92, 111
015-0663-00	SPO2 vinyl sensor, adult, single use	24/bx	58, 68, 78, 92, 111
015-0664-00	SPO2 foam sensor, pediatric, single use	24/bx	58, 68, 78, 92, 111
015-0665-00	SPO2 vinyl sensor, pediatric single use	24/bx	58, 68, 78, 92, 111
015-0666-00	SPO2 finger sensor, adult, single use	24/bx	58, 68, 78, 92, 111
015-0667-00	SPO2 cloth sensor, neonate, single use	24/bx	58, 68, 78, 111
015-0683-00	CO2 sidestream Nomoline adapter w/sampling line, single use	25/bx	70
015-0701-00	CO2 canister w/LoFloSorb, green to violet	6/cs	106
015-0702-00	CO2 canister w/Spherasorb, white to violet	6/cs	106
015-0715-00	Reusable Pouch, 90227	ea	97
015-0716-00	Reusable Shoulder Strap	ea	97
015-1826-00	ABP Cuff, Spacelabs Soft, self-applied, sm adult, 18-26 cm, no fitting	20/bx	97
015-2635-00	ABP Cuff, Spacelabs Soft, self-applied, adult, 26-35 cm, no fitting	20/bx	97
015-3242-00	ABP Cuff, Spacelabs Soft, self-applied, lg adult, 32-42 cm, no fitting	20/bx	97
015-4250-00	ABP Cuff, Spacelabs Soft, self-applied, xl adult, 42-50 cm, no fitting	20/bx	97
016-0077-05Q	ABP Cuff, Spacelabs Cloth, self-applied, lg adult, 24-32 cm	ea	96
016-0080-00	Nylon web adjustable waist belt, D-ring adjustable fastener	ea	98
016-0109-04Q	ABP Cuff, Spacelabs Cloth, self-applied, xl adult, 38-50 cm	ea	96
016-0262-00	Nylon web adjustable shoulder strap w/velcro fasteners	ea	98
016-0263-00	Vinyl carrying pouch wth Velcro closure, 90207	ea	97
016-0340-00	Carrying pouch, cloth, 90217	ea	97
016-0410-00	Cloth carrying pouch with Velcro closure, 90207	ea	97

Order Code	Description	UOM	Page(s)
016-0428-00	Filter cartridge	ea	76
016-1320-01	ABP Cuff, UltraCheck Nylon, self-applied, child, 13 - 20 cm	ea	96
016-1826-01	ABP Cuff, UltraCheck Nylon, self-applied, small adult, 18 - 26 cm	ea	96
016-2635-01	ABP Cuff, UltraCheck Nylon, self-applied, adult, 26 - 35 cm	ea	96
016-3242-01	ABP Cuff, UltraCheck Nylon, self-applied, lg adult, 32 - 42 cm	ea	96
019-0539-00	Oral/axillary temperature probe, 91330	ea	95
019-0539-00	Replacement Temperature Probe, Oral/Axillary	ea	95
0315BOX	Nikotab Eelectrode	500/bx	104
050-1120-00	O2 sensor Kit	ea	76
073-0267-00	Patient Diary, English, (Aria, EVO, LCF) CE0123	ea	102
103-0059-00	TruLink ECG cable adapter, 17-pin to 6-pin adapter	ea	54
103-0234-00	CO2 sidestream Nomoline adapter, mult-use	25/bx	70, 106
1103FGA-24	Unifusor Plus, 500mL w/aneroid gauge and thumbwheel valve	24/bx	39
1103SGA-24	Unifusor Plus, 500mL w/aneroid gauge and stopcock valve	24/bx	39
1103TGA-24	Unifusor Plus, 500mL w/piston gauge and thumbwheel valve	24/bx	39
1103X	Unifusor Plus, 500mL, cuff only	ea	40
1103X-05	Unifusor Plus, 500mL, cuff only	5/bx	40
1103XB-05	Unifusor Plus, 500mL, cuff only	5/bx	39
1104FGA-24	Unifusor Plus, 1000mL w/aneroid gauge and thumbwheel valve	24/bx	39
1104SGA-24	Unifusor Plus, 1000mL w/aneroid gauge and stopcock valve	24/bx	39
1104TGA-24	Unifusor Plus, 1000mL w/piston gauge and thumbwheel valve	24/bx	39
1104X	Unifusor Plus, 1000mL, cuff only	ea	40
1104X-05	Unifusor Plus, 1000mL, cuff only	5/bx	40
1104XB-05	Unifusor Plus, 1000mL, cuff only	5/bx	39
1105FGA-10	Unifusor Plus, 3000mL w/aneroid gauge and thumbwheel valve	10/bx	39
1105SGA-10	Unifusor Plus, 3000mL w/aneroid gauge and stopcock valve	10/bx	39
1105TGA-10	Unifusor Plus, 3000mL w/piston gauge and thumbwheel valve	10/bx	39
1105X	Unifusor Plus, 3000mL, cuff only	ea	40
1105X-05	Unifusor Plus, 3000mL, cuff only	5/bx	40
1105XB-05	Unifusor Plus, 3000mL, cuff only	5/bx	39
1113SGA-24	Unifusor Plus, 500mL w/piston gauge and stopcock valve	24/bx	39
1114SGA-24	Unifusor Plus, 1000mL w/piston gauge and stopcock valve	24/bx	39
1115SGA-10	Unifusor Plus, 3000mL w/piston gauge and stopcock valve	10/bx	39
118-0256-01	Esophageal/rectal probe, reusable	ea	84
13-9083	Razor with blade cover	ea	103
162-0046-01	Anesthetic agent sampling line, 90515/91518	50/bx	75
162-0047-01	Anesthetic agent sampling line, 90515/91518	50/bx	75
166-7085-00	CO2 sampling line, sidestream, 2m/6.5 ft	25/bx	71, 106
166-7121-00	CO2 sampling line, sidestream, 3m/10 ft	50/bx	71, 106
19-7118	Lifecard Compact USB Flashcard Reader	ea	102
19-7130	Laces for pouch set	ea	101
19-7131	Cordlock for pouch set	ea	101
19-7509	6-Lead wire set, snap, 61 cm/24 in	ea	99
19-7510	7-lead wire set, 61 cm/24 in	ea	99
200-1508-001	Pouch & straps for 485	ea	102
23-0026	Download cable	ea	105
23-0051	Patient cable, 2 electrodes, 1 ch	ea	105
23-0059	Patient cable, 3 electrodes, 2 ch	ea	105
23-0245	Programing card & key, CCVS20	ea	105
23-0252	Programing card & key, CCST80	ea	105
23-0255	Lanyard clip & assembly	ea	105
23189-106	4 electrode hook-up kit	ea	99
23189-106-401	Aria hook-up kit with case, pouch and clear tape electrodes	20/cs	99
23189-110-401	Hook-up Kit (no pouch)	20/cs	99
25-0004	10 Electrode patient cable, 4MM, IEC	ea	104

Order Code	Description	UOM	Page(s)
25-0005	10 Electrode patient cable, 4MM, AHA	ea	104
25-0010	Limb Clamps	4/st	104
25-0020	10 Electrode patient cable, 4MM, AHA	ea	104
25-0027	10 Electrode patient cable, short with banana connector	ea	104
25-0028	10 Electrode patient cable, short, AHA	ea	104
25450-187	indv color lead wire, snap, white, 14 in	ea	99
25452-519	7-Lead wire set, snap, 46 cm/18 in, white/black/red/green/orange/blue/brown	ea	99
25452-525	4-Lead wire set, snap, three 14 in leads and one 7 in white lead	ea	99
25452-529	4-Lead wire set, snap, 102 cm/40 in, Aria and evo belt clip config	ea	99
25452-535	4-Lead wire set, snap, three 14 in leads and one 14 in white lead	ea	99
27249-7-401	Aria disposable pouch	50/bx	101
27565-5	ECG Electrode, foam	5/pk	100
27565-7	ECG Electrode, foam	7/pk	100
27565-9	ECG Electrode, clear tape	4/pk	100
27565-9CASE	ECG Electrode, clear tape	150/cs	100
27685-401	Disposable pouch, LC12	50/pk	102
2803FGA-24	Unifusor II, 500mL w/aneroid gauge and thumbwheel valve	24/bx	38
2803SGA-24	Unifusor II, 500mL w/aneroid gauge and stopcock valve	24/bx	38
2803X-24	Unifusor II, 500mL, cuff only	24/bx	38
2804FGA-24	Unifusor II, 1000mL w/aneroid gauge and thumbwheel valve	24/bx	38
2804SGA-24	Unifusor II, 1000mL w/aneroid gauge and stopcock valve	24/bx	38
2804X-24	Unifusor II, 1000mL, cuff only	24/bx	38
2805FGA-10	Unifusor II, 3000mL w/aneroid gauge and thumbwheel valve	10/bx	38
2805SGA-10	Unifusor II, 3000mL w/aneroid gauge and stopcock valve	10/bx	38
2805X-10	Unifusor II, 3000mL, cuff only	10/bx	38
2903SGA-24	Unifusor II, 500mL w/piston gauge and stopcock valve	24/bx	38
2903TGA-24	Unifusor II, 500mL w/piston gauge and thumbwheel valve	24/bx	38
2904SGA-24	Unifusor II, 1000mL w/piston gauge and stopcock valve	24/bx	38
2904TGA-24	Unifusor II, 1000mL w/piston gauge and thumbwheel valve	24/bx	38
2905SGA-10	Unifusor II, 3000mL w/piston gauge and stopcock valve	10/bx	38
2905TGA-10	Unifusor II, 3000mL w/piston gauge and thumbwheel valve	10/bx	38
305240-001	Roll thermal paper, single-channel printing	10/bx	89
306655-003	Cardiac output cable, Baxter/Edwards/Abbott	ea	88
306655-004	Cardiac output cable, Becton Dickinson/Datex-Ohmeda	ea	88
307438-001	Roll thermal paper, single- or dual-channel printing	10/pk	89
307568-001	Z-fold thermal paper, 2- or 4- channel printing	10/bx	114
344-0056-00	Garment clip	ea	98
369014-003	Nafion tubing, single use	ea	76
369062-001	ECG electrode, adult, tape, single use	bx	50, 57, 63, 67, 110
369063-001	ECG electrode, transparent tape, pediatric/neonatal, single use	bx	51, 64
376-0561-01	Pediatric Patient flow sensor	ea	106
378-0050-00	Water trap, single use	ea	76
392003-001	Skin prep pad	100/bx	103
41-0053EU	12 lead hook-up set, IEC	ea	104
41-0053US	12-lead hook-up set, AHA	ea	104
41-0080	CD12USB stress tests hook-up kit	ea	103
46-0063	Plastic memory card cases	ea	102
46-0092	Fabric pouch	ea	101
46-0408	Fabric pouch set, 2 pouches with laces, cordlock & instructions	ea	101
46-0413	LCF opener tool	ea	102
46-0435	Patient diary, 7 day, German	ea	102
46-0459	Patient diary, 7 day, English	ea	102
46-0512	Disposable pouch, lifecard cf	100/pk	102
46-0573	Patient diary lcf, eng,ita, fra, spa, pol	ea	102
50-0042	Evo holder for waist belt	ea	101

Order Code	Description	UOM	Page(s)
50-0048	4 electrode hook-up kit	ea	99
50-0049	4 electrode hook-up kit with disposable chest pouch	ea	99
50-0050	Disposable chest pouch for evo recorder	50/pk	101
50-0190	4 electrode hook-up kit	20/cs	99
50-0191	4 electrode hook-up kit with chest pouch	20/cs	99
50-0220	LCF / evo pen	ea	102
685-0037-00	ECG pre-wired neonatal electrode, DIN unshielded, 3-lead, single use	100/bx	51, 64
690-0003-01	Nellcor OxiMax DS-100A finger sensor, adult, reusable	ea	83, 95, 113
690-0004-01	Nellcor OXI-A/N sensor, adult/neonatal, reusable	ea	83, 95
690-0008-00	PRESS,TRANSD,PAT-MOUNT KIT,SPU,152.5CM/60IN,5/BX	5/bx	87, 114
690-0009-00	PRESS,TRANSD,POL-MOUNT KIT,SPU,152.5CM/60IN,5/BX	5/bx	87, 114
690-0010-00	Transducer kit, pole mount, 182.9 cm/72 in	5/bx	87, 114
690-0011-00	Pressure Transducer, neonatal	10/bx	87
690-0012-00	Transducer mounting background	ea	87, 114
690-0018-00	Diaphragm dome, Datex-Ohmeda TA101OND	12/bx	87, 113
690-0020-00	Pressure kit, for reusable transducer, 122 cm/48 in	5/bx	87, 114
690-0021-00	INVPB interface cable for pressure transducers, Datex-Ohmeda/Becton Dickinson	ea	86, 113
690-0023-00	Injectate probe, in-line, reusable, Datex-Ohmeda SP4045	ea	88
690-0024-00	Injectate probe, in-line, reusable, Datex-Ohmeda SP5030	ea	88
690-0028-01	Skin temperature probe, small, YSI729	ea	84
690-0031-00	CritiKit SP4572 CO in-line room temperature injectate kit, closed-loop	5/bx	88
690-0032-00	CritiKit SP4500 CO in-line ice temperature injectate kit, closed-loop	4/bx	88
690-0033-00	Injectate probe, in-line, reusable, Datex-Ohmeda SP4042	ea	88
690-0034-00	Pressure transducer holder, reusable	ea	87, 114
690-0039-01	Nellcor OXI-P/I sensor, pediatric/infant, reusable	ea	83, 95
690-0043-01	Injectate probe, in-line, Datex-Ohmeda	25/bx	88
690-0215-01	Nellcor DS-Y multi-site sensor, universal, reusable	ea	83, 95
690-0247-00	Masimo LNCS DC-I sensor, adult, reusable	ea	79, 93, 112
690-0248-00	Masimo LNCS DC-IP sensor, pediatric, reusable	ea	79, 93, 112
690-0249-00	Masimo LNCS TC-I tip clip sensor, reusable	ea	79, 93, 112
690-0250-00	Masimo LNCS TF-I forehead sensor, adult, reusable	ea	79, 93, 112
690-0251-00	Masimo LNCS Adtx adhesive sensor, adult, single use	20/bx	80, 93, 112
690-0252-00	Masimo LCNS Pdtx adhesive sensor, pediatric, single use	20/bx	80, 112
690-0256-00	SpO2 Adapter Cable, Masimo, LNCS, 305 cm/10 ft, Elance	ea	112
690-0258-00	Masimo LNCS Neo-L adhesive sensor, neonatal, single use	20/bx	80
690-0259-00	Masimo LNCS Inf-L adhesive sensor, infant, single use	20/bx	80
690-0261-01	Skin/Surface Temp probe, reusable, YSI-400, 305cm/10ft, Elance	ea	113
690-0262-00	ECG electrode, transparent tape, general/stress, adult	30/bx	50, 57, 63, 67, 110
690-0267-00	Masimo RD SET DCI, Adult Reusable Sensor	EA	81, 94
690-0270-00	Masimo RD SET DCI-P, Pediatric Reusable Sensor	EA	81
690-0271-00	SpO2 Adapter Cable, Nellcor, 305 cm/10 ft, Elance	ea	113
690-0272-00	Esophageal/rectal probe, reusable	ea	84
690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO2 sensor	ea	81
690-0283-00	Masimo RD SET Adt, Adult SpO2 adhesive Sensor	20/bx	81, 94
690-0284-00	Masimo RD SET Pdt, Pediatric SpO2 adhesive Sensor	20/bx	81, 94
690-0285-00	Masimo RD SET Inf, Infant SpO2 adhesive Sensor	20/bx	81, 94
690-0286-00	Masimo RD SET Neo, Neonatal/Adult SpO2 adhesive Sensor	20/bx	81, 94
690-0287-00	Masimo RD SET NeoPt, Neonatal SpO2 adhesive Sensor	20/bx	81, 94
690-0289-00	ECG electrode, transparent tape, adult, single use	30/bx	50, 57, 63, 67, 110
690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	30/bx	51, 64
700-0006-00	ECG Lead Wire Set, TruLink, shielded, 3L pinch, AAMI, 61 cm/24 in	1 Set	48, 110
700-0006-01	ECG Lead Wire Set, TruLink, shielded, 3L pinch, IEC, 61 cm/24 in	ea	61, 110
700-0006-02	ECG Lead Wire Set, TruLink, shielded, 3L pinch, AAMI, 102 cm/40 in	1 Set	48, 110
700-0006-03	ECG Lead Wire Set, TruLink, shielded, 3L pinch, IEC, 102 cm/40 in	ea	61, 110
700-0006-08	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, AAMI, 61 cm/24 in	1 Set	48, 110

Order Code	Description	UOM	Page(s)
700-0006-09	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, IEC, 61 cm/24 in	ea	61, 110
700-0006-10	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, AAMI, 102 cm/40 in	1 Set	48, 110
700-0006-11	ECG Lead Wire Set, TruLink, shielded, 5L, pinch, IEC, 102 cm/40 in	ea	61, 110
700-0006-16	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, AAMI, 61 cm/24 in	1 Set	48
700-0006-17	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, IEC, 61 cm/24 in	1 Set	61
700-0006-18	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, AAMI, 102 cm/40 in	1 Set	48
700-0006-19	ECG Lead Wire Set, TruLink, shielded, 10L, pinch, IEC, 102 cm/40 in	1 Set	61
700-0006-30	ECG Combiner Disp Lead Wire Set, 3L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-31	ECG Combiner Disp Lead Wire Set, 4V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-32	ECG Combiner Disp Lead Wire Set, 5L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-33	ECG Combiner Disp Lead Wire Set, 5V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41, 60
700-0006-34	ECG Combiner Disp Lead Wire Set, 6L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0006-35	ECG Combiner Lead Wire Set, 3L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-36	ECG Combiner Lead Wire Set, 4V, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-37	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-38	ECG Combiner Lead Wire Set, 5V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41, 60
700-0006-39	ECG Combiner Lead Wire Set, 6L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0006-41	ECG Combiner Lead Wire Set, 3L, multi-pinch, IEC, 74 cm/29 in	1 Set	60
700-0006-42	ECG Combiner Disp Lead Wire Set, 3L, multi-pinch, IEC, 127 cm/50 in	25/bx	60
700-0006-43	ECG Combiner Lead Wire Set, 5L, multi-pinch, IEC, 74 cm/29 in	1 Set	60
700-0006-44	ECG Combiner Disp Lead Wire Set, 5L, multi-pinch, IEC, 127 cm/50 in, 25/bx	1 Set	60
700-0007-00	ECG Lead Wire Set, TruLink, shielded, 3L, snap, AAMI, 61 cm/24 in	1 Set	48, 109
700-0007-01	ECG Lead Wire Set, TruLink, shielded, 3L, snap, IEC, 61 cm/24 in	1 Set	61, 109
700-0007-02	ECG Lead Wire Set, TruLink, shielded, 3L, snap, AAMI, 102 cm/40 in	1 Set	48, 109
700-0007-03	ECG Lead Wire Set, TruLink, shielded, 3L, snap, IEC, 102 cm/40 in	1 Set	61, 109
700-0007-08	ECG Lead Wire Set, TruLink, shielded, 5L, snap, AAMI, 61 cm/24 in	1 Set	48, 110
700-0007-09	ECG Lead Wire Set, TruLink, shielded, 5L, snap, IEC, 61 cm/24 in	1 Set	61, 110
700-0007-10	ECG Lead Wire Set, TruLink, shielded, 5L, snap, AAMI, 102 cm/40 in	1 Set	48, 110
700-0007-11	ECG Lead Wire Set, TruLink, shielded, 5L, snap, IEC, 102 cm/40 in	1 Set	61, 110
700-0007-16	ECG Lead Wire Set, TruLink, shielded, 10L, snap, AAMI, 61 cm/24 in	1 Set	48
700-0007-17	ECG Lead Wire Set, shielded, lmb/chst, snap, IEC, 61 cm/24 in	1 Set	61
700-0007-18	ECG Lead Wire Set, TruLink, shielded, 10L, snap, AAMI, 102 cm/40 in	1 Set	48
700-0007-19	TruLink ECG leadwire set, shielded, 10-Lead, limb/chest, snap, 102 cm/40 in, IEC	1 Set	61
700-0007-30	ECG Combiner Disp Lead Wire Set, 3L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-31	ECG Combiner Disp Lead Wire Set, 4V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-32	ECG Combiner Disp Lead Wire Set, 5L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-33	ECG Combiner Disp Lead Wire Set, 5V, multi-pinch, AAMI, 127 cm/50 in	25/bx	41, 60
700-0007-34	ECG Combiner Disp Lead Wire Set, 6L, multi-pinch, AAMI, 127 cm/50 in	25/bx	41
700-0007-35	ECG Combiner Lead Wire Set, 3L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-36	ECG Combiner Lead Wire Set, 4V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-37	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-38	ECG Combiner Lead Wire Set, 5V-Leads, multi-pinch, AAMI, 74 cm/29 in	1 Set	41, 60
700-0007-39	ECG Combiner Lead Wire Set, 6L, multi-pinch, AAMI, 74 cm/29 in	1 Set	41
700-0007-41	ECG Combiner Lead Wire Set, 3L, Snap, IEC, 74 cm/29 in	1 Set	60
700-0007-42	ECG Combiner Disp Lead Wire Set, 3L, Snap, IEC, 127 cm/50 in, 25/bx	1 Set	60
700-0007-44	ECG Combiner Lead Wire Set, 5L, Snap, IEC, 74 cm/29 in	1 Set	60
700-0007-45	ECG Combiner Disp Lead Wire Set, 5L, Snap, IEC, 127 cm/50 in, 25/bx	1 Set	60
700-0008-00	ECG cable, shielded, 10-Lead, AAMI	ea	50
700-0008-01	ECG cable, shielded, 10-Lead, 305cm/10ft, IEC	ea	63
700-0008-06	ECG cable, shielded, 5-Lead, AAMI	ea	50, 109
700-0008-07	ECG cable, shielded, 5-lead, 305cm/10ft, IEC	ea	63, 109
700-0008-08	ECG cable, shielded, 3-Lead, AAMI	ea	50, 109
700-0008-09	ECG cable, shielded, 3-lead, 305cm/10ft, IEC	ea	63, 109
700-0008-12	ECG neonatal cable, DIN unshielded, 3-Lead, AAMI	ea	51
700-0008-13	ECG neonatal cable, DIN unshielded, 3-Lead, 305cm/10ft, IEC	ea	64

Order Code	Description	UOM	Page(s)
700-0008-22	ECG lead wire set w/integrated cable, shielded, 5L, snap, AAMI	ea	109
700-0008-23	ECG lead wire set w/integrated cable, shielded, 5L, snap, IEC	ea	109
700-0008-24	ECG lead wire set w/integrated cable, shielded, 3L, snap, AAMI	ea	109
700-0008-25	ECG lead wire set w/integrated cable, shielded, 3L, snap, IEC	ea	109
700-0008-30	ECG Combiner Cable, 3-Lead, Datascope, AAMI	ea	44
700-0008-31	ECG Combiner Cable, 5-Lead, Datascope, AAMI	ea	44
700-0008-33	ECG Combiner Cable, 3-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-34	ECG Combiner Cable, 5-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-35	ECG Combiner Cable, 6-Lead, Datex-Ohmeda, AAMI	ea	44
700-0008-36	ECG Combiner Cable, 3-Lead, GE, AAMI	ea	44
700-0008-37	ECG Combiner Cable, 5-Lead, GE, AAMI	ea	44
700-0008-38	ECG Combiner Cable, 5+5-Lead, GE, AAMI	ea	44
700-0008-39	ECG Combiner Cable, 6-Lead, GE, AAMI	ea	44
700-0008-40	ECG Combiner Cable, 6+4-Lead, GE, AAMI	ea	45
700-0008-41	ECG Combiner Cable, 3-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-42	ECG Combiner Cable, 5-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-43	ECG Combiner Cable, 6-Lead, Medtronic LifePak 90 angle, AAMI	ea	45
700-0008-44	ECG Combiner Cable, 3-Lead, Nihon Kohden, AAMI	ea	45
700-0008-45	ECG Combiner Cable, 5-Lead, Nihon Kohden, AAMI	ea	45
700-0008-46	ECG Combiner Cable, 6-Lead, Nihon Kohden, AAMI	ea	45
700-0008-47	ECG Combiner Cable, 3-Lead, Philips, AAMI	ea	45
700-0008-48	ECG Combiner Cable, 5+5-Lead, Philips, AAMI	ea	45
700-0008-49	ECG Combiner Cable, 5-Lead, Philips, AAMI	ea	45
700-0008-50	ECG Combiner Cable, 6+4-Lead, Philips, AAMI	ea	45
700-0008-51	ECG Combiner Cable, 6-Lead, Philips, AAMI	ea	45
700-0008-52	ECG Combiner Cable, 3-Lead, Zoll 90 angle, AAMI	ea	45
700-0008-53	ECG Combiner Cable, 5-Lead, Zoll 90 angle, AAMI	ea	45
700-0008-54	ECG Combiner Cable, 6-Lead, Zoll 90 angle, AAMI	ea	46
700-0008-55	ECG Combiner Cable, 3-Lead Spacelabs monitors, AAMI	ea	43
700-0008-56	ECG Combiner Cable, 5+5-Lead Spacelabs monitors, AAMI	ea	43
700-0008-57	ECG Combiner Cable, 5-Lead Spacelabs monitors, AAMI	ea	43
700-0008-58	ECG Combiner Cable, 6+4-Lead, Spacelabs monitors, AAMI	ea	43
700-0008-59	ECG Combiner Cable, 6-Lead, Spacelabs monitors, AAMI	ea	43
700-0008-60	ECG Combiner adapter, 5-Lead, Datascope Telemetry , AAMI	ea	46
700-0008-61	ECG Combiner adapter, 3-Lead, Datascope , AAMI	ea	46
700-0008-62	ECG Combiner adapter, 5-Lead, Datascope , AAMI	ea	46
700-0008-63	ECG Combiner adapter, 3-Lead, Datex-Ohmeda , AAMI	ea	46
700-0008-64	ECG Combiner adapter, 5-Lead, Datex-Ohmeda , AAMI	ea	46
700-0008-65	ECG Combiner adapter, 5-Lead, GE Apex Pro FH Telemetry , AAMI	ea	46
700-0008-66	ECG Combiner adapter, 6-Lead, GE Apex Pro FH Telemetry , AAMI	ea	46
700-0008-67	ECG Combiner adapter, 3-Lead, GE , AAMI	ea	46
700-0008-68	ECG Combiner adapter, 5-Lead, GE , AAMI	ea	46
700-0008-69	ECG Combiner adapter, 6-Lead, GE , AAMI	ea	47
700-0008-70	ECG Combiner adapter, 3-Lead, Nihon Kohden , AAMI	ea	47
700-0008-71	ECG Combiner adapter, 6-Lead, Nihon Kohden , AAMI	ea	47
700-0008-72	ECG Combiner adapter, 3/5-Lead, Philips Telemetry , AAMI	ea	47
700-0008-73	ECG Combiner adapter, 3-Lead, Philips , AAMI	ea	47
700-0008-74	ECG Combiner adapter, 5-Lead, Philips , AAMI	ea	47
700-0008-75	ECG Combiner adapter, 3-Lead, Siemens Draeger , AAMI	ea	47
700-0008-76	ECG Combiner adapter, 5-Lead, Siemens Draeger , AAMI	ea	47
700-0008-78	ECG Combiner adapter, 3-Lead Spacelabs Monitor compatible, AAMI	ea	43
700-0008-79	ECG Combiner adapter, 5-Lead Spacelabs Monitor compatible, AAMI	ea	43
700-0008-80	ECG Combiner adapter, 5-Lead Spacelabs Telemetry compatible, AAMI	ea	43,56
700-0008-81	ECG Combiner Cable, 3-Lead Spacelabs monitors, IEC	ea	60
700-0008-83	ECG Combiner Cable, 5+5-Lead Spacelabs monitors, IEC	ea	60

Order Code	Description	UOM	Page(s)
700-0008-84	ECG Combiner Cable, 5-Lead Spacelabs monitors, IEC	ea	60
700-0008-88	ECG Combine adapter, Aria Telemetry, 96281, AAMI	ea	43,56
700-0010-00	ECG Combiner Lead Wire Set, 3L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-01	ECG Combiner Lead Wire Set, 5L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-02	ECG Combiner Lead Wire Set, 5L , multi-pinch, AAMI, 130 cm/51 in	1 Set	41
700-0010-03	ECG Combiner Radiolucent Lead Wire Set, multi-pinch, AAMI, 74 cm/29 in	1 Set	42
700-0010-04	ECG Combiner Radiolucent Lead Wire Set, 5L , multi-pinch, AAMI, 74 cm/29 in	1 Set	42
700-0010-05	ECG Combiner Radiolucent Lead Wire Set, 3L multi-pinch, AAMI, 122 cm/48 in	1 Set	42
700-0010-06	ECG Combiner Radiolucent Lead Wire Set, 5L multi-pinch, AAMI, 122 cm/48 in	1 Set	42
700-0010-07	ECG Combiner Radiolucent Lead Wire Set, 3L multi-pinch, AAMI, 183 cm/72 in	1 Set	42
700-0010-08	ECG Combiner Radiolucent Lead Wire Set, 5L multi-pinch, AAMI, 183 cm/72 in	1 Set	42
700-0010-09	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, black	5/pk	55
700-0010-10	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, brown	5/pk	55
700-0010-11	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, green	5/pk	55
700-0010-12	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, red	5/pk	55
700-0010-13	ECG Lead Wires, DIN unshielded, snap, AAMI, 74 cm/29 in, white	5/pk	55
700-0010-14	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI,102 cm/ 40 in	1 Set	42
700-0010-15	ECG Combiner Lead Wire Set, 5L, snap, AAMI, 61 cm/24 in	1 Set	42
700-0010-16	ECG Combiner Lead Wire Set, 5L, multi-pinch, AAMI, 61 cm/24 in	1 Set	42
700-0014-00	TruLink SpO2 Sensor Adapter Cable for Tele transmitter	ea	58, 68
700-0015-00	TruLink ABP Interface Cable, 90343/91343/96217, 60 cm/24 in	ea	59, 69
700-0027-00	CO2 sampling line, sidestream, 91517	ea	88
700-0028-00	INVBP interface adapter cable, dual pressure	ea	86
700-0029-00	Novamatrix sensor adapter cable	ea	78
700-0030-00	Trulink SPO2 cable	ea	77, 83, 92
700-0031-00	Dual temperature adapter cable, 30 cm/12 in	ea	84
700-0287-00	Trulink SPO2 Cable, elance, 305 cm/10 ft	ea	111
700-0290-00	SPO2 Adapter Cable, elance, Nellcor OxiMax, 305 cm / 10 ft	ea	113
700-0292-01	INVBP integrated adapter cable, Edwards Lifesciences, dual pressure	ea	86
700-0293-01	INVBP integrated adapter cable, Edwards Lifesciences, single pressure	ea	86
700-0294-01	INVBP integrated adapter cable, Abbott Labs/Hospira/ICU Medical, dual pressure	ea	86
700-0295-01	Single pressure integrated cable	ea	86
700-0316-00	Patient cable, 3-Channel, 6-lead, 22 in	ea	100
700-0317-00	Patient cable, 3-Channel, 3-lead, Belt Clip, 41 in	ea	100
700-0318-00	Patient cable, 3-channel	ea	100
700-0319-00	Patient cable, 3-Channel, 3-lead, Neck Pouch, 26 in	ea	100
700-0320-00	Patient cable, 2-Channel	ea	100
700-0321-00	Patient cable, 2-Channel, short	ea	100
700-0322-00	Patient cable, 3-Channel, 7-lead, 23 in	ea	100
700-0323-00	Patient cable, 2-Channel, 4-lead, Belt Clip, long	ea	100
700-0324-00	Patient cable, 2-Channel, 4-lead, Neck Pouch, short	ea	100
700-0325-00	Lifecard 12 varios yoke	ea	101
700-0326-00	Lifecard 12 varios yoke, including 12-lead patient cable, IEC	ea	101
700-0327-00	Lifecard 12 varios yoke, including 12-lead patient cable, AHA	ea	101
700-0328-00	Patient cable, Lifecard 12, IEC	ea	101
700-0329-00	Patient cable, Lifecard 12, AHA	ea	101
700-0330-00	Patient cable,long, Lifecard 12, IEC	ea	101
700-0331-00	Patient cable, long, Lifecard 12, AHA	ea	101
700-0791-00	INVBP interface adapter cable, single pressure	ea	86
700-0792-00	Nellcor OxiMax sensor adapter cable	ea	83, 95
700-0906-01	Masimo SET SPO2 sensor adapter cables	ea	79, 93
700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible	ea	84
701-0122-00	Temperature probe covers	200/bx	95
701-0122-00	Disposable Temperature Probe Covers	ea	95
701-4000-00	Skin temperture sensor	20/bx	85

Order Code	Description	UOM	Page(s)
701-4009-00	Esophageal/stethoscope temperature sensor, 9FR	20/bx	85
701-4012-00	Esophageal/stethoscope temperature sensor, 12FR	20/bx	85
701-4018-00	Esophageal/stethoscope temperature sensor, 18FR	20/bx	85
701-4024-00	Esophageal/stethoscope temperature sensor, 24 FR	20/bx	85
701-4109-00	Esophageal/rectal temperature probe	20/bx	85
701-4112-00	Esophageal/stethoscope temperature sensor, 12 FR	20/bx	85
701-4214-00	Foley catheter w/temperature probe, 14FR	20/bx	85
701-4216-00	Foley catheter w/temperature probe, 16FR	20/bx	85
701-4218-00	Foley catheter w/temperature probe, 18FR	20/bx	85
701-4300-00	Tympanic temperature sensor, adult	20/bx	85
704-0001-00	CO2 mainstream sensor, 91517	ea	74
704-0002-00	CO2 airway adapter, mainstream, 91517, reusable, adult	ea	74
704-0003-00	CO2 airway adapter, mainstream, 91517, reusable, neonate	ea	74
704-0004-00	CO2 airway adapter, mainstream, 91517, single use, neonate	24/bx	74
704-0005-00	CO2 airway adapter, mainstream, 91517, single use, adult	24/bx	74
704-0006-00	Cable holding clip	5/pk	74
704-0010-00	CO2 nasal cannula, 91517, adult	10/bx	72
704-0011-00	CO2 nasal cannula, 91517, infant	10/bx	72
704-0012-00	CO2 nasal cannula, 91517, pediatric	10/bx	72
704-0014-00	CO2 nasal cannula w/O2 delivery, 91517, adult	10/bx	72
704-0015-00	CO2 nasal cannula w/O2 delivery, 91517, infant/neonatal	10/bx	72
704-0016-00	CO2 nasal cannula w/O2 delivery, 91517, pediatric	10/bx	72
704-0017-00	CO2 oral/nasal cannula, 91517, adult	10/bx	72
704-0019-00	CO2 oral/nasal cannula, 91517, pediatric	10/bx	72
704-0020-00	CO2 oral/nasal cannula w/O2 delivery, 91517, adult	10/bx	72
704-0022-00	CO2 oral/nasal cannula w/O2 delivery, 91517, pediatric	10/bx	72
704-0023-00	Airway adapter + sampling line, sidestream, 91517, adult/pediatric	10/bx	73
704-0024-00	Airway adapter + sampling line, sidestream, 91517, pediatric/infant	10/bx	73
704-0025-00	Airway adapter + sampling line w/nafion, sidestream, 91517, adult/ped	10/bx	73
704-0026-00	Airway adapter + sampling line w/nafion, sidestream, 91517, ped/infant	10/bx	73
704-0027-00	CO2 sampling line, sidestream, 91517	10/bx	73
704-0028-00	CO2 sampling line w/nafion, sidestream, 91517	10/bx	73
704-0167-00	FilterLine H Set, Adult/Pediatric	25/pk	114
704-0168-00	FilerLine Set, Adult/Pediatric	25/pk	114
704-0170-00	Smart CapnoLine, Pediatric	25/pk	114
704-0171-00	Smart CapnoLine Plus, Adult	25/pk	114
704-0173-00	Airway adapter, mainstream, 92517, adult/pediatric	25/bx	71
704-0174-00	Airway adapter, mainstream, 92517, infant	10/bx	71
704-0175-00	Airway T-adapter, sidestream, FLL, adult/pediatric	25/bx	71
704-0176-00	CO2 nasal cannula, MLL, adult, 2.1m/7 ft	25/cs	70
704-0177-00	CO2 nasal cannula, MLL, pediatric, 2.1m/7 ft	25/cs	70
704-0178-00	CO2 nasal cannula, infant, 2.1m/7 ft	25/cs	70
704-0179-00	CO2 nasal cannula w/O2 delivery, MLL, infant, 2.1m/7 ft	25/cs	70
704-0180-00	CO2 nasal cannula w/O2 delivery, MLL, pediatric, 2.1m/7 ft	25/cs	70
704-0181-00	CO2 nasal cannula w/O2 delivery, MLL, adult, 2.1m/7 ft	25/cs	70
704-0182-00	CO2 oral/nasal cannula w/O2, delivery, MLL, pediatric, 2.1m/7 ft	25/cs	70
704-0183-00	CO2 oral/nasal cannula w/O2 delivery, MLL, adult, 2.1m/7 ft	25/cs	70
704-0184-00	CO2 nasal cannula, MML, adult, 3m/10 ft	25/cs	70
704-0185-00	CO2 nasal cannula, MML, pediatric, 3m/10 ft	25/cs	70
704-0186-00	CO2 nasal cannula w/O2 delivery, MLL, adult, 3m/10 ft	25/cs	70
704-0187-00	CO2 oral/nasal cannula w/O2 delivery, MLL, adult, 3m/10 ft	25/cs	70
704-0189-00	CO2 oral/nasal cannula w/O2 delivery, MLL, pediatric, 3m/10 ft	25/cs	70
704-0191-00	CO2 sampling line w/T-adapter, MML, adult, 2.5 m/6 ft	25/bx	71
704-0192-00	CO2 sampling line w/T-adapter, MML, adult, 3 m/10 ft	25/bx	71
704-0193-00	CO2 emergency nasal cannula, adult	25/cs	71

Order Code	Description	UOM	Page(s)
706-0182-00	Disposable cuff cover for reusable ABP cuffs, standard	50/bx	96
706-0183-00	Disposable cuff cover for reusable ABP cuffs, large	50/bx	96
706-0189-00	Disposable pouch for evo recorder	100/pk	101
706-0190-00	ECG Lead Wire Grouper, 5 Lead, 96281, IEC	5/bg	66
706-0191-00	ECG Lead Wire Grouper, 5 Lead, 96281, AAMI	5/bg	56
706-0200-00	Disposable Pouch, 90227	50/cs	97
712-1039-00	Quick Release Metal Fittings, female bayonet	3/bg	97
714-0018-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9, 91, 108
714-0018-12	BP hose, dual tube to single tube, adult	ea	9
714-0019-01	BP hose, single tube, neonatal	ea	5, 8, 11
714-0019-12	BP hose, dual tube to single tube, neonatal	ea	11
714-0020-00	BP cuff, Spacelabs Nylon, 1T, child, 12-19 cm, MHP	ea	6
714-0021-00	BP cuff, Spacelabs Nylon, 1T, sm adult, 17-25 cm, MHP	ea	6
714-0023-00	BP cuff, Spacelabs Nylon, 1T, adult, 23-33 cm, MHP	ea	6
714-0025-00	BP cuff, Spacelabs Nylon, 1T, lg adult, 31-40 cm, MHP	ea	6
714-0027-00	BP cuff, Spacelabs Nylon, 1T, thigh, 38-50 cm, MHP	ea	6
714-0038-00	BP Slip Fit Connector, 3mm ID Adult NIBP Hose	ea	9
714-0043-00	BP TruLink dual hose conversion kit, adult hose assembly	ea	10
714-0045-00	BP Slip Fit Connector, 6mm ID Adult NIBP Hose	ea	9
714-0046-00	BP Slip Fit Connector, 5mm ID Adult NIBP Hose	ea	9
714-0047-00	BP Slip Fit Connector, 4mm ID Adult NIBP Hose	ea	9
714-0052-50	BP hose fitting, plastic, neonatal, 50/pk	50/pk	11
714-0053-10	BP hose fitting, neonatal	10/pk	11
714-0054-50	BP hose fitting, plastic, infant, 50/pk	50/pk	11
714-0056-00	BP Y adapter, adult, 3 mm/1/8 in, barb	5/pk	10
714-0057-00	BP Y adapter, dual to single hose, adult	5/pk	10
714-0057-10	BP hose hose connector, non-luer	10/pk	11
714-0060-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9, 91
714-0061-00	BP hose, single tube, adult	ea	1, 2, 3, 4, 6, 7, 9
714-0062-00	BP hose, single tube, neonatal	ea	5, 8, 11
714-0070-00	BP cuff, Spacelabs Nylon, 1T, child long, 12-19 cm, MHP	ea	6
714-0071-00	BP cuff, Spacelabs Nylon, 1T, sm adult long, 17-25 cm, MHP	ea	6
714-0073-00	BP cuff, Spacelabs Nylon, 1T, adult long, 23-33 cm, MHP	ea	6
714-0075-00	BP cuff, Spacelabs Nylon, 1T, lg adult long, 31-40 cm, MHP	ea	6
714-1020-00	BP cuff, Spacelabs Vinyl, 1T, child, 12-19 cm, HP	5/pk	7
714-1021-00	BP cuff, Spacelabs Vinyl, 1T, small adult, 17-25 cm, HP	5/pk	7
714-1023-00	BP cuff, Spacelabs Vinyl, 1T, adult, 23-33 cm, HP	5/pk	7
714-1025-00	BP cuff, Spacelabs Vinyl, 1T, lg adult, 31-40 cm, HP	5/pk	7
714-1027-00	BP cuff, Spacelabs Vinyl, 1T, thigh, 38-50 cm, HP	5/pk	7
714-1028-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 1, 3-6 cm, M316	5/pk	8
714-1029-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 2, 4-8 cm, M316	5/pk	8
714-1030-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 3, 6-11cm, M316	5/pk	8
714-1032-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 4, 7-13 cm, M316	5/pk	8
714-1033-01	BP cuff, Spacelabs Soft Neonatal, 1T, size 5, 8-15 cm, M316	5/pk	8
714-1073-00	BP cuff, Spacelabs Vinyl, 1T, adult long, 23-33 cm, HP	5/pk	7
714-1075-00	BP cuff, Spacelabs Vinyl, 1T, lg adult long, 31-40 cm, HP	5/pk	7
714-4001-00	BP cuff, UltraCheck Nylon, 1T, child, 13-20cm,MHP	ea	108
714-4002-00	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26cm,MHP	ea	108
714-4003-00	BP cuff, UltraCheck Nylon, 1T, adult, 26-35cm,MHP	ea	108
714-4004-00	BP cuff, UltraCheck Nylon, 1T, adult long,29-38cm,MHP	ea	108
714-4005-00	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42cm,MHP	ea	108
714-4006-00	BP cuff, UltraCheck Nylon, 1T, lg adult long, 35-44cm,MHP	ea	108
714-4007-00	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50cm,MHP	ea	108
715-1121-10	BP cuff, TruLink Nylon, 1T, child, 13-20 cm, HP	ea	1, 91
715-1122-10	BP cuff, TruLink Vinyl, 1T, child, 13-20 cm, HP	5/bg	3, 91

Order Code	Description	UOM	Page(s)
715-1123-10	BP cuff, TruLink Soft, 1T, child, 13-20 cm, HP	5/bg	2, 91
715-1131-10	BP cuff, TruLink Nylon, 1T, sm adult, 18-26 cm, HP	ea	1, 91
715-1132-10	BP cuff, TruLink Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	3, 91
715-1133-10	BP cuff, TruLink Soft, 1T, sm adult, 18-26 cm, HP	5/bg	2, 91
715-1135-10	BP cuff, TruLink Yellow Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	4
715-1141-10	BP cuff, TruLink Nylon, 1T, adult, 26-36 cm, HP	ea	1, 91
715-1142-10	BP cuff, TruLink Vinyl, 1T, adult, 26-35 cm, HP	5/bg	3, 91
715-1143-10	BP cuff, TruLink Soft, 1T, adult, 26-35 cm, HP	5/bg	2, 91
715-1145-10	BP cuff, TruLink Yellow Vinyl, 1T, adult, 26-35 cm, HP	5/bg	4
715-1151-10	BP cuff, TruLink Nylon, 1T, adult long, 23-33 cm, HP	ea	1, 91
715-1152-10	BP cuff, TruLink Vinyl, 1T, adult long, 23-33 cm, HP	5/bg	3, 91
715-1153-10	BP cuff, TruLink Soft, 1T, adult long, 23 - 33 cm, HP	5/bg	2, 91
715-1161-10	BP cuff, TruLink Nylon, 1T, lg adult, 32-42 cm, HP	ea	1, 91
715-1162-10	BP cuff, TruLink Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	3, 91
715-1163-10	BP cuff, TruLink Soft, 1T, lg adult, 32-42 cm, HP	5/bg	2, 91
715-1165-10	BP cuff, TruLink Yellow Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	4
715-1171-10	BP cuff, TruLink Nylon, 1T, lg adult long, 35-44 cm, HP	ea	1, 91
715-1172-10	BP cuff, TruLink Vinyl, 1T, lg adult long, 35-44 cm, HP	5/bg	3, 91
715-1173-10	BP cuff, TruLink Soft, 1T, lg adult long, 35-44 cm, HP	5/bg	2, 91
715-1181-10	BP cuff, TruLink Nylon, 1T, thigh, 42-50 cm, HP	ea	1, 91
715-1182-10	BP cuff, TruLink Vinyl, 1T, thigh, 42-50 cm, HP	5/bg	3, 91
715-1183-10	BP cuff, TruLink Soft, 1T, thigh, 42-50 cm, HP	5/bg	2, 91
715-2114-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	5
715-2124-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	5
715-2134-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	5
715-2144-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	5
715-2154-10	BP cuff, TruLink Non-PVC Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	5
717-0071-00	ECG Electrode, AMBU	5/pk	100
717-0072-00	ECG electrodes, foam with snap, adult	25/pk	104
717-0408-00	Exhaust scavenging line w/female luer connection on both ends (90518)	ea	77
717-0447-00	D-FEND Water Trap Cartridge	10/pk	76
717-0448-00	D-FEND+ Water Trap Cartridge	10/pk	76
719-0001-10	BISX Pediatric Sensor	25/bx	90
719-0002-10	BISX Quarto Sensor	25/bx	90
719-0003-10	BISX Extend Sensor	25/bx	90
719-0004-01	BISX Patient Interface cable	ea	90
719-0006-20	BISx Pod	ea	90
719-0007-01	BISX cable	ea	90
803FGA-24	Unifusor, 500mL w/aneroid gauge and thumbwheel valve	24/bx	37
803SGA-24	Unifusor, 500mL w/aneroid gauge and stopcock valve	24/bx	37
804FGA-24	Unifusor, 1000mL w/aneroid gauge and thumbwheel valve	24/bx	37
804SGA-24	Unifusor, 1000mL w/aneroid gauge and stopcock valve	24/bx	37
805FGA-10	Unifusor, 3000mL w/aneroid gauge and thumbwheel valve	10/bx	37
805SGA-10	Unifusor, 3000mL w/aneroid gauge and stopcock valve	10/bx	37
899F-04	Aneroid gauge and thumbwheel valve assy	4/bx	38,39
899S-04	Aneroid gauge and stopcock valve assy	4/bx	38,39
903SGA-24	Unifusor, 500mL w/piston gauge and stopcock valve	24/bx	37
903TGA-24	Unifusor, 500mL w/piston gauge and thumbwheel valve	24/bx	37
904SGA-24	Unifusor, 1000mL w/piston gauge and stopcock valve	24/bx	37
904TGA-24	Unifusor, 1000mL w/piston gauge and thumbwheel valve	24/bx	37
905SGA-10	Unifusor, 3000mL w/piston gauge and stopcock valve	10/bx	37
905TGA-10	Unifusor, 3000mL w/piston gauge and thumbwheel valve	10/bx	37
98-0002-50	Microphone pads, 5 PK	5/pk	103
98-0006-00	K-Sound microphone	ea	103
98-0061-01	Cuff, Orbit-K, small adult	ea	103

Order Code	Description	UOM	Page(s)
999S-04	Piston gauge and stopcock valve assembly	4/bx	38,39
999T-04	Piston gauge and thumbwheel valve assembly	4/bx	38,39
ADP-182YBPF	BP Y adapter, 182YBPF	ea	35
ADP-182YRX4	BP Y adapter, 182YRX4	ea	35
ADP-5322YBPF	BP Y adapter, 5322YBPF	ea	35
ADP-5322YRX4	BP Y adapter, 5322YRX4	ea	35
ADP-DM2YRX4	BP Y adapter, Metal Female HP to 2 DM	ea	35
ADP-MDM2YMHP	BY Y adapter, Metal DM Connector and Metal HP Fitting	ea	35
BRDT3854-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, no fitting	5/bx	25
BRDT3854DM-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, DM	5/bx	25
BRDT3854FL-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, FL	5/pk	25
BRDT3854HP-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, HP	5/bx	25
BRDT3854MB-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, manual	5/bx	25
BRDT3854ML-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, ML	5/pk	25
BRDT3854MQ-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQ	5/bx	25
BRDT3854MQMF-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQMF	5/bx	25
BRDT3854MQML-05	BP cuff, SoftCheck Curve, 2T, XL-adult, 38-54 cm, MQML	5/pk	25
BRST3854-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, DM	5/bx	25
BRST3854DM-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, FL	5/bx	25
BRST3854FL-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, HP	5/bx	25
BRST3854HP-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, ML	5/bx	25
BRST3854ML-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQ	5/bx	25
BRST3854MQ-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQML	5/bx	25
BRST3854MQML-05	BP cuff, SoftCheck Curve, 1T, XL-adult, 38-54 cm, MQML	5/bx	25
BRUD3854	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, no fitting	ea	18
BRUD3854DM	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, DM	ea	18
BRUD3854HP	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, HP	ea	18
BRUD3854MB	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, manual	ea	18
BRUD3854ML	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, ML	ea	18
BRUD3854MQ	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, MQ	ea	18
BRUD3854MQMF	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, MQMF	ea	18
BRUD3854MQML	BP cuff, Ultracheck Curve, 2T, XL-adult, 38-54 cm, MQML	ea	18
BRUS3854	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, no fitting	ea	18
BRUS3854DM	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, DM	ea	18
BRUS3854FL	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, FL	ea	18
BRUS3854HP	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, HP	ea	18
BRUS3854MB	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, manual	ea	18
BRUS3854MHP	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, MHP	ea	18
BRUS3854ML	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, ML	ea	18
BRUS3854MQ	BP cuff, Ultracheck Curve, 1T, XL-adult, 38-54 cm, MQ	ea	18
CFC256MB	Flashcard LifeCard CF, 256MB	ea	102
CFC90MB	Flashcard Memory, 90MByte	ea	102
DT0814-05	BP cuff, SoftCheck Fabric, 2T infant, 8-14 cm, no fitting	5/bg	24
DT0814DM-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, DM	5/bg	22
DT0814DM-20	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, DM	20/bx	22
DT0814MB-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, manual	5/bg	24
DT0814MQ-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQ	5/bg	23
DT0814MQ-20	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQ	20/bx	23
DT0814MQMF-05	BP cuff, SoftCheck Fabric, 2T, infant, 8-14 cm, MQMF	5/bg	24
DT1320-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, no fitting	5/bg	24
DT1320DM-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, DM	5/bg	22
DT1320DM-20	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, DM	20/bx	22
DT1320MB-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, manual	5/bg	24
DT1320MQ-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQ	5/bg	23
DT1320MQ-20	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQ	20/bx	23

Order Code	Description	UOM	Page(s)
DT1320MQMF-05	BP cuff, SoftCheck Fabric, 2T, child, 13-20 cm, MQMF	5/bg	24
DT1826-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, no fitting	5/bg	24
DT1826DM-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, DM	5/bg	22
DT1826DM-20	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, DM	20/bx	22
DT1826MB-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, manual	5/bg	24
DT1826MQ-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQ	5/bg	23
DT1826MQ-20	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQ	20/bx	23
DT1826MQMF-05	BP cuff, SoftCheck Fabric, 2T, sm adult, 18-26 cm, MQMF	5/bg	24
DT2635-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, no fitting	5/bg	24
DT2635DM-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, DM	5/bg	22
DT2635DM-20	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, DM	20/bx	22
DT2635MB-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, manual	5/bg	24
DT2635MQ-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQ	5/bg	23
DT2635MQ-20	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQ	20/bx	23
DT2635MQMF-05	BP cuff, SoftCheck Fabric, 2T, adult, 26-35 cm, MQMF	5/bg	24
DT2938-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, no fitting	5/bg	24
DT2938DM-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, DM	5/bg	22
DT2938DM-20	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, DM	20/bx	22
DT2938MB-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, manual	5/bg	24
DT2938MQ-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQ	5/bg	23
DT2938MQ-20	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQ	20/bx	23
DT2938MQMF-05	BP cuff, SoftCheck Fabric, 2T, adult long, 29-38 cm, MQMF	5/bg	24
DT3242-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, no fitting	5/bg	24
DT3242DM-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, DM	5/bg	22
DT3242DM-20	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, DM	20/bx	22
DT3242MB-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, manual	5/bg	24
DT3242MQ-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQ	5/bg	23
DT3242MQ-20	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQ	20/bx	23
DT3242MQMF-05	BP cuff, SoftCheck Fabric, 2T, lg adult, 32-42 cm, MQMF	5/bg	24
DT3544-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, no fitting	5/bg	24
DT3544DM-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, DM	5/bg	22
DT3544DM-20	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, DM	20/bx	22
DT3544MB-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, manual	5/bg	24
DT3544MQ-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQ	5/bg	23
DT3544MQ-20	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQ	20/bx	23
DT3544MQMF-05	BP cuff, SoftCheck Fabric, 2T, lg adult long, 35-44 cm, MQMF	5/bg	24
DT4250-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, no fitting	5/bg	24
DT4250DM-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, DM	5/bg	22
DT4250DM-20	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, DM	20/bx	22
DT4250MB-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, manual	5/bg	24
DT4250MQ-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQ	5/bg	23
DT4250MQ-20	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQ	20/bx	23
DT4250MQMF-05	BP cuff, SoftCheck Fabric, 2T, thigh, 42-50 cm, MQMF	5/bg	24
M1050-4	Reusable Y patient circuit, adult	ea	107
MS1R-106902-A0	ECG cable, 4MM, banana connector, IEC	ea	104
MS1R-107048-A0	ECG cable, 4MM, banana connector, AHA	ea	104
MS1R-107371-A0	Recording paper, Z-fold, 210 MM X 295 MM	ea	104
MS1R-107451-A0	Recording paper, Z-fold, 215MM x 280MM	ea	104
MS1R-107560-A0	Recording paper, Z-fold, 110MM*140MM	ea	104
MS1R-107581-A0	ECG cable, snap connector, IEC	ea	104
MS1R-107582-A0	ECG cable, snap connector, AHA	ea	104
MS1R-107583-A0	ECG cable, 4MM, graber connection, IEC	ea	104
MS1R-107584-A0	ECG cable, 4MM, graber connection, AHA	ea	104
MS1R-32461	Recording paper, Roll, 210 MM x 30M	ea	104
MS1R-32462	Recording paper, Z-Fold, 210MM x 295MM	ea	104

Order Code	Description	UOM	Page(s)
MS1R-32463	Recording paper, Roll, 215 MM x 30M	ea	104
MS4R-111897	Paper roller, plastic, SL12	ea	104
NN1DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	31
NN1ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	31
NN1ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	31
NN1ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	31
NN2DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	31
NN2ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	31
NN2ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	31
NN2ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	31
NN3DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	31
NN3ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	31
NN3ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	31
NN3ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	31
NN4DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	31
NN4ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	31
NN4ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	31
NN4ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	31
NN5DT-10	BP cuff, SoftCheck Fabric Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	31
NN5ST-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	31
NN5ST-HP-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	31
NN5ST-M316-10	BP cuff, SoftCheck Fabric Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	31
PK0001	Fabric waist belt	ea	101
PNN1DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	34
PNN1ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	34
PNN1ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	34
PNN1ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	34
PNN2DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	34
PNN2ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	34
PNN2ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	34
PNN2ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	34
PNN3DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	34
PNN3ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	34
PNN3ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	34
PNN3ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	34
PNN4DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	34
PNN4ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	34
PNN4ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	34
PNN4ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	34
PNN5DT-10	BP cuff, SoftCheck Free Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	34
PNN5ST-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	34
PNN5ST-HP-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	34
PNN5ST-M316-10	BP cuff, SoftCheck Free Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	34
ST0814-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, no fitting	5/bg	21
ST0814DM-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, DM	5/bg	20
ST0814FL-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, FL	5/bg	21
ST0814HP-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, HP	5/bg	19
ST0814HP-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, HP	20/bx	19
ST0814ML-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, ML	5/bg	19
ST0814ML-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, ML	20/bx	19
ST0814MQ-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQ	5/bg	20
ST0814MQ-20	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQ	20/bx	20
ST0814MQML-05	BP cuff, SoftCheck Fabric, 1T, infant, 8-14 cm, MQML	5/bg	21
ST1320-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, no fitting	5/bg	21
ST1320DM-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, DM	5/bg	20

Order Code	Description	UOM	Page(s)
ST1320FL-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, FL	5/bg	21
ST1320HP-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, HP	5/bg	19
ST1320HP-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, HP	20/bx	19
ST1320ML-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, ML	5/bg	19
ST1320ML-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, ML	20/bx	19
ST1320MQ-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQ	5/bg	20
ST1320MQ-20	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQ	20/bx	20
ST1320MQML-05	BP cuff, SoftCheck Fabric, 1T, child, 13-20 cm, MQML	5/bg	21
ST1826-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, no fitting	5/bg	21
ST1826DM-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, DM	5/bg	20
ST1826FL-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, FL	5/bg	21
ST1826HP-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, HP	5/bg	19
ST1826HP-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, HP	20/bx	19
ST1826ML-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, ML	5/bg	19
ST1826ML-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, ML	20/bx	19
ST1826MQ-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQ	5/bg	20
ST1826MQ-20	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQ	20/bx	20
ST1826MQML-05	BP cuff, SoftCheck Fabric, 1T, sm adult, 18-26 cm, MQML	5/bg	21
ST2111	2.0 L Rebreathing bag, 22mm neck	ea	107
ST2635-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, no fitting	5/bg	21
ST2635DM-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, DM	5/bg	20
ST2635FL-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, FL	5/bg	21
ST2635HP-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, HP	5/bg	19
ST2635HP-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, HP	20/bx	19
ST2635ML-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, ML	5/bg	19
ST2635ML-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, ML	20/bx	19
ST2635MQ-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQ	5/bg	20
ST2635MQ-20	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQ	20/bx	20
ST2635MQML-05	BP cuff, SoftCheck Fabric, 1T, adult, 25-35 cm, MQML	5/bg	21
ST2670	Malleable face mask, size 4	ea	107
ST2938-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, no fitting	5/bg	21
ST2938DM-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, DM	5/bg	20
ST2938FL-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, FL	5/bg	21
ST2938HP-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, HP	5/bg	19
ST2938HP-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, HP	20/bx	19
ST2938ML-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, ML	5/bg	19
ST2938ML-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, ML	20/bx	19
ST2938MQ-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQ	5/bg	20
ST2938MQ-20	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQ	20/bx	20
ST2938MQML-05	BP cuff, SoftCheck Fabric, 1T, adult long, 29-38 cm, MQML	5/bg	21
ST3242-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, no fitting	5/bg	21
ST3242DM-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, DM	5/bg	20
ST3242FL-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, FL	5/bg	21
ST3242HP-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, HP	5/bg	19
ST3242HP-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, HP	20/bx	19
ST3242ML-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, ML	5/bg	19
ST3242ML-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, ML	20/bx	19
ST3242MQ-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQ	5/bg	20
ST3242MQ-20	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQ	20/bx	20
ST3242MQML-05	BP cuff, SoftCheck Fabric, 1T, lg adult, 32-42 cm, MQML	5/bg	21
ST3544-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, no fitting	5/bg	21
ST3544DM-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, DM	5/bg	20
ST3544FL-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, FL	5/bg	21
ST3544HP-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, HP	5/bg	19
ST3544HP-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, HP	20/bx	19

Order Code	Description	UOM	Page(s)
ST3544ML-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, ML	5/bg	19
ST3544ML-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, ML	20/bx	19
ST3544MQ-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQ	5/bg	20
ST3544MQ-20	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQ	20/bx	20
ST3544MQML-05	BP cuff, SoftCheck Fabric, 1T, lg adult long, 35-44 cm, MQML	5/bg	21
ST4250-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, no fitting	5/bg	21
ST4250DM-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, DM	5/bg	20
ST4250FL-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, FL	5/bg	21
ST4250HP-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, HP	5/bg	19
ST4250HP-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, HP	20/bx	19
ST4250ML-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, ML	5/bg	19
ST4250ML-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, ML	20/bx	19
ST4250MQ-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQ	5/bg	20
ST4250MQ-20	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQ	20/bx	20
ST4250MQML-05	BP cuff, SoftCheck Fabric, 1T, thigh, 42-50 cm, MQML	5/bg	21
ST7557	Magill breathing circuit w/face mask	ea	107
UD0814-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, no fitting	5/bg	17
UD0814DM-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, DM	5/bg	16
UD0814HP-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, HP	5/bg	16
UD0814MB-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, manual	5/bg	17
UD0814MQ-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, MQ	5/bg	16
UD0814MQMF-05	BP cuff, UltraCheck Nylon, 2T, infant, 8-14 cm, MQMF	5/bg	17
UD1320-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, no fitting	5/bg	17
UD1320DM-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, DM	5/bg	16
UD1320HP-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, HP	5/bg	16
UD1320MB-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, manual	5/bg	17
UD1320MQ-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, MQ	5/bg	16
UD1320MQMF-05	BP cuff, UltraCheck Nylon, 2T, child, 13-20 cm, MQMF	5/bg	17
UD1826-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, no fitting	5/bg	17
UD1826DM-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, DM	5/bg	16
UD1826HP-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, HP	5/bg	16
UD1826MB-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, manual	5/bg	17
UD1826MQ-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, MQ	5/bg	16
UD1826MQMF-05	BP cuff, UltraCheck Nylon, 2T, sm adult, 18-26 cm, MQMF	5/bg	17
UD2635-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, no fitting	5/bg	17
UD2635DM-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, DM	5/bg	16
UD2635HP-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, HP	5/bg	16
UD2635MB-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, manual	5/bg	17
UD2635MQ-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, MQ	5/bg	16
UD2635MQMF-05	BP cuff, UltraCheck Nylon, 2T, adult, 26-35 cm, 2T, MQMF	5/bg	17
UD2938-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, no fitting	5/bg	17
UD2938DM-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, DM	5/bg	16
UD2938HP-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, HP	5/bg	16
UD2938MB-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, manual	5/bg	17
UD2938MQ-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, MQ	5/bg	16
UD2938MQMF-05	BP cuff, UltraCheck Nylon, 2T, adult long, 29-39 cm, 2T, MQMF	5/bg	17
UD3242-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, no fitting	5/bg	17
UD3242DM-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, DM	5/bg	16
UD3242HP-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, HP	5/bg	16
UD3242MB-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, manual	5/bg	17
UD3242MQ-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, MQ	5/bg	16
UD3242MQMF-05	BP cuff, UltraCheck Nylon, 2T, lg adult, 32-42 cm, MQMF	5/bg	17
UD3544-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, no fitting	5/bg	17
UD3544DM-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, DM	5/bg	16
UD3544HP-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, HP	5/bg	16

Order Code	Description	UOM	Page(s)
UD3544MB-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, manual	5/bg	17
UD3544MQ-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, MQ	5/bg	16
UD3544MQMF-05	BP cuff, UltraCheck Nylon, 2T, lg adult long, 35-44 cm, MQMF	5/bg	17
UD4250-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, no fitting	5/bg	17
UD4250DM-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, DM	5/bg	16
UD4250HP-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, HP	5/bg	16
UD4250MB-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50, manual	5/bg	17
UD4250MQ-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, MQ	5/bg	16
UD4250MQMF-05	BP cuff, UltraCheck Nylon, 2T, thigh, 42-50 cm, MQMF	5/bg	17
US0814-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, no fitting	5/bg	15
US0814DM-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, DM	5/bg	13
US0814FL-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, FL	5/bg	14
US0814HP	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, HP	ea	12
US0814HP-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, HP	5/bg	12
US0814MB	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, manual	ea	15
US0814MHP-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MHP	5/bg	12
US0814ML-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, ML	5/bg	13
US0814MML-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MML	5/bg	13
US0814MQ-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MQ	5/bg	14
US0814MQML-05	BP cuff, UltraCheck Nylon, 1T, infant, 8-14 cm, MQML	5/bg	14
US1320-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, no fitting	5/bg	15
US1320DM-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, DM	5/bg	13
US1320FL-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, FL	5/bg	14
US1320HP	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, HP	ea	12
US1320HP-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, HP	5/bg	12
US1320MB	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, manual with sphyg	ea	15
US1320MHP-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MHP	5/bg	12
US1320ML-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, ML	5/bg	13
US1320MML-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MML	5/bg	13
US1320MQ-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MQ	5/bg	14
US1320MQML-05	BP cuff, UltraCheck Nylon, 1T, child, 13-20 cm, MQML	5/bg	14
US1826-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, no fitting	5/bg	15
US1826DM-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, DM	5/bg	13
US1826FL-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, FL	5/bg	14
US1826HP	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, HP	ea	12
US1826HP-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, HP	5/bg	12
US1826MB	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, manual	ea	15
US1826MHP-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MHP	5/bg	12
US1826ML-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, ML	5/bg	13
US1826MML-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MML	5/bg	13
US1826MQ-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MQ	5/bg	14
US1826MQML-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, 1T, MQML	5/bg	14
US2635-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, no fitting	5/bg	15
US2635DM-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, DM	5/bg	13
US2635FL-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, FL	5/bg	14
US2635HP	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, HP	ea	12
US2635HP-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, HP	5/bg	12
US2635MB	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, manual with sphyg	ea	15
US2635MHP-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MHP	5/bg	12
US2635ML-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, ML	5/bg	13
US2635MML-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MML	5/bg	13
US2635MQ-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MQ	5/bg	14
US2635MQML-05	BP cuff, UltraCheck Nylon, 1T, adult, 26-35 cm, MQML	5/bg	14
US2938-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, no fitting	5/bg	15
US2938DM-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, DM	5/bg	13

Order Code	Description	UOM	Page(s)
US2938FL-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, FL	5/bg	14
US2938HP	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, HP	ea	12
US2938HP-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, HP	5/bg	12
US2938MB	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, manual with sphyg	ea	15
US2938MHP-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, 1T, MHP	5/bg	12
US2938ML-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, ML	5/bg	13
US2938MML-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, MML	5/bg	13
US2938MQ-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, MQ	5/bg	14
US2938MQML-05	BP cuff, UltraCheck Nylon, 1T, adult long, 29-38 cm, MQML	5/bg	14
US3242-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, no fitting	5/bg	15
US3242DM-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, DM	5/bg	13
US3242FL-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, FL	5/bg	14
US3242HP	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, HP	ea	12
US3242HP-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, HP	5/bg	12
US3242MB	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, manual with sphyg	ea	15
US3242MHP-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MHP	5/bg	12
US3242ML-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, ML	5/bg	13
US3242MML-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MML	5/bg	13
US3242MQ-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MQ	5/bg	14
US3242MQML-05	BP cuff, UltraCheck Nylon, 1T, lg adult, 32-42 cm, MQML	5/bg	14
US3544-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, no fitting	5/bg	15
US3544DM-05	BP cuff, UltraCheck Nylon, 1T, sm adult, 18-26 cm, DM	5/bg	13
US3544FL-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, FL	5/bg	14
US3544HP	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, HP	ea	12
US3544HP-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, HP	5/bg	12
US3544MB	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, manual	ea	15
US3544MHP-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 35-44 cm, MHP	5/bg	12
US3544ML-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, ML	5/bg	13
US3544MML-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, MML	5/bg	13
US3544MQ-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, MQ	5/bg	14
US3544MQML-05	BP cuff, UltraCheck Nylon, 1T, lg adult long, 29-39 cm, MQML	5/bg	14
US4250-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, no fitting	5/bg	15
US4250DM-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, DM	5/bg	13
US4250FL-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, FL	5/bg	14
US4250HP	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, HP	ea	12
US4250HP-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, HP	5/bg	12
US4250MB	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, manual	ea	15
US4250MHP-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MHP	5/bg	12
US4250ML-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, ML	5/bg	13
US4250MML-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MML	5/bg	13
US4250MQ-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MQ	5/bg	14
US4250MQML-05	BP cuff, UltraCheck Nylon, 1T, thigh, 42-50 cm, MQML	5/bg	14
VDT0814-05	BP cuff, SoftCheck Vinyl, 2T, infant, 8-14 cm, no fitting	5/bg	29
VDT0814DM-05	BP cuff, SoftCheck Vinyl, 2T, infant, 8-14 cm, DM	5/bg	29
VDT1320-05	BP cuff, SoftCheck Vinyl, 2T, child, 13-20 cm, no fitting	5/bg	29
VDT1320DM-05	BP cuff, SoftCheck Vinyl, 2T, child, 13-20 cm, DM	5/bg	29
VDT1826-05	BP cuff, SoftCheck Vinyl, 2T, sm adult, 18-26 cm, no fitting	5/bg	29
VDT1826DM-05	BP cuff, SoftCheck Vinyl, 2T, sm adult, 18-26 cm, DM	5/bg	29
VDT2635-05	BP cuff, SoftCheck Vinyl, 2T, adult, 26-35 cm, no fitting	5/bg	29
VDT2635DM-05	BP cuff, SoftCheck Vinyl, 2T, adult, 26-35 cm, DM	5/bg	29
VDT2938-05	BP cuff, SoftCheck Vinyl, 2T, adult long, 29-38 cm, no fitting	5/bg	29
VDT2938DM-05	BP cuff, SoftCheck Vinyl, 2T, adult long, 29-38 cm, DM	5/bg	29
VDT3242-05	BP cuff, SoftCheck Vinyl, 2T, lg adult, 32-42 cm, no fitting	5/bg	29
VDT3242DM-05	BP cuff, SoftCheck Vinyl, 2T, lg adult, 32-42 cm, DM	5/bg	29
VDT4250-05	BP cuff, SoftCheck Vinyl, 2T, thigh, 42-50 cm, no fitting	5/bg	29

Order Code	Description	UOM	Page(s)
VDT4250DM-05	BP cuff, SoftCheck Vinyl, 2T, thigh, 42-50 cm, DM	5/bg	29
VNN1DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 1, 3-6 cm, ML	10/bx	33
VNN1ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, ML	10/bx	32
VNN1ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, HP	10/bx	32
VNN1ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, M316	10/bx	32
VNN1ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, MQ	10/bx	32
VNN1ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 1, 3-6 cm, SQ	10/bx	32
VNN2DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 2, 4-8 cm, ML	10/bx	33
VNN2ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, ML	10/bx	32
VNN2ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, HP	10/bx	32
VNN2ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, M316	10/bx	32
VNN2ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, MQ	10/bx	32
VNN2ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 2, 4-8 cm, SQ	10/bx	32
VNN3DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 3, 6-11 cm, ML	10/bx	33
VNN3ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, ML	10/bx	32
VNN3ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, HP	10/bx	32
VNN3ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, M316	10/bx	32
VNN3ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, MQ	10/bx	32
VNN3ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 3, 6-11 cm, SQ	10/bx	32
VNN4DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 4, 7-13 cm, ML	10/bx	33
VNN4ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, ML	10/bx	32
VNN4ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, HP	10/bx	32
VNN4ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, M316	10/bx	32
VNN4ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, MQ	10/bx	32
VNN4ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 4, 7-13 cm, SQ	10/bx	32
VNN5DT-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 5, 8-15 cm, ML	10/bx	33
VNN5ST-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, ML	10/bx	32
VNN5ST-HP-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, HP	10/bx	32
VNN5ST-M316-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, M316	10/bx	32
VNN5ST-MQ-10	BP cuff, SoftCheck Vinyl Neonatal, 1T, size 5, 8-15 cm, MQ	10/bx	32
VNN5ST-SQ-10	BP cuff, SoftCheck Vinyl Neonatal, 2T, size 5, 8-15 cm, SQ	10/bx	32
VST0814HP-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, HP	5/bg	26
VST0814HP-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, HP	20/bx	26
VST0814ML-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, ML	5/bg	27
VST0814ML-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, ML	20/bx	27
VST0814MQ-05	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQ	5/bg	28
VST0814MQ-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQ	20/bx	28
VST0814MQML-20	BP cuff, SoftCheck Vinyl, 1T, infant, 8-14 cm, MQML	20/bx	28
VST1320-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, no fitting	5/bg	28
VST1320HP-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, HP	5/bg	26
VST1320HP-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, HP	20/bx	26
VST1320ML-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, ML	5/bg	27
VST1320ML-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, ML	20/bx	27
VST1320MQ-05	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQ	5/bg	28
VST1320MQ-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQ	20/bx	28
VST1320MQML-20	BP cuff, SoftCheck Vinyl, 1T, child, 13-20 cm, MQML	20/bx	28
VST1826-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, no fitting	5/bg	28
VST1826HP-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	26
VST1826HP-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, HP	20/bx	26
VST1826ML-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, ML	5/bg	27
VST1826ML-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, ML	20/bx	27
VST1826MQ-05	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, MQ	5/bg	28
VST1826MQML-20	BP cuff, SoftCheck Vinyl, 1T, sm adult, 18-26 cm, MQML	20/bx	28
VST2635-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, no fitting	5/bg	28
VST2635HP-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, HP	5/bg	26

Order Code	Description	UOM	Page(s)
VST2635HP-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, HP	20/bx	26
VST2635ML-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, ML	5/bg	27
VST2635ML-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, ML	20/bx	27
VST2635MQ-05	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, MQ	5/bg	28
VST2635MQ-20	BP cuff, SoftCheck Vinyl, 1T, adult, 26-35 cm, MQ	20/bx	28
VST2635MQML-20	BP cuff, SoftCheck Vinyl, 1T, adult, 25-35 cm, MQML	20/bx	28
VST2938HP-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, HP	5/bg	26
VST2938HP-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, HP	20/bx	26
VST2938ML-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, ML	5/bg	27
VST2938ML-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, ML	20/bx	27
VST2938MQ-05	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQ	5/bg	28
VST2938MQ-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQ	20/bx	28
VST2938MQML-20	BP cuff, SoftCheck Vinyl, 1T, adult long, 29-38 cm, MQML	20/bx	28
VST3242-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, no fitting	5/bg	28
VST3242HP-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	26
VST3242HP-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, HP	20/bx	26
VST3242ML-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, ML	5/bg	27
VST3242ML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, ML	20/bx	27
VST3242MQ-05	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQ	5/bg	28
VST3242MQ-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQ	20/bx	28
VST3242MQML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult, 32-42 cm, MQML	20/bx	28
VST3544HP-05	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, HP	5/bg	26
VST3544HP-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, HP	20/bx	26
VST3544ML-05	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, ML	5/bg	27
VST3544ML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, ML	20/bx	27
VST3544MQML-20	BP cuff, SoftCheck Vinyl, 1T, lg adult long, 35-44 cm, MQML	20/bx	28
VST4250-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, no fitting	5/bg	28
VST4250HP-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, HP	5/bg	26
VST4250HP-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, HP	20/bx	26
VST4250ML-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, ML	5/bg	27
VST4250ML-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, ML	20/bx	27
VST4250MQ-05	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQ	5/bg	28
VST4250MQ-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQ	20/bx	28
VST4250MQML-20	BP cuff, SoftCheck Vinyl, 1T, thigh, 42-50 cm, MQML	20/bx	28
YVST1826HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, sm adult, 18-26 cm, HP	5/bg	30
YVST2635HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, adult, 26-35 cm, HP	5/bg	30
YVST3242HP-05	BP cuff, SoftCheck Yellow Vinyl, 1T, lg adult, 32-42 cm, HP	5/bg	30
ZFT-BPF230	BP fitting, BPF230	ea	35
ZFT-BPF240	BP fitting, BP240	ea	35
ZFT-CS901	BP fitting, CS901	ea	36
ZFT-DM-20	BP fitting, DM	20/bx	36
ZFT-DM2YRX1	BP Y adapter, DM2YRX1	ea	35
ZFT-DM532-20	BP fitting, DM532	20/bx	36
ZFT-DMML-20	BP fitting, DMML	20/bx	36
ZFT-FL532-20	BP fitting, FL532	20/bx	36
ZFT-FLL240-20	BP fitting, FLL240	20/bx	36
ZFT-HP-20	BP fitting, HP	20/bx	35
ZFT-HP532-20	BP fitting, HP532	20/bx	35
ZFT-HPDMAD-05	BP Y adapter, HPDMAD	5/bg	35
ZFT-ML2YFL	BP Y adapter, ML2YFL	ea	35
ZFT-ML532-20	BP fitting, ML532	20/bx	36
ZFT-MLL532-20	Male Locking Luer	20/bx	36
ZFT-MQ-20	BP fitting, MQ	20/bx	36
ZFT-MQ2YRX1	BP Y adapter, MQ2YRX1	ea	35
ZFT-MQMF2YFL	Y Adapter, 5/32 Female Luer to 1 MQ & 1 MQML	ea	36

Order Code	Description	UOM	Page(s)
ZFT-MQMF2YHP	BP Y adapter, MQMF2YHP	ea	36
ZFT-MQMF2YRX1	BP Y adapter, MQMF2YRX1	ea	35
ZFT-MQML-20	BP fitting, MQML	20/bx	36
ZFT-MQML2YRX1	BP Y adapter, MQML2YRX1	ea	35
ZFT-MQML5-20	BP fitting, MQML5	20/bx	36
ZFT-MQML5-20	Male Twist-Lock Fitting for Reusable Cuffs	ea	36
ZFT-QC8030	BP fitting, QC8030	ea	35
ZFT-QC8030	Male Metal Bayonet, 5/32" barb	ea	35
ZFT-QC8220	BP fitting, QC8220	ea	35
ZFT-QC8220	Female Metal Bayonet Connector, 1/8" barb	ea	35
ZFT-QC8230	BP fitting, QC8230	ea	35
ZFT-QC8230	Female Metal Bayonet Connector, 5/32" barb	ea	35
ZFT-RX1DM	BP adapter, RX1DM	ea	35
ZFT-RX1DM	5" Tube w/Metal Female Bayonet and Screw Fitting	ea	35
ZFT-RX1FL	BP adapter, RX1FL	ea	35
ZFT-RX1FL	5" Tube w/Metal Female Bayonet and Female Luer	ea	35
ZFT-RX1ML	5 in Tube w/Metal Female Bayonet and Male Luer	ea	35
ZFT-RX1ML	5" Tube w/Metal Female Bayonet and Male Luer	ea	35
ZFT-RX1MQ	BP adapter, RX1MQ	ea	35
ZFT-RX1MQ	5" Tube w/Metal Female Bayonet and Twist-Lock Fitting	ea	35
ZFT-RX1MQML	BP adapter, RX1MQML	ea	35
ZFT-RX1MQML	5" Tube w/Metal Female Bayonet and Male Twist-Lock Fitting	ea	35
ZGA-BPGBHB	BP cuff Inflator w/gauge, BPGBHB	ea	36
ZGA-BPGBHB	BP Gauge w/Black 2 oz. Handbulb	ea	36
ZGA-BPHBRX1	BP cuff Inflator w/gauge, BPHBA1241	ea	36
ZGA-BPHBRX1	Aneroid Gauge BPC Inflator, 8 in. tubing w/Metal Female Bayonet	ea	36
ZHB-4FBPADPRX1	BP manual bulb adapter, 4FBPADPRX1	ea	36
ZHB-4FBPADPRX1	Four foot single tube Manual Bulb BP Cuff Adapter w/Metal Female Bayonet	ea	36
ZHB-BPCOILADPRX1	BP manual bulb adapter, BPCOILADPRX1	ea	35
ZHB-BPCOILADPRX1	Coiled Tube Manual Bulb Adapter w/metal female bayonet fitting	ea	35
ZHB-HBTHM	Hand bulb, THM	ea	36
ZHB-HBTHM	Black 2oz Handbulb w/Metal Thumbwheel	ea	36

SUPPLIES & ACCESSORIES

NON-INVASIVE BLOOD PRESSURE CUFFS



TRULINK® NYLON REUSABLE BLOOD PRESSURE CUFFS



TruLink Nylon Reusable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Color coded for accurate sizing
- Manufactured without PVC/DEHP
- Provides a fluid barrier
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TruLink Soft Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Soft Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Vinyl Disposable Blood Pressure Cuffs are available in seven adult/pediatric sizes for accurate measurements.

- Comfortable soft vinyl
- Easy to clean with PDI cloths
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

715-1122-10	Single Tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single Tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single Tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single Tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single Tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single Tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single Tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



TruLink Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in three adult/pediatric sizes.

- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Yellow Vinyl Disposable BP Cuffs

Single Tube



Bayonet (HP)

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220

715-1135-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1145-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1165-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

TRULINK® NEONATAL NON-PVC BLOOD PRESSURE CUFFS



TruLink Neonatal Non-PVC Blood Pressure Cuffs are available in five neonatal sizes for accurate measurements.

- Color coded for accurate sizing
- Genuine low profile Velcro® closures
- Manufactured without natural latex
- Fabric fold is more comfortable than a weld



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

TruLink Neonatal Non-PVC BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



Mini Male Luer
(ML316)

715-2114-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
715-2124-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
715-2134-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
715-2144-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
715-2154-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box

Single Tube Hose



714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube, 91220	Neonatal	275 cm/9 ft	1 Each

SPACELABS NYLON REUSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Nylon Reusable Blood Pressure Cuffs are available in both regular and long sizes.

- Long size cuffs for children, small adults, adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Nylon Reusable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



**Metal Bayonet
(MHP)**

714-0020-00	Single tube, MHP	Child	12 - 19 cm	1 Each
714-0070-00	Single tube, MHP	Child Long	12 - 19 cm	1 Each
714-0021-00	Single tube, MHP	Small Adult	17 - 25 cm	1 Each
714-0071-00	Single tube, MHP	Small Adult Long	17 - 25 cm	1 Each
714-0023-00	Single tube, MHP	Adult	23 - 33 cm	1 Each
714-0073-00	Single tube, MHP	Adult Long	23 - 33 cm	1 Each
714-0025-00	Single tube, MHP	Large Adult	31 - 40 cm	1 Each
714-0075-00	Single tube, MHP	Large Adult Long	31 - 40 cm	1 Each
714-0027-00	Single tube, MHP	Thigh	38 - 50 cm	1 Each

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS VINYL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs. The Vinyl Disposable Blood Pressure Cuffs are available in seven sizes for accurate measurement.

- Long size cuffs for adults and large adults
- Long cuffs feature additional length allowing ¼ extra turn for greater security and elimination of “pop-offs”



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Vinyl Disposable BP Cuffs

Single Tube

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



Bayonet (HP)

714-1020-00	Single tube, HP	Child	12 - 19 cm	5/Pack
714-1021-00	Single tube, HP	Small Adult	17 - 25 cm	5/Pack
714-1023-00	Single tube, HP	Adult	23 - 33 cm	5/Pack
714-1073-00	Single tube, HP	Adult Long	23 - 33 cm	5/Pack
714-1025-00	Single tube, HP	Large Adult	31 - 40 cm	5/Pack
714-1075-00	Single tube, HP	Large Adult Long	31 - 40 cm	5/Pack
714-1027-00	Single tube, HP	Thigh	38 - 50 cm	5/Pack

Single Tube Hose



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

SPACELABS SOFT CUFF NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



Spacelabs brand cuffs provide a comprehensive solution for your NIBP needs.



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Spacelabs Soft Cuff Neonatal BP Cuffs

Single Tube

For use with Command Modules 90496/91496 and mCare Monitor 91220



Mini Male Luer
(ML316)

714-1028-01	Single tube, 3/16" ML	Size 1	3 - 6 cm	10/Pack
714-1029-01	Single tube, 3/16" ML	Size 2	4 - 8 cm	10/Pack
714-1030-01	Single tube, 3/16" ML	Size 3	6 - 11 cm	10/Pack
714-1032-01	Single tube, 3/16" ML	Size 4	7 - 13 cm	10/Pack
714-1033-01	Single tube, 3/16" ML	Size 5	8 - 15 cm	10/Pack

Single Tube Hose



714-0019-01	Single tube, 90496/91496	Neonatal	275 cm/9 ft	1 Each
714-0062-00	Single tube hose, 91220	Neonatal	275 cm/9 ft	1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Non-Invasive Blood Pressure Hose

Single Tube Hose

For use with Spacelabs Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



714-0018-00	90496/91496/91330	Adult	275 cm/9 ft	1 Each
714-0060-00	90496/91496/91330	Adult	366 cm/12 ft	1 Each
714-0061-00	91220	Adult	275 cm/9 ft	1 Each

Double Tube to Single Tube Hose

For use with Command Modules 90496/91496



714-0018-12	90496/91496	Adult	275 cm/9 ft	1 Each
-------------	-------------	-------	-------------	--------

Hose Connectors

Slip Fit Connector

Connectors fit Spacelabs standard NIBP tubing. Adapts hoses for use with quick connect fittings.



714-0038-00	3mm ID Adult NIBP hoses			1 Each
714-0047-00	4mm ID Adult NIBP hoses			1 Each
714-0046-00	5mm ID Adult NIBP hoses			1 Each
714-0045-00	6mm ID Adult NIBP hoses			1 Each

NIBP accessories for use with Spacelabs and TruLink branded cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

NIBP Cuff Conversion Accessories

Hose Conversion Kit



714-0043-00	Dual hose conversion kit Use with TruLink and Spacelabs cuffs and dual hose blood pressure monitors incorporating thread type fittings	Adult		1 Each
-------------	---	-------	--	--------

Y-Adapters



714-0056-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single-tube 3mm ID hose and monitor/modules	Adult		5/Pack
-------------	--	-------	--	--------








714-0057-00	NIBP dual to single hose Y-adapter Connects dual-tube NIBP cuffs to single tube NIBP hoses and module/monitors	Adult		5/Pack
-------------	---	-------	--	--------

NIBP accessories for use with Spacelabs and TruLink branded neonatal cuffs.

- Hose
- Fittings/Adapters
- Conversion Accessories



Order Code	Description	Size	Range	Qty	
Neonatal BP Cuff Accessories					
Single Tube Hose					
<i>For use with Command Modules 90496/91496 and mCare Monitor 91220</i>					
	714-0019-01	90496/91496	Neonatal	275 cm/9 ft	1 Each
	714-0062-00	91220	Neonatal	275 cm/9 ft	1 Each
Double Tube Cuff to Single Tube Hose					
<i>For use with Command Modules 90496/91496</i>					
	714-0019-12	90496/91496	Neonatal	275 cm/9 ft	1 Each
Metal Cuff Fitting					
<i>For use with NIBP hose 714-0019-01</i>					
	714-0053-10	Hose fitting, metal	Neonatal		10/Pack
Plastic Cuff Fittings					
<i>For use with NIBP Hoses 714-0019-01 & 714-0019-12</i>					
	714-0052-50	Cuff fitting, plastic	Neonatal		50/Pack
	714-0054-50	Cuff fitting, plastic	Infant		50/Pack
Hose Connectors					
	714-0057-10	NIBP non-luer hose metal connector Connects neonatal cuffs to non-Spacelabs monitors with a 3mm ID hose	Neonatal		10/Pack

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Bayonet (HP)

US0814HP	Single tube, HP	Infant	8 - 14 cm	Each
US0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
US1320HP	Single tube, HP	Child	13 - 20 cm	Each
US1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
US1826HP	Single tube, HP	Small Adult	18 - 26 cm	Each
US1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
US2635HP	Single tube, HP	Adult	26 - 35 cm	Each
US2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
US2938HP	Single tube, HP	Adult Long	29 - 38 cm	Each
US2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
US3242HP	Single tube, HP	Large Adult	32 - 42 cm	Each
US3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
US3544HP	Single tube, HP	Large Adult Long	35 - 44 cm	Each
US3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
US4250HP	Single tube, HP	Thigh	42 - 50 cm	Each
US4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag



Metal Bayonet (MHP)

US0814MHP-05	Single tube, MHP	Infant	8 - 14 cm	5/Bag
US1320MHP-05	Single tube, MHP	Child	13 - 20 cm	5/Bag
US1826MHP-05	Single tube, MHP	Small Adult	18 - 26 cm	5/Bag
US2635MHP-05	Single tube, MHP	Adult	26 - 35 cm	5/Bag
US2938MHP-05	Single tube, MHP	Adult Long	29 - 38 cm	5/Bag
US3242MHP-05	Single tube, MHP	Large Adult	32 - 42 cm	5/Bag
US3544MHP-05	Single tube, MHP	Large Adult Long	35 - 44 cm	5/Bag
US4250MHP-05	Single tube, MHP	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Screw (DM)

US0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
US1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
US1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
US2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
US2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
US3544DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
US4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag



Male Luer (ML)

US0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag
US1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag
US1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag
US2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
US2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
US3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
US3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
US4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag



Male Locking Luer (MML)

US0814MML-05	Single tube, MML	Infant	8 - 14 cm	5/Bag
US1320MML-05	Single tube, MML	Child	13 - 20 cm	5/Bag
US1826MML-05	Single tube, MML	Small Adult	18 - 26 cm	5/Bag
US2635MML-05	Single tube, MML	Adult	26 - 35 cm	5/Bag
US2938MML-05	Single tube, MML	Adult Long	29 - 38 cm	5/Bag
US3242MML-05	Single tube, MML	Large Adult	32 - 42 cm	5/Bag
US3544MML-05	Single tube, MML	Large Adult Long	35 - 44 cm	5/Bag
US4250MML-05	Single tube, MML	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



Twist Lock (MQ)

US0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
US1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
US1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
US2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
US2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
US3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
US3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
US4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag



Female Luer (FL)

US0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
US1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
US1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
US2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
US2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
US3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
US3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
US4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag



Male Twist Lock (MQML)

US0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
US1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
US1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
US2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
US2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
US3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
US3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
US4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Single Tube



No Fitting

US0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
US1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
US1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
US2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
US2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
US3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
US3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
US4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag



US0814MB	Single tube, MB w/gauge	Infant	8 - 14 cm	1 Each
US1320MB	Single tube, MB w/gauge	Child	13 - 20 cm	1 Each
US1826MB	Single tube, MB w/gauge	Small Adult	18 - 26 cm	1 Each
US2635MB	Single tube, MB w/gauge	Adult	26 - 35 cm	1 Each
US2938MB	Single tube, MB w/gauge	Adult Long	29 - 38 cm	1 Each
US3242MB	Single tube, MB w/gauge	Large Adult	32 - 42 cm	1 Each
US3544MB	Single tube, MB w/gauge	Large Adult Long	35 - 44 cm	1 Each
US4250MB	Single tube, MB w/gauge	Thigh	42 - 50 cm	1 Each

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Double Tube



Screw (DM)

UD0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
UD1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
UD1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
UD2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
UD2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
UD3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
UD3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
UD4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag



Bayonet (HP)

UD0814HP-05	Double tube, HP	Infant	8 - 14 cm	5/Bag
UD1320HP-05	Double tube, HP	Child	13 - 20 cm	5/Bag
UD1826HP-05	Double tube, HP	Small Adult	18 - 26 cm	5/Bag
UD2635HP-05	Double tube, HP	Adult	26 - 35 cm	5/Bag
UD2938HP-05	Double tube, HP	Adult Long	29 - 38 cm	5/Bag
UD3242HP-05	Double tube, HP	Large Adult	32 - 42 cm	5/Bag
UD3544HP-05	Double tube, HP	Large Adult Long	35 - 44 cm	5/Bag



Twist Lock (MQ)

UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag
UD0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag
UD1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag
UD1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag
UD2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag
UD2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag
UD3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag
UD3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® NYLON REUSABLE BLOOD PRESSURE CUFFS



UltraCheck Nylon Reusable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Fluid barrier
- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Nylon Reusable BP Cuffs

Double Tube



Male & Female MQ (MQMF)

UD0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
UD1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
UD1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
UD2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
UD2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
UD3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
UD3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
UD4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag



No Fitting

UD0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
UD1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag
UD1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
UD2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
UD2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
UD3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
UD3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
UD4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag



Manual Bulb (MB)

UD0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
UD1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag
UD1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
UD2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
UD2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
UD3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
UD3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
UD4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

ULTRACHECK® CURVE™ REUSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Manufactured without DEHP
- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

UltraCheck Curve Reusable BP Cuffs

Single Tube



BRUS3854	Single tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854DM	Single tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854FL	Single tube, FL	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854HP	Single tube, HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MHP	Single tube, metal HP	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854ML	Single tube, ML	Extra Large Adult	38 - 54 cm	1 Each
BRUS3854MQ	Single tube, MQ	Extra Large Adult	38 - 54 cm	1 Each

Double Tube



BRUD3854	Double tube, no fitting	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854DM	Double tube, DM	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854HP	Double tube, HP	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854ML	Double tube, ML	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQ	Double tube, MQ	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQMF	Double tube, MQMF	Extra Large Adult	38 - 54 cm	1 Each
BRUD3854MQML	Double tube, MQML	Extra Large Adult	38 - 54 cm	1 Each

Manual Cuffs



BRUS3854MB	Single tube, manual w/ gauge	Extra Large Adult	38 - 54 cm	1 Each
------------	------------------------------	-------------------	------------	--------



BRUD3854MB	Double tube, manual	Extra Large Adult	38 - 54 cm	1 Each
------------	---------------------	-------------------	------------	--------

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Bayonet (HP)

ST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag
ST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box
ST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag
ST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box
ST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
ST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box
ST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
ST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box
ST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
ST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box
ST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
ST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box
ST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
ST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box
ST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag
ST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box



Male Luer (ML)

ST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag
ST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	5/Bag
ST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box
ST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag
ST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box
ST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag
ST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box
ST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag
ST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

ST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
ST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
ST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
ST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
ST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
ST1826MQ-20	Single tube, MQ	Small Adult	18 - 26 cm	20/Box
ST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
ST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
ST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	5/Bag
ST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
ST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
ST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
ST3544MQ-05	Single tube, MQ	Large Adult Long	35 - 44 cm	5/Bag
ST3544MQ-20	Single tube, MQ	Large Adult Long	35 - 44 cm	20/Box
ST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
ST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Screw (DM)

ST0814DM-05	Single tube, DM	Infant	8 - 14 cm	5/Bag
ST1320DM-05	Single tube, DM	Child	13 - 20 cm	5/Bag
ST1826DM-05	Single tube, DM	Small Adult	18 - 26 cm	5/Bag
ST2635DM-05	Single tube, DM	Adult	26 - 35 cm	5/Bag
ST2938DM-05	Single tube, DM	Adult Long	29 - 38 cm	5/Bag
ST3242DM-05	Single tube, DM	Large Adult	32 - 42 cm	5/Bag
ST3544DM-05	Single tube, DM	Large Adult Long	35 - 44 cm	5/Bag
ST4250DM-05	Single tube, DM	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Single Tube



Female Luer (FL)

ST0814FL-05	Single tube, FL	Infant	8 - 14 cm	5/Bag
ST1320FL-05	Single tube, FL	Child	13 - 20 cm	5/Bag
ST1826FL-05	Single tube, FL	Small Adult	18 - 26 cm	5/Bag
ST2635FL-05	Single tube, FL	Adult	26 - 35 cm	5/Bag
ST2938FL-05	Single tube, FL	Adult Long	29 - 38 cm	5/Bag
ST3242FL-05	Single tube, FL	Large Adult	32 - 42 cm	5/Bag
ST3544FL-05	Single tube, FL	Large Adult Long	35 - 44 cm	5/Bag
ST4250FL-05	Single tube, FL	Thigh	42 - 50 cm	5/Bag



Male Twist Lock (MQML)

ST0814MQML-05	Single tube, MQML	Infant	8 - 14 cm	5/Bag
ST1320MQML-05	Single tube, MQML	Child	13 - 20 cm	5/Bag
ST1826MQML-05	Single tube, MQML	Small Adult	18 - 26 cm	5/Bag
ST2635MQML-05	Single tube, MQML	Adult	26 - 35 cm	5/Bag
ST2938MQML-05	Single tube, MQML	Adult Long	29 - 38 cm	5/Bag
ST3242MQML-05	Single tube, MQML	Large Adult	32 - 42 cm	5/Bag
ST3544MQML-05	Single tube, MQML	Large Adult Long	35 - 44 cm	5/Bag
ST4250MQML-05	Single tube, MQML	Thigh	42 - 50 cm	5/Bag



No Fitting

ST0814-05	Single tube, no fitting	Infant	8 - 14 cm	5/Bag
ST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
ST1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
ST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
ST2938-05	Single tube, no fitting	Adult Long	29 - 38 cm	5/Bag
ST3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
ST3544-05	Single tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
ST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs



Screw (DM)

Double Tube				
DT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
DT0814DM-20	Double tube, DM	Infant	8 - 14 cm	20/Box
DT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
DT1320DM-20	Double tube, DM	Child	13 - 20 cm	20/Box
DT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
DT1826DM-20	Double tube, DM	Small Adult	18 - 26 cm	20/Box
DT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
DT2635DM-20	Double tube, DM	Adult	26 - 35 cm	20/Box
DT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
DT2938DM-20	Double tube, DM	Adult Long	29 - 38 cm	20/Box
DT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
DT3242DM-20	Double tube, DM	Large Adult	32 - 42 cm	20/Box
DT3544DM-05	Double tube, DM	Large Adult Long	35 - 44 cm	5/Bag
DT3544DM-20	Double tube, DM	Large Adult Long	35 - 44 cm	20/Box
DT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
DT4250DM-20	Double tube, DM	Thigh	42 - 50 cm	20/Box

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs



Twist Lock (MQ)

Double Tube					
DT0814MQ-05	Double tube, MQ	Infant	8 - 14 cm	5/Bag	
DT0814MQ-20	Double tube, MQ	Infant	8 - 14 cm	20/Box	
DT1320MQ-05	Double tube, MQ	Child	13 - 20 cm	5/Bag	
DT1320MQ-20	Double tube, MQ	Child	13 - 20 cm	20/Box	
DT1826MQ-05	Double tube, MQ	Small Adult	18 - 26 cm	5/Bag	
DT1826MQ-20	Double tube, MQ	Small Adult	18 - 26 cm	20/Box	
DT2635MQ-05	Double tube, MQ	Adult	26 - 35 cm	5/Bag	
DT2635MQ-20	Double tube, MQ	Adult	26 - 35 cm	20/Box	
DT2938MQ-05	Double tube, MQ	Adult Long	29 - 38 cm	5/Bag	
DT2938MQ-20	Double tube, MQ	Adult Long	29 - 38 cm	20/Box	
DT3242MQ-05	Double tube, MQ	Large Adult	32 - 42 cm	5/Bag	
DT3242MQ-20	Double tube, MQ	Large Adult	32 - 42 cm	20/Box	
DT3544MQ-05	Double tube, MQ	Large Adult Long	35 - 44 cm	5/Bag	
DT3544MQ-20	Double tube, MQ	Large Adult Long	35 - 44 cm	20/Box	
DT4250MQ-05	Double tube, MQ	Thigh	42 - 50 cm	5/Bag	
DT4250MQ-20	Double tube, MQ	Thigh	42 - 50 cm	20/Box	

SOFTCHECK® FABRIC DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Fabric Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Soft fabric helps prevent perspiration under cuff
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Disposable BP Cuffs

Double Tube



DT0814MQMF-05	Double tube, MQMF	Infant	8 - 14 cm	5/Bag
DT1320MQMF-05	Double tube, MQMF	Child	13 - 20 cm	5/Bag
DT1826MQMF-05	Double tube, MQMF	Small Adult	18 - 26 cm	5/Bag
DT2635MQMF-05	Double tube, MQMF	Adult	26 - 35 cm	5/Bag
DT2938MQMF-05	Double tube, MQMF	Adult Long	29 - 38 cm	5/Bag
DT3242MQMF-05	Double tube, MQMF	Large Adult	32 - 42 cm	5/Bag
DT3544MQMF-05	Double tube, MQMF	Large Adult Long	35 - 44 cm	5/Bag
DT4250MQMF-05	Double tube, MQMF	Thigh	42 - 50 cm	5/Bag



DT0814-05	Double tube, no fitting	Infant	8 - 14 cm	5/Bag
DT1320-05	Double tube, no fitting	Child	13 - 20 cm	5/Bag
DT1826-05	Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
DT2635-05	Double tube, no fitting	Adult	26 - 35 cm	5/Bag
DT2938-05	Double tube, no fitting	Adult Long	29 - 38 cm	5/Bag
DT3242-05	Double tube, no fitting	Large Adult	32 - 42 cm	5/Bag
DT3544-05	Double tube, no fitting	Large Adult Long	35 - 44 cm	5/Bag
DT4250-05	Double tube, no fitting	Thigh	42 - 50 cm	5/Bag



DT0814MB-05	Double tube, MB	Infant	8 - 14 cm	5/Bag
DT1320MB-05	Double tube, MB	Child	13 - 20 cm	5/Bag
DT1826MB-05	Double tube, MB	Small Adult	18 - 26 cm	5/Bag
DT2635MB-05	Double tube, MB	Adult	26 - 35 cm	5/Bag
DT2938MB-05	Double tube, MB	Adult Long	29 - 38 cm	5/Bag
DT3242MB-05	Double tube, MB	Large Adult	32 - 42 cm	5/Bag
DT3544MB-05	Double tube, MB	Large Adult Long	35 - 44 cm	5/Bag
DT4250MB-05	Double tube, MB	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® CURVE™ DISPOSABLE BLOOD PRESSURE CUFFS



Worldwide obesity rates are climbing. This disease is strongly associated with increased mortality and morbidity, cardiovascular disease, diabetes, and many other illnesses. Blood pressure is the number one vital sign in detecting these issues and accurate blood pressure is critical. The Curve Cuffs are designed to fit a bariatric patient's arm correctly. Their conical shape fits snugly against the arm and properly compresses the artery.

- Will not slide during inflation
- Genuine Velcro® closures
- Manufactured without natural latex
- Adapters allow for compatibility with monitors worldwide



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Curve Disposable BP Cuffs

Single Tube

BRST3854HP-05	Single tube, HP	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854DM-05	Single tube, DM	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854FL-05	Single tube, FL	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854ML-05	Single tube, ML	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854MQ-05	Single tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854MQML-05	Single tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag
BRST3854-05	Single tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag



Double Tube

BRDT3854HP-05	Double tube, HP	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854DM-05	Double tube, DM	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854FL-05	Double tube, FL	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854ML-05	Double tube, ML	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQ-05	Double tube, MQ	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQMF-05	Double tube, MQMF	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854MQML-05	Double tube, MQML	Extra Large Adult	38 - 54 cm	5/Bag
BRDT3854-05	Double tube, no fitting	Extra Large Adult	38 - 54 cm	5/Bag



Manual Cuff

BRDT3854MB-05	Double tube, MB	Extra Large Adult	38 - 54 cm	5/Bag
---------------	-----------------	-------------------	------------	-------



SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs



Bayonet (HP)

Single Tube					
VST0814HP-05	Single tube, HP	Infant	8 - 14 cm	5/Bag	
VST0814HP-20	Single tube, HP	Infant	8 - 14 cm	20/Box	
VST1320HP-05	Single tube, HP	Child	13 - 20 cm	5/Bag	
VST1320HP-20	Single tube, HP	Child	13 - 20 cm	20/Box	
VST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag	
VST1826HP-20	Single tube, HP	Small Adult	18 - 26 cm	20/Box	
VST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag	
VST2635HP-20	Single tube, HP	Adult	26 - 35 cm	20/Box	
VST2938HP-05	Single tube, HP	Adult Long	29 - 38 cm	5/Bag	
VST2938HP-20	Single tube, HP	Adult Long	29 - 38 cm	20/Box	
VST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag	
VST3242HP-20	Single tube, HP	Large Adult	32 - 42 cm	20/Box	
VST3544HP-05	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag	
VST3544HP-20	Single tube, HP	Large Adult Long	35 - 44 cm	20/Box	
VST4250HP-05	Single tube, HP	Thigh	42 - 50 cm	5/Bag	
VST4250HP-20	Single tube, HP	Thigh	42 - 50 cm	20/Box	

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube					
VST0814ML-05	Single tube, ML	Infant	8 - 14 cm	5/Bag	
VST0814ML-20	Single tube, ML	Infant	8 - 14 cm	20/Box	
VST1320ML-05	Single tube, ML	Child	13 - 20 cm	5/Bag	
VST1320ML-20	Single tube, ML	Child	13 - 20 cm	20/Box	
VST1826ML-05	Single tube, ML	Small Adult	18 - 26 cm	5/Bag	
VST1826ML-20	Single tube, ML	Small Adult	18 - 26 cm	20/Box	
VST2635ML-05	Single tube, ML	Adult	26 - 35 cm	5/Bag	
VST2635ML-20	Single tube, ML	Adult	26 - 35 cm	20/Box	
VST2938ML-05	Single tube, ML	Adult Long	29 - 38 cm	20/Box	
VST2938ML-20	Single tube, ML	Adult Long	29 - 38 cm	20/Box	
VST3242ML-05	Single tube, ML	Large Adult	32 - 42 cm	5/Bag	
VST3242ML-20	Single tube, ML	Large Adult	32 - 42 cm	20/Box	
VST3544ML-05	Single tube, ML	Large Adult Long	35 - 44 cm	5/Bag	
VST3544ML-20	Single tube, ML	Large Adult Long	35 - 44 cm	20/Box	
VST4250ML-05	Single tube, ML	Thigh	42 - 50 cm	5/Bag	
VST4250ML-20	Single tube, ML	Thigh	42 - 50 cm	20/Box	



Male Luer (ML)

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Single Tube



Twist Lock (MQ)

VST0814MQ-05	Single tube, MQ	Infant	8 - 14 cm	5/Bag
VST0814MQ-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQ-05	Single tube, MQ	Child	13 - 20 cm	5/Bag
VST1320MQ-20	Single tube, MQ	Child	13 - 20 cm	20/Box
VST1826MQ-05	Single tube, MQ	Small Adult	18 - 26 cm	5/Bag
VST2635MQ-05	Single tube, MQ	Adult	26 - 35 cm	5/Bag
VST2635MQ-20	Single tube, MQ	Adult	26 - 35 cm	20/Box
VST2938MQ-05	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST2938MQ-20	Single tube, MQ	Adult Long	29 - 38 cm	20/Box
VST3242MQ-05	Single tube, MQ	Large Adult	32 - 42 cm	5/Bag
VST3242MQ-20	Single tube, MQ	Large Adult	32 - 42 cm	20/Box
VST4250MQ-05	Single tube, MQ	Thigh	42 - 50 cm	5/Bag
VST4250MQ-20	Single tube, MQ	Thigh	42 - 50 cm	20/Box



Male Twist Lock (MQML)

VST0814MQML-20	Single tube, MQ	Infant	8 - 14 cm	20/Box
VST1320MQML-20	Single tube, MQML	Child	13 - 20 cm	20/Box
VST1826MQML-20	Single tube, MQML	Small Adult	18 - 26 cm	20/Box
VST2635MQML-20	Single tube, MQML	Adult	26 - 35 cm	20/Box
VST2938MQML-20	Single tube, MQML	Adult Long	29 - 38 cm	20/Box
VST3242MQML-20	Single tube, MQML	Large Adult	32 - 42 cm	20/Box
VST3544MQML-20	Single tube, MQML	Thigh	35 - 44 cm	20/Box
VST4250MQML-20	Single tube, MQML	Thigh	42 - 50 cm	20/Box



No Fitting

VST1320-05	Single tube, no fitting	Child	13 - 20 cm	5/Bag
VST1826-05	Single tube, no fitting	Small Adult	18 - 26 cm	5/Bag
VST2635-05	Single tube, no fitting	Adult	26 - 35 cm	5/Bag
VST3242-05	Single tube, no fitting	Large Adult	32 - 42 cm	5/Bag
VST4250-05	Single tube, no fitting	Thigh	42 - 50 cm	5/Bag

SOFTCHECK® VINYL DISPOSABLE BLOOD PRESSURE CUFFS





SoftCheck Vinyl Disposable Blood Pressure Cuffs are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Disposable BP Cuffs

Double Tube					
 Screw (DM)	VDT0814DM-05	Double tube, DM	Infant	8 - 14 cm	5/Bag
	VDT1320DM-05	Double tube, DM	Child	13 - 20 cm	5/Bag
	VDT1826DM-05	Double tube, DM	Small Adult	18 - 26 cm	5/Bag
	VDT2635DM-05	Double tube, DM	Adult	26 - 35 cm	5/Bag
	VDT2938DM-05	Double tube, DM	Adult Long	29 - 38 cm	5/Bag
	VDT3242DM-05	Double tube, DM	Large Adult	32 - 42 cm	5/Bag
	VDT4250DM-05	Double tube, DM	Thigh	42 - 50 cm	5/Bag
	 No Fitting	VDT0814-05	Double tube, no fitting	Infant	8 - 14 cm
VDT1320-05		Double tube, no fitting	Child	13 - 20 cm	20/Box
VDT1826-05		Double tube, no fitting	Small Adult	18 - 26 cm	5/Bag
VDT2635-05		Double tube, no fitting	Adult	26 - 35 cm	20/Box
VDT2938-05		Double tube, no fitting	Adult Long	29 - 38 cm	20/Box
VDT3242-05		Double tube, no fitting	Large Adult	32 - 42 cm	20/Box
VDT4250-05		Double tube, no fitting	Thigh	42 - 50 cm	20/Box

SOFTCHECK® YELLOW VINYL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Yellow Vinyl Disposable Blood Pressure Cuffs intended for use in isolation areas are available in multiple sizes along with a variety of adapters allowing you to provide standardization across your facility.

- Comfortable soft vinyl
- Color coded for accurate sizing
- Genuine Velcro® closures
- Manufactured without natural latex
- Fabric fold at elbow is more comfortable than a weld
- Standardization - compatible with most brands of automated blood pressure monitors and wall or portable manometers



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Yellow Vinyl Disposable BP Cuffs



Bayonet (HP)

Single Tube				
YVST1826HP-05	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
YVST2635HP-05	Single tube, HP	Adult	26 - 35 cm	5/Bag
YVST3242HP-05	Single tube, HP	Large Adult	32 - 42 cm	5/Bag

SOFTCHECK® FABRIC NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.

- Meets AHA / AAMI standards
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Fabric Neonatal Disposable BP Cuffs

Single Tube



Mini Male Luer
(ML316)

NN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
NN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
NN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
NN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
NN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box



Male Luer (ML)

NN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
NN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
NN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
NN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
NN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Bayonet (HP)

NN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
NN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
NN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
NN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
NN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



Male Luer (ML)

NN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
NN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
NN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
NN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
NN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.

- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs

Single Tube



Mini Male Luer (ML316)

VNN1ST-M316-10	Single tube, 3/16" ML	Size 1	3-6 cm	10/Box
VNN2ST-M316-10	Single tube, 3/16" ML	Size 2	4-8 cm	10/Box
VNN3ST-M316-10	Single tube, 3/16" ML	Size 3	6-11 cm	10/Box
VNN4ST-M316-10	Single tube, 3/16" ML	Size 4	7-13 cm	10/Box
VNN5ST-M316-10	Single tube, 3/16" ML	Size 5	8-15 cm	10/Box



SnapQuik (SQ)

VNN1ST-SQ-10	Single tube, SQ	Size 1	3-6 cm	10/Box
VNN2ST-SQ-10	Single tube, SQ	Size 2	4-8 cm	10/Box
VNN3ST-SQ-10	Single tube, SQ	Size 3	6-11 cm	10/Box
VNN4ST-SQ-10	Single tube, SQ	Size 4	7-13 cm	10/Box
VNN5ST-SQ-10	Single tube, SQ	Size 5	8-15 cm	10/Box



Male Luer (ML)

VNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
VNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
VNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
VNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
VNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Bayonet (HP)

VNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
VNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
VNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
VNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
VNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box



Twist Lock (MQ)

VNN1ST-MQ-10	Single tube, MQ	Size 1	3-6 cm	10/Box
VNN2ST-MQ-10	Single tube, MQ	Size 2	4-8 cm	10/Box
VNN3ST-MQ-10	Single tube, MQ	Size 3	6-11 cm	10/Box
VNN4ST-MQ-10	Single tube, MQ	Size 4	7-13 cm	10/Box
VNN5ST-MQ-10	Single tube, MQ	Size 5	8-15 cm	10/Box

SOFTCHECK® VINYL NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.

- Comfortable soft vinyl, created with babies' skin integrity in mind
- Easy to clean
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Vinyl Neonatal Disposable BP Cuffs



Male Luer (ML)

Double Tube					
VNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box	
VNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box	
VNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box	
VNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box	
VNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box	

SOFTCHECK FREE® NEONATAL DISPOSABLE BLOOD PRESSURE CUFFS



SoftCheck Free Neonatal Disposable Blood Pressure Cuffs are some of the world's finest quality blood pressure cuffs. In a NICU environment, this means better care for babies as well as greater effectiveness for practitioners.

- Created with extremely soft material
- Color coded for accurate sizing
- 5 sizes for accurate measurements
- Genuine low profile Velcro® closures



Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

SoftCheck Free Non-PVC Neonatal Disposable BP Cuffs

Single tube



Mini Male Luer (ML316)

PNN1ST-M316-10	Single tube, M316	Size 1	3-6 cm	10/Box
PNN2ST-M316-10	Single tube, M316	Size 2	4-8 cm	10/Box
PNN3ST-M316-10	Single tube, M316	Size 3	6-11 cm	10/Box
PNN4ST-M316-10	Single tube, M316	Size 4	7-13 cm	10/Box
PNN5ST-M316-10	Single tube, M316	Size 5	8-15 cm	10/Box



Male Luer (ML)

PNN1ST-10	Single tube, ML	Size 1	3-6 cm	10/Box
PNN2ST-10	Single tube, ML	Size 2	4-8 cm	10/Box
PNN3ST-10	Single tube, ML	Size 3	6-11 cm	10/Box
PNN4ST-10	Single tube, ML	Size 4	7-13 cm	10/Box
PNN5ST-10	Single tube, ML	Size 5	8-15 cm	10/Box



Bayonet (HP)

PNN1ST-HP-10	Single tube, HP	Size 1	3-6 cm	10/Box
PNN2ST-HP-10	Single tube, HP	Size 2	4-8 cm	10/Box
PNN3ST-HP-10	Single tube, HP	Size 3	6-11 cm	10/Box
PNN4ST-HP-10	Single tube, HP	Size 4	7-13 cm	10/Box
PNN5ST-HP-10	Single tube, HP	Size 5	8-15 cm	10/Box

Double Tube



Male Luer (ML)

PNN1DT-10	Double tube, ML	Size 1	3-6 cm	10/Box
PNN2DT-10	Double tube, ML	Size 2	4-8 cm	10/Box
PNN3DT-10	Double tube, ML	Size 3	6-11 cm	10/Box
PNN4DT-10	Double tube, ML	Size 4	7-13 cm	10/Box
PNN5DT-10	Double tube, ML	Size 5	8-15 cm	10/Box

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
------------	-------------	-----

Blood Pressure Cuff Fittings & Adapters

Used in conversions to HP Cuffs

ZFT-HP-20	Male Bayonet Connector for Disposable Cuffs	20/Pack
ZFT-HP532-20	Male Bayonet Connector for Reusable Cuffs	20/Pack
ZFT-QC8030	Male Metal Bayonet Connector, 5/32" barb	1 Each
ZFT-QC8220	Female Metal Bayonet Connector, 1/8" barb	1 Each
ZFT-QC8230	Female Metal Bayonet Connector, 5/32" barb	1 Each
ZFT-BPF230	Female Plastic Bayonet Connector, 1/8" barb	1 Each
ZFT-BPF240	Female Plastic Bayonet Connector, 5/32" barb	1 Each
ZHB-BPCOILADPRX1	Coiled Tube Manual Bulb Adapter w/metal female bayonet fitting	1 Each
ADP-182YBPF	Y Adapter w/1/8" barb fittings & plastic female bayonet connector	1 Each
ADP-182YRX4	Y Adapter w/1/8" barb fittings & metal female bayonet connector	1 Each
ADP-5322YBPF	Y Adapter w/5/32" barb fittings & plastic female bayonet connector	1 Each
ADP-5322YRX4	Y Adapter w/5/32" barb fittings & metal female bayonet connector	1 Each
ADP-DM2YRX4	Y Adapter, Metal Female HP to 2 DM	1 Each
ZFT-HPDMAD-05	Y Adapter /Male Bayonet to 2 Male Screw	5/Pack
ZFT-ML2YFL	Y Adapter w/Female Luer to 2 Male Luers	1 Each

Used in conversions to HP Cuffs-Straight

ZFT-RX1DM	5" Tube w/Metal Female Bayonet and Screw Fitting	1 Each
ZFT-RX1FL	5" Tube w/Metal Female Bayonet and Female Luer	1 Each
ZFT-RX1ML	5" Tube w/Metal Female Bayonet and Male Luer	1 Each
ZFT-RX1MQ	5" Tube w/Metal Female Bayonet and Twist-Lock Fitting	1 Each
ZFT-RX1MQML	5" Tube w/Metal Female Bayonet and Male Twist-Lock Fitting	1 Each

Used in conversions to HP Cuffs-Ys

ZFT-DM2YRX1	Y Adapter, Metal Female Bayonet to 2 Screw fittings	1 Each
ZFT-MQ2YRX1	Y Adapter w/Metal Female Bayonet to 2 Twist Lock fittings	1 Each
ZFT-MQML2YRX1	Y Adapter w/Metal Female Bayonet to 2 Male Twist Lock fittings	1 Each
ZFT-MQMF2YRX1	Y Adapter w/Metal Female Bayonet to one female Twist Lock and one male Twist Lock fitting	1 Each

Eighty-five percent of all blood pressure cuffs are made with the bayonet (HP) fitting. Spacelabs Healthcare can convert your monitors and sphygmomanometers to this fitting to save you time, energy and costs. Plus, your patient can move from Emergency, to the OR, to PACU, to the general floor with one cuff that attaches to every monitor in the hospital.



Order Code	Description	Qty
------------	-------------	-----

Blood Pressure Cuff Fittings & Adapters

Miscellaneous		
ZFT-CS901	Brass Male Screw Fitting	1 Each
ZFT-DM-20	Female Screw Connector for Disposable Cuffs	20/Pack
ZFT-DM532-20	Female Screw Connector for Reusable Cuffs	20/Pack
ZFT-DMML-20	Nylon Male Screw Fitting	20/Pack
ZFT-FL532-20	Female Luer for Disposable Cuffs	20/Pack
ZFT-FLL240-20	Female Luer for Reusable Cuffs	20/Pack
ZFT-ML532-20	Male Slip Luer for Disposable and Reusable Cuffs	20/Pack
ZFT-MLL532-20	Male Locking Luer	20/Pack
ZFT-MQ-20	Female Twist-Lock Fitting for Disposable and Reusable Cuffs	20/Pack
ZFT-MQMF2YHP	Y Adapter - Twist Lock and Male Twist Lock to Male Bayonet	1 Each
ZFT-MQML-20	Male Twist-Lock Fitting for Disposable Cuffs	20/Pack
ZFT-MQML5-20	Male Twist-Lock Fitting for Reusable Cuffs	20/Pack
ZGA-BPGBHB	BP Gauge w/Black 2 oz. Handbulb	1 Each
ZGA-BPHBRX1	Aneroid Gauge BPC Inflator, 8 in. tubing w/Metal Female Bayonet	1 Each
ZHB-4FBPADPRX1	4' single tube Manual Bulb BP Cuff Adapter w/Metal Female Bayonet	1 Each
ZHB-HBTHM	Black 2oz Handbulb w/Metal Thumbwheel	1 Each

SUPPLIES & ACCESSORIES

PRESSURE INFUSORS



Unifusor is designed for rapid IV infusion and invasive pressure monitoring. The Unifusor line is the clear choice in pressure infusors.

3 Models: Unifusor, Unifusor II, Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook



Order Code	Description	Range	Qty
------------	-------------	-------	-----

Unifusor Disposable Pressure Infusor



Disposable Infusor with aneroid gauge - sold as one piece/complete unit

803FGA-24	Unifusor with thumbwheel valve & aneroid gauge	500mL	24/Box
804FGA-24	Unifusor with thumbwheel valve & aneroid gauge	1000mL	24/Box
805FGA-10	Unifusor with thumbwheel valve & aneroid gauge	3000mL	10/Box



803SGA-24	Unifusor with stopcock valve & aneroid gauge	500mL	24/Box
804SGA-24	Unifusor with stopcock valve & aneroid gauge	1000mL	24/Box
805SGA-10	Unifusor with stopcock valve & aneroid gauge	3000mL	10/Box



Disposable Infusor with piston gauge - sold as one piece/complete unit

903TGA-24	Unifusor with thumbwheel valve & piston gauge	500mL	24/Box
904TGA-24	Unifusor with thumbwheel valve & piston gauge	1000mL	24/Box
905TGA-10	Unifusor with thumbwheel valve & piston gauge	3000mL	10/Box



903SGA-24	Unifusor with stopcock valve & piston gauge	500mL	24/Box
904SGA-24	Unifusor with stopcock valve & piston gauge	1000mL	24/Box
905SGA-10	Unifusor with stopcock valve & piston gauge	3000mL	10/Box

Unifusor is designed for rapid IV infusion and invasive pressure monitoring. The Unifusor line is the clear choice in pressure infusors.

3 Models: Unifusor, Unifusor II, Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook



Order Code	Description	Range	Qty
------------	-------------	-------	-----

Unifusor II Disposable Pressure Infusor with Reusable Inflation System



Disposable Infusor - white disposable bag with reusable inflation system

2803FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	500mL	24/Box
2804FGA-24	Unifusor II with aneroid gauge & thumbwheel valve	1000mL	24/Box
2805FGA-10	Unifusor II with aneroid gauge & thumbwheel valve	3000mL	10/Box



2803SGA-24	Unifusor II with aneroid gauge & stopcock valve	500mL	24/Box
2804SGA-24	Unifusor II with aneroid gauge & stopcock valve	1000mL	24/Box
2805SGA-10	Unifusor II with aneroid gauge & stopcock valve	3000mL	10/Box



2903TGA-24	Unifusor II with piston gauge & thumbwheel valve	500mL	24/Box
2904TGA-24	Unifusor II with piston gauge & thumbwheel valve	1000mL	24/Box
2905TGA-10	Unifusor II with piston gauge & thumbwheel valve	3000mL	10/Box



2903SGA-24	Unifusor II with piston gauge & stopcock valve	500mL	24/Box
2904SGA-24	Unifusor II with piston gauge & stopcock valve	1000mL	24/Box
2905SGA-10	Unifusor II with piston gauge & stopcock valve	3000mL	10/Box

Disposable Replacement Parts




2803X-24	Unifusor II disposable white infusion bag	500mL	24/Box
2804X-24	Unifusor II disposable white infusion bag	1000mL	24/Box
2805X-10	Unifusor II disposable white infusion bag	500mL	10/Box

Reusable Replacement Parts



899F-04	Aneroid gauge & thumbwheel valve assembly	4/Pk
899S-04	Aneroid gauge & stopcock valve assembly	4/Pk



999T-04	Piston gauge & thumbwheel valve assembly	4/Pk
999S-04	Piston gauge & stopcock valve assembly	4/Pk

Unifusor is designed for rapid IV infusion and invasive pressure monitoring. The Unifusor line is the clear choice in pressure infusors.

3 Models: Unifusor, Unifusor II, Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook



Order Code	Description	Range	Qty
------------	-------------	-------	-----

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with white hose

1103FGA-24	Unifusor Plus with aneroid gauge & thumbwheel valve	500mL	24/Box
1104FGA-24	Unifusor Plus with aneroid gauge & thumbwheel valve	1000mL	24/Box
1105FGA-10	Unifusor Plus with aneroid gauge & thumbwheel valve	3000mL	10/Box



1103SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	500mL	24/Box
1104SGA-24	Unifusor Plus with aneroid gauge & stopcock valve	1000mL	24/Box
1105SGA-10	Unifusor Plus with aneroid gauge & stopcock valve	3000mL	10/Box



1103TGA-24	Unifusor Plus with piston gauge & thumbwheel valve	500mL	24/Box
1104TGA-24	Unifusor Plus with piston gauge & thumbwheel valve	1000mL	24/Box
1105TGA-10	Unifusor Plus with piston gauge & thumbwheel valve	3000mL	10/Box



1113SGA-24	Unifusor Plus with piston gauge & stopcock valve	500mL	24/Box
1114SGA-24	Unifusor Plus with piston gauge & stopcock valve	1000mL	24/Box
1115SGA-10	Unifusor Plus with piston gauge & stopcock valve	3000mL	10/Box


Reusable replacement parts



1103XB-05	Unifusor Plus blue reusable bag with white hose	500mL	5/Box
1104XB-05	Unifusor Plus blue reusable bag with white hose	1000mL	5/Box
1105XB-05	Unifusor Plus blue reusable bag with white hose	3000mL	5/Box



899F-04	Aneroid gauge & thumbwheel valve assembly		4/Pk
899S-04	Aneroid gauge & stopcock valve assembly		4/Pk



999T-04	Piston gauge & thumbwheel valve assembly		4/Pk
999S-04	Piston gauge & stopcock valve assembly		4/Pk

Unifusor is designed for rapid IV infusion and invasive pressure monitoring. The Unifusor line is the clear choice in pressure infusors.

3 Models: Unifusor, Unifusor II, Unifusor Plus

- Manufactured without natural latex
- 100% leak tested during manufacturing process
- Reinforced translucent sleeve (see type of fluid and level through sleeve)
- IV pole loop
- Fluid bag hook



Order Code	Description	Range	Qty
------------	-------------	-------	-----

Unifusor Plus Reusable Pressure Infusor



Reusable Infusor - blue bag with black hose

1103	Unifusor Plus with trigger / gauge assembly	500mL	1 Each
1104	Unifusor Plus with trigger / gauge assembly	1000mL	1 Each
1105	Unifusor Plus with trigger / gauge assembly	3000mL	1 Each

Reusable replacement parts



1103X	Unifusor Plus blue reusable bag with black hose	500mL	1 Each
1103X-05	Unifusor Plus blue reusable bag with black hose	500mL	5/Box
1104X	Unifusor Plus blue reusable bag with black hose	1000mL	1 Each
1104X-05	Unifusor Plus blue reusable bag with black hose	1000mL	5/Box
1105X	Unifusor Plus blue reusable bag with black hose	3000mL	1 Each
1105X-05	Unifusor Plus blue reusable bag with black hose	3000mL	5/Box



1199	Trigger/ gauge assembly only		1 Each
------	------------------------------	--	--------

SUPPLIES & ACCESSORIES
ELECTROCARDIOGRAPHY
AAMI Color Coding



ECG COMBINER LEAD WIRE SYSTEM

AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter

Disposable, DIN Shielded, Multi-Pinch, AAMI



700-0006-30	3 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-31	4 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets
700-0006-32	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-33	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets
700-0006-34	6 Lead	Multi-Pinch	127cm /50 in	25 Sets

Disposable, DIN Shielded, Snap, AAMI



700-0007-30	3 Lead	Snap	127cm /50 in	25 Sets
700-0007-31	4 Lead (V-Lead)	Snap	127cm /50 in	25 Sets
700-0007-32	5 Lead	Snap	127cm /50 in	25 Sets
700-0007-33	5 Lead (V-Lead)	Snap	127cm /50 in	25 Sets
700-0007-34	6 Lead	Snap	127cm /50 in	25 Sets

Reusable, DIN Shielded, Multi-Pinch, AAMI



700-0006-35	3 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-36	4 Lead (V-Lead)	Multi-Pinch	74cm /29 in	1 Set
700-0006-37	5 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-38	5 Lead (V-Lead)	Multi-Pinch	74cm /29 in	1 Set
700-0006-39	6 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0010-00	3 Lead	Multi-Pinch	130cm/51 in	1 Set
700-0010-01	5 Lead	Multi-Pinch	130cm/51 in	1 Set
700-0010-02	5 Lead (V-Lead)	Multi-Pinch	130cm/51 in	1 Set

Reusable, DIN Shielded, Snap, AAMI



700-0007-35	3 Lead	Snap	74cm /29 in	1 Set
700-0007-36	4 Lead (V-Lead)	Snap	74cm /29 in	1 Set
700-0007-37	5 Lead	Snap	74cm /29 in	1 Set
700-0007-38	5 Lead (V-Lead)	Snap	74cm /29 in	1 Set
700-0007-39	6 Lead	Snap	74cm /29 in	1 Set

ECG COMBINER LEAD WIRE SYSTEM

AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Radiolucent Lead Wire Sets - requires ECG Combiner Cable or Adapter

Radiolucent, Reusable, DIN Shielded, Multi-Pinch, AAMI					
700-0010-03	3 Lead, Radiolucent	Multi-Pinch	74cm /29 in	1 Set	
700-0010-04	5 Lead, Radiolucent	Multi-Pinch	74cm /29 in	1 Set	
700-0010-05	3 Lead, Radiolucent	Multi-Pinch	122cm/48 in	1 Set	
700-0010-06	5 Lead, Radiolucent	Multi-Pinch	122cm/48 in	1 Set	
700-0010-07	3 Lead, Radiolucent	Multi-Pinch	183cm/72 in	1 Set	
700-0010-08	5 Lead, Radiolucent	Multi-Pinch	183cm/72 in	1 Set	



ECG Combiner Lead Wire Sets with Spacelabs Connection - for use with Command Module 90496/91496 and cable 700

Shielded, Spacelabs Connection, Multi-Pinch, Reusable, AAMI					
700-0010-16	5 Lead	Spacelabs Connection, Multi-Pinch	60 cm/24 in	1 Set	
700-0010-14	5 Lead	Spacelabs Connection, Multi-Pinch	102 cm/40 in	1 Set	



Shielded, Spacelabs Connection, Snap, Reusable, AAMI					
700-0010-15	5 Lead	Spacelabs Connection, Snap	60 cm/24 in	1 Set	



ECG COMBINER LEAD WIRE SYSTEM AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for use with Spacelabs command module 90496/91496				
700-0008-55	3 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-57	5 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-59	6 Lead	Spacelabs	305 cm/10 ft	1 Cable
700-0008-56	10 Lead (5+5)	Spacelabs	305 cm/10 ft	1 Cable
700-0008-58	10 Lead (6+4)	Spacelabs	305 cm/10 ft	1 Cable

ECG Combiner Adapters



Spacelabs Monitors, AAMI			
700-0008-78	3 Lead	Spacelabs	1 Adapter
700-0008-79	5 Lead	Spacelabs	1 Adapter
Spacelabs Telemetry, AAMI			
700-0008-80	5 Lead	Ultraview Telemetry (9034X/ 9134X)	1 Adapter
700-0008-88	5 Lead	Aria Telemetry (96281)	1 Adapter

ECG COMBINER LEAD WIRE SYSTEM AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization



Other manufacturers					
700-0008-30	3 Lead	Datascope	305 cm/10 ft	1 Cable	
700-0008-31	5 Lead	Datascope	305 cm/10 ft	1 Cable	
700-0008-33	3 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable	
700-0008-34	5 Lead	Datex-Ohmeda	305 cm/10 ft	1 Cable	
700-0008-35	6 lead	Datex-Ohmeda	305 cm/10 ft	1 Cable	
700-0008-36	3 Lead	GE	305 cm/10 ft	1 Cable	
700-0008-37	5 Lead	GE	305 cm/10 ft	1 Cable	
700-0008-38	10 Lead (5+5)	GE	305 cm/10 ft	1 Cable	
700-0008-39	6 Lead	GE	305 cm/10 ft	1 Cable	

ECG COMBINER LEAD WIRE SYSTEM

AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Cables for Standardization



Other manufacturers					
700-0008-40	10 Lead (6+4)	GE	305 cm/10 ft	1 Cable	
700-0008-47	3 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-48	10 Lead (5+5)	Philips	305 cm/10 ft	1 Cable	
700-0008-49	5 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-50	10 Lead (6+4)	Philips	305 cm/10 ft	1 Cable	
700-0008-51	6 Lead	Philips	305 cm/10 ft	1 Cable	
700-0008-41	3 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-42	5 Lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-43	6 lead	Medtronic LifePak	305 cm/10 ft	1 Cable	
700-0008-44	3 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-45	5 Lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-46	6 lead	Nihon Kohden	305 cm/10 ft	1 Cable	
700-0008-52	3 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-53	5 Lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	
700-0008-54	6 lead	Zoll, 90° angle	305 cm/10 ft	1 Cable	

ECG COMBINER LEAD WIRE SYSTEM AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.









Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers

	700-0008-60	5 Lead	Datascope Telemetry	1 Adapter
	700-0008-61	3 Lead	Datascope	1 Adapter
	700-0008-62	5 Lead	Datascope	1 Adapter
	700-0008-63	3 Lead	Datex-Ohmeda	1 Adapter
	700-0008-64	5 Lead	Datex-Ohmeda	1 Adapter
	700-0008-65	5 Lead	GE Apex Pro FH Telemetry	1 Adapter
	700-0008-66	6 Lead	GE Apex Pro FH Telemetry	1 Adapter
	700-0008-67	3 Lead	GE	1 Adapter
	700-0008-68	5 Lead	GE	1 Adapter

ECG COMBINER LEAD WIRE SYSTEM AAMI COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.

Additional cables, adapters and lead wires are also available to standardize most other manufacturer monitors and telemetry transmitters, allowing your hospital to use one set of lead wires per patient throughout the continuum of care.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Combiner Adapters for Standardization

Other manufacturers				
	700-0008-69	6 Lead	GE Monitors	1 Adapter
	700-0008-70	3 Lead	Nihon Kohden	1 Adapter
	700-0008-71	6 Lead	Nihon Kohden	1 Adapter
	700-0008-72	3/5 Lead	Philips Telemetry	1 Adapter
	700-0008-73	3 Lead	Philips	1 Adapter
	700-0008-74	5 Lead	Philips	1 Adapter
	700-0008-75	3 Lead	Siemens Drager	1 Adapter
	700-0008-76	5 Lead	Siemens Drager	1 Adapter

Shielded ECG lead wires and cables for use with Spacelabs
Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, Shielded, Reusable, AAMI



Snap

Snap Sets

700-0007-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0007-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0007-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0007-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0007-16*	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0007-18*	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220




Pinch

Pinch Sets

700-0006-00	3 Lead, Shielded	61 cm/24 in	1 Set
700-0006-02	3 Lead, Shielded	102 cm/40 in	1 Set
700-0006-08	5 Lead, Shielded	61 cm/24 in	1 Set
700-0006-10	5 Lead, Shielded	102 cm/40 in	1 Set
700-0006-16	10 Lead, limb/chest, Shielded	61 cm/24 in	1 Set
700-0006-18	10 Lead, limb/chest, Shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220

Shielded ECG lead wires and cables for use with Spacelabs
 Command Modules 90496 & 91496 and mCare Monitor 91220.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wires, Individual Colors, Shielded, Reusable, AAMI



Snap

24" Snap Individual Color

012-0498-00	Red	61 cm/24 in	3/Pack
012-0498-04	Black	61 cm/24 in	3/Pack
012-0498-34	Black	61 cm/24 in	5/Pack
012-0498-05	Brown	61 cm/24 in	3/Pack
012-0498-03	White	61 cm/24 in	3/Pack
012-0498-01	Green	61 cm/24 in	3/Pack
012-0498-11	Yellow/Brown	61 cm/24 in	3/Pack
012-0498-12	Green/Brown	61 cm/24 in	3/Pack
012-0498-13	Blue/Brown	61 cm/24 in	3/Pack
012-0498-14	Orange/Brown	61 cm/24 in	3/Pack
012-0498-15	Violet/Brown	61 cm/24 in	3/Pack



Snap

40" Snap Individual Color

012-0498-06	Red	102 cm/40 in	3/Pack
012-0498-09	Black	102 cm/40 in	3/Pack
012-0498-39	Black	102 cm/40 in	5/Pack
012-0498-27	Brown	102 cm/40 in	3/Pack
012-0498-08	White	102 cm/40 in	3/Pack
012-0498-07	Green	102 cm/40 in	3/Pack

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Cables, Shielded, AAMI

With 17 PIN Adapter



700-0008-08	3 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-06	5 Lead, Shielded	305 cm/10 ft	1 Cable
700-0008-00	10 Lead, Shielded	305 cm/10 ft	1 Cable

ECG Accessories

Electrodes



369062-001	ECG electrode, transparent tape, adult, single use	Adult	100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG electrode, transparent tape, general/stress, adult, single use	Adult	30/Box

Shielded ECG lead wires and cables for use with Spacelabs
 Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Neonatal ECG Supplies & Accessories

DIN Unshielded, Neonatal ECG Cable



700-0008-12	3 Lead		1 Cable
-------------	--------	--	---------

DIN Unshielded, Neonatal Electrodes



685-0037-00	Pre-wired neonatal electrodes, Pre-gelled, silver/silver chloride, single use	3-Lead	100 Sets/Box
-------------	---	--------	-----------------

Electrode



369063-001	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	100/Box
------------	---	---------	---------



690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	30/Box
-------------	---	---------	--------

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap

Snap Sets			
012-0283-00	3 Lead	46 cm/18 in	1 Set
012-0283-11	3 Lead	61 cm/24 in	1 Set
012-0283-02	3 Lead	102 cm/40 in	1 Set
012-0285-00	5 Lead	46 cm/18 in	1 Set
012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0285-02	5 Lead	102 cm/40 in	1 Set



Pinch

Pinch Sets			
012-0284-00	3 Lead	46 cm/18 in	1 Set
012-0284-11	3 Lead	61 cm/24 in	1 Set
012-0284-02	3 Lead	102 cm/40 in	1 Set
012-0286-00	5 Lead	46 cm/18 in	1 Set
012-0286-11	5 Lead	61 cm/24 in	1 Set
012-0286-02	5 Lead	102 cm/40 in	1 Set

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

Snap Individual Color			
012-0295-00	Red	46 cm/18 in	3/Pack
012-0295-10	Black	46 cm/18 in	3/Pack
012-0295-20	Brown	46 cm/18 in	3/Pack
012-0295-30	White	46 cm/18 in	3/Pack
012-0295-40	Green	46 cm/18 in	3/Pack
012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack
012-0295-02	Red	102 cm/40 in	3/Pack
012-0295-22	Black	102 cm/40 in	3/Pack
012-0295-42	Brown	102 cm/40 in	3/Pack
012-0295-12	White	102 cm/40 in	3/Pack
012-0295-32	Green	102 cm/40 in	3/Pack

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch


Pinch Individual Color			
012-0296-00	Red	46 cm/18 in	3/Pack
012-0296-20	Black	46 cm/18 in	3/Pack
012-0296-10	White	46 cm/18 in	3/Pack
012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-02	Red	102 cm/40 in	3/Pack
012-0296-22	Black	102 cm/40 in	3/Pack
012-0296-42	Brown	102 cm/40 in	3/Pack
012-0296-12	White	102 cm/40 in	3/Pack
012-0296-32	Green	102 cm/40 in	3/Pack

ECG Adapter and Cable, AAMI



Required for use with Multiparameter Module 90470		
103-0059-00	17-pin to 6-pin adapter, allows the ECG cables designed for Command Modules 90496/91496 to work on the Multiparameter Module 90470	1 Each

SUPPLIES & ACCESSORIES
TELEMETRY
AAMI Color Coding



AHA color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, AAMI



Snap

Snap Sets

012-0285-11	5 Lead	61 cm/24 in	1 Set
012-0605-10	5 Lead	64 cm/25.2 in	1 Set



Pinch

Pinch Sets

012-0286-11	5 Lead	61 cm/24 in	1 Set
-------------	--------	-------------	-------

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Snap

Snap Individual Color

012-0295-06	Red	61 cm/24 in	3/Pack
012-0295-26	Black	61 cm/24 in	3/Pack
012-0295-46	Brown	61 cm/24 in	3/Pack
012-0295-16	White	61 cm/24 in	3/Pack
012-0295-36	Green	61 cm/24 in	3/Pack

Snap, Individual Color

For use with Spacelabs Ultraview Telemetry 9034X/9134X only



Snap

700-0010-12	Red	74 cm/29 in	5/Pack
700-0010-09	Black	74 cm/29 in	5/Pack
700-0010-10	Brown	74 cm/29 in	5/Pack
700-0010-13	White	74 cm/29 in	5/Pack
700-0010-11	Green	74 cm/29 in	5/Pack

AHA color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, AAMI



Pinch Individual Color			
012-0296-04	Red	61 cm/24 in	3/Pack
012-0296-14	White	61 cm/24 in	3/Pack
012-0296-24	Black	61 cm/24 in	3/Pack
012-0296-34	Green	61 cm/24 in	3/Pack
012-0296-44	Brown	61 cm/24 in	3/Pack

Telemetry Transmitter Adapters



For use with ECG Combiner Lead Wires		
700-0008-80	Ultraview Telemetry Combiner Adapter (9034X/ 9134X)	1 Adapter



700-0008-88	Aria Telemetry Combiner Adapter (96281)	1 Adapter
-------------	---	-----------



For use with DIN Unshielded Lead Wires		
706-0191-00	Aria Telemetry Grouper (96281)	5/Bag

AHA color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Telemetry Accessories

Disposable Pouches and Electrodes



015-0500-02	Ultraview Telemetry pouch (90347)	100/Box	100/Box
-------------	-----------------------------------	---------	---------



369062-001	ECG electrode, transparent tape, adult, single use	Adult	100/Box
------------	--	-------	---------










690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
-------------	--	-------	--------



690-0262-00	ECG transparent electrode, general/stress	Adult	Adult	30/Box
-------------	---	-------	-------	--------

AHA color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
SpO₂ Cable & Sensors			
SpO₂ Cable			
 700-0014-00	TruLink SpO ₂ Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO ₂ sensors	60 cm/24 in	1 Cable
Reusable SpO₂ Sensors			
 015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
Disposable SpO₂ Sensors			
 015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
 015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
 015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
 015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
 015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box

AHA color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
ABP Cable			
ABP Cable			
700-0015-00	TruLink ABP Interface Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281)	60 cm/24 in	1 Cable

SUPPLIES & ACCESSORIES

ELECTROCARDIOGRAPHY

IEC Color Coding




ECG COMBINER LEAD WIRE SYSTEM

IEC COLOR CODED



The ECG Combiner System can standardize Spacelabs monitors and Ultraview & Aria Telemetry transmitters to use one set of lead wires.



Order Code Description Length Qty

ECG Combiner Lead Wire Sets - requires ECG Combiner Cable or Adapter



Disposable, DIN Shielded, Multi-Pinch, IEC

700-0006-42	3 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-44	5 Lead	Multi-Pinch	127cm /50 in	25 Sets
700-0006-33	5 Lead (V-Lead)	Multi-Pinch	127cm /50 in	25 Sets



Disposable, DIN Shielded, Snap, IEC

700-0007-42	3 Lead	Snap	127cm /50 in	25 Sets
700-0007-45	5 Lead	Snap	127cm /50 in	25 Sets
700-0007-33	5 Lead (V-Lead)	Snap	127cm /50 in	25 Sets



Reusable, DIN Shielded, Multi-Pinch, IEC

700-0006-41	3 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-43	5 Lead	Multi-Pinch	74cm /29 in	1 Set
700-0006-38	5 Lead (V-Lead)	Multi-Pinch	74cm /29 in	1 Set



Reusable, DIN Shielded, Snap, IEC

700-0007-41	3 Lead	Snap	74cm /29 in	1 Set
700-0007-44	5 Lead	Snap	74cm /29 in	1 Set
700-0007-38	5 Lead (V-Lead)	Snap	74cm /29 in	1 Set

ECG Combiner Cables



Spacelabs Connect Hardwired Cables for monitors using Command Module (90496/91496)

700-0008-81	3 Lead	IEC	305 cm/10 ft	1 Cable
700-0008-83	5 Lead	IEC	305 cm/10 ft	1 Cable
700-0008-84	10 Lead (5+5)	IEC	305 cm/10 ft	1 Cable

Shielded ECG lead wires and cables for use with Spacelabs
Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, Shielded, Reusable, IEC



Snap

Snap Sets

700-0007-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0007-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0007-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0007-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0007-17*	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0007-19*	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220



Pinch

Pinch Sets

700-0006-01	3 Lead, shielded	61 cm/24 in	1 Set
700-0006-03	3 Lead, shielded	102 cm/40 in	1 Set
700-0006-09	5 Lead, shielded	61 cm/24 in	1 Set
700-0006-11	5 Lead, shielded	102 cm/40 in	1 Set
700-0006-17	10 Lead, limb/chest, shielded	61 cm/24 in	1 Set
700-0006-19	10 Lead, limb/chest, shielded	102 cm/40 in	1 Set

*10 Lead ECG Lead Wire Sets not compatible with mCare Monitor 91220

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wires, Individual Colors, Shielded, Reusable, IEC



Snap


24" Snap Individual Color			
012-0698-00	Red	61 cm/24 in	3/Pack
012-0698-01	Green	61 cm/24 in	3/Pack
012-0698-03	White	61 cm/24 in	5/Pack
012-0698-04	Black	61 cm/24 in	3/Pack
012-0498-02	Yellow	61 cm/24 in	3/Pack
012-0498-21	Yellow/White	61 cm/24 in	3/Pack
012-0498-22	Green/White	61 cm/24 in	3/Pack
012-0498-24	Black/White	61 cm/24 in	3/Pack
012-0498-25	Violet/White	61 cm/24 in	3/Pack



Snap

40" Snap Individual Color			
012-0498-26	Yellow	102 cm/40 in	3/Pack
012-0698-06	Red	102 cm/40 in	3/Pack
012-0698-07	Green	102 cm/40 in	5/Pack
012-0698-08	White	102 cm/40 in	3/Pack
012-0698-09	Black	102 cm/40 in	3/Pack

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded

Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Cables, Shielded, IEC

With 17 PIN Adapter



700-0008-09	3 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-07	5 Lead, shielded	305 cm/10 ft	1 Cable
700-0008-01	10 Lead, shielded	305 cm/10 ft	1 Cable

ECG Accessories

Electrodes



369062-001	ECG electrode, transparent tape, adult, single use	Adult	100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG electrode, transparent tape, general/stress	Adult	30/Box

Shielded ECG lead wires and cables for use with Spacelabs Command Modules 90496 & 91496 and mCare Monitor 91220.



Shielded



Order Code	Description	Length	Qty
------------	-------------	--------	-----

Neonatal ECG Supplies & Accessories



DIN Unshielded, Neonatal ECG Cable

700-0008-13	3 Lead		1 Cable
-------------	--------	--	---------



DIN Unshielded, Neonatal Electrodes

685-0037-00	Pre-wired neonatal electrodes, single use Pre-gelled. Silver/Silver Chloride	3-Lead	100 Sets/Box
-------------	---	--------	--------------



Electrode

369063-001	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	100/Box
690-0290-00	ECG electrode, transparent tape, pediatric/neonatal, single use	Ped/Neo	30/Box

DIN Unshielded ECG lead wires and cables for use with Spacelabs Multiparameter Module 90470 ; available in 3- and 5-lead configurations.



Order Code	Description	Length	Qty
------------	-------------	--------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC



Snap

Snap Sets

012-0326-10	3 Lead	61 cm/24 in	1 Set
012-0326-01	3 Lead	102 cm/40 in	1 Set
012-0328-10	5 Lead	61 cm/24 in	1 Set
012-0328-01	5 Lead	102 cm/40 in	1 Set



Pinch

Pinch Sets

012-0325-10	3 Lead	61 cm/24 in	1 Set
012-0325-01	3 Lead	102 cm/40 in	1 Set
012-0327-10	5 Lead	61 cm/24 in	1 Set
012-0327-01	5 Lead	102 cm/40 in	1 Set


ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, IEC



Snap Individual Color

012-0295-07	Red	61 cm/24 in	3/Pack
012-0295-17	White	61 cm/24 in	3/Pack
012-0295-27	Black	61 cm/24 in	3/Pack
012-0295-38	Green	61 cm/24 in	3/Pack
012-0295-48	Yellow	61 cm/24 in	3/Pack
012-0295-05	Red	102 cm/40 in	3/Pack
012-0295-15	White	102 cm/40 in	3/Pack
012-0295-25	Black	102 cm/40 in	3/Pack
012-0295-35	Green	102 cm/40 in	3/Pack
012-0295-52	Yellow	102 cm/40 in	3/Pack

SUPPLIES & ACCESSORIES
TELEMETRY
IEC Color Coding



IEC color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

ECG Lead Wire Sets, DIN Unshielded, Reusable, IEC



Snap

Snap Sets

012-0328-10	5 Lead	61 cm/24 in	1 Set
-------------	--------	-------------	-------



Pinch

Pinch Sets

012-0327-10	5 Lead	61 cm/24 in	1 Set
-------------	--------	-------------	-------

ECG Lead Wires, Individual Colors, DIN Unshielded, Reusable, IEC



Snap

Snap Individual Color

012-0295-07	Red	61 cm/24 in	3/Pack
012-0295-17	White	61 cm/24 in	3/Pack
012-0295-27	Black	61 cm/24 in	3/Pack
012-0295-38	Green	61 cm/24 in	3/Pack
012-0295-48	Yellow	61 cm/24 in	3/Pack

For use with DIN Unshielded Lead Wires

706-0190-00	Aria Telemetry Grouper (96281), IEC		5/Bag
-------------	-------------------------------------	--	-------

IEC color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Telemetry Accessories








Disposable Pouches and Electrodes



015-0500-02	Ultraview Telemetry pouch (90347)		100/Box
369062-001	ECG electrode, transparent tape, adult, single use	Adult	100/Box
690-0289-00	ECG electrode, transparent tape, adult, single use	Adult	30/Box
690-0262-00	ECG transparent electrode, general/stress	Adult	Adult 30/Box

IEC color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
SpO₂ Cable & Sensors			
SpO₂ Cable			
 015-0014-00	TruLink SpO ₂ Sensor Adapter Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96281) and TruLink SpO ₂ sensors	60 cm/24 in	1 Cable
Reusable SpO₂ Sensors			
 015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
Disposable SpO₂ Sensors			
 015-0662-00	TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
 015-0663-00	TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
 015-0664-00	TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
 015-0665-00	TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
 015-0667-00	TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box


IEC color coded lead wires for use with Spacelabs Telemetry Transmitter models: 90341, 90343, 90347, 91341, 91343, 91347 and 96281.



Order Code	Description	Length/Size	Qty
ABP Cable			
ABP Cable			
700-0015-00	TruLink ABP Interface Cable Compatible with SL multi-parameter telemetry transmitters (90343/91343/96217)	60 cm/24 in	1 Cable

SUPPLIES & ACCESSORIES

CAPNOGRAPHY & MULTIGAS



The Nomoline™ multi-patient use adapter for sidestream monitoring offers a cost effective solution for both intubated and non-intubated patients. The sophisticated technology removes water from the sample line and reduces potential for water intrusion and cross contamination. Utilizing the Nomoline multi-patient use adapter in conjunction with single use CO₂ cannulas, we are able to offer a capnography solution at a much lower cost.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Sidestream Monitoring Supplies

Multi-patient use Nomoline Adapter - required for sidestream monitoring



103-0234-00	Sidestream Nomoline adapter with female luer lock. To be used with single patient use cannulas , sidestream sampling lines and T-adapters. T-adapters require a sampling line for use.	Multi-use	25/Box
-------------	--	-----------	--------



015-0683-00	Sidestream Nomoline adapter with sampling line	Single Use	25/Box
-------------	--	------------	--------

CO2 Nasal Cannula with male luer connector



704-0176-00	Adult	2 m/7 ft	25/Box
704-0177-00	Pediatric	2 m/7 ft	25/Box
704-0178-00	Infant	2 m/7 ft	25/Box
704-0184-00	Adult	3 m/10 ft	25/Box
704-0185-00	Pediatric	3 m/10 ft	25/Box



CO2 Nasal Cannula w/O2 Delivery and male luer connector

704-0181-00	Adult	2 m/7 ft	25/Box
704-0180-00	Pediatric	2 m/7 ft	25/Box
704-0179-00	Infant	2 m/7 ft	25/Box
704-0186-00	Adult	3 m/10 ft	25/Box



CO2 Oral/Nasal Cannula w/O2 Delivery and male luer connector

704-0183-00	Adult	2 m/7 ft	25/Box
704-0182-00	Pediatric	2 m/7 ft	25/Box
704-0187-00	Adult	3 m/10 ft	25/Box
704-0189-00	Pediatric	3 m/10 ft	25/Box

The Nomoline™ multi-patient use adapter for sidestream monitoring offers a cost effective solution for both intubated and non-intubated patients. The sophisticated technology removes water from the sample line and reduces potential for water intrusion and cross contamination. Utilizing the Nomoline multi-patient use adapter in conjunction with single use CO₂ cannulas, we are able to offer a capnography solution at a much lower cost.



(92516) Capno Pod (92517) Capnography Module (92518) Multi-Gas Module

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Sidestream Monitoring Supplies - Non-intubated patients



Emergency CO2 Nasal Cannula			
704-0193-00	Adult	2 m/7 ft	25/Box



Sampling Lines & Airway Adapters, single use			
704-0191-00	Adult CO2 Sample Line w/ T-Adapter	2 m/7 ft	25/Box
704-0192-00	Adult CO2 Sample Line w/ T-Adapter	3 m/10 ft	25/Box



166-7085-00	Sidestream sampling line with male luer lock on both ends	2 m/6.5 ft	25/Box
166-7121-00	Sidestream sampling line with male luer lock on both ends	3 m/10 ft	50/Box



704-0175-00	Sidestream airway adapter/T-Adapter. Requires a sampling line for use.		25/Box
-------------	---	--	--------

Mainstream Monitoring Supplies for use only with the Capnography Module 92517 on intubated patients



Sensor & Airway Adapters - required for mainstream monitoring			
010-1980-00	Mainstream sensor, reusable		1 Each



704-0173-00	Mainstream airway adapter, single use	Adult/Pediatric	25/Box
-------------	---------------------------------------	-----------------	--------







704-0174-00	Mainstream airway adapter, single use	Infant	10/Box
-------------	---------------------------------------	--------	--------

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.









(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream, Non-intubated, single patient use			
	CO₂ Nasal Cannula		
704-0010-00	Adult		10/Box
704-0012-00	Pediatric		10/Box
704-0011-00	Infant/Neonatal		10/Box
	CO₂ Nasal Cannula w/O₂ Delivery		
704-0014-00	Adult		10/Box
704-0016-00	Pediatric		10/Box
704-0015-00	Infant/Neonatal		10/Box
	CO₂ Oral/Nasal Cannula		
704-0017-00	Adult		10/Box
704-0019-00	Pediatric		10/Box
	CO₂ Oral/Nasal Cannula w/O₂ Delivery		
704-0020-00	Adult		10/Box
704-0022-00	Pediatric		10/Box

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
Sidestream, Intubated, single patient use			
Sidestream - short term up to 12 hours (patient non-humidified)			
	704-0023-00 Airway adapter + sampling line, ET tube > 4.0 mm	Adult/Ped	10/Box
	704-0024-00 Airway adapter + sampling line, ET tube ≤ 4.0 mm	Ped/Infant	10/Box
	704-0027-00 CO ₂ sampling line w/male luer <i>Recommended for heat moisture exchange filters</i>		10/Box
Sidestream - long term up to 120 hours (patient humidified)			
	704-0025-00 Airway adapter + sampling line with Nafion®, ET tube > 4.0 mm	Adult/Ped	10/Box
	704-0026-00 Airway adapter + sampling line with Nafion	Ped/Infant	10/Box
	704-0028-00 CO ₂ sampling line w/male luer and Nafion <i>Recommended for heat moisture exchange filters</i>		10/Box

The Dual Capnography Module (91517) provides both sidestream and mainstream EtCO₂ capabilities in a single module. Supplies available for all patient types.



(91517)
Dual Capnography Module

Order Code	Description	Size	Qty
------------	-------------	------	-----

Mainstream, for Intubated Patients Only

Mainstream Cable with Sensor



704-0001-00	CO ₂ sensor, reusable		1 Each
-------------	----------------------------------	--	--------



704-0006-00	Sensor cable holding clip		5/Pack
-------------	---------------------------	--	--------

Mainstream Airway Adapters, for use with CO₂ Sensor 704-0001-00



704-0005-00	CO ₂ airway adapter, single use	Adult	24/Box
-------------	--	-------	--------



704-0004-00	CO ₂ airway adapter, single use	Neonate	24/Box
-------------	--	---------	--------



704-0002-00	CO ₂ airway adapter, reusable	Adult	1 Each
-------------	--	-------	--------



704-0003-00	CO ₂ airway adapter, reusable	Neonate	1 Each
-------------	--	---------	--------

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)




(91518)

Multigas Analyzers

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Multigas Supplies & Accessories

Anesthetic Agent Samples Lines



015-0312-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------



015-0313-01	Anesthetic agent sample line with gas sampling tee and airway adapter, single use, Sampling tee: 22mm OD/15 mm ID x 15mm OD Sampling line: 1.22 mm/0.05 inches ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0046-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	305 cm/10 ft.	50/Box
-------------	--	---------------	--------



162-0047-01	Anesthetic agent sample line with a male luer connector on each end, single use, 1.22 mm/.048 in ID tubing (90518/91518)	213 cm/7 ft.	50/Box
-------------	--	--------------	--------

For use with Ultraview SL™ Multigas Analyzers (90518 & 91518).



(90518)



(91518)

Multigas Analyzers

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Multigas Supplies & Accessories

Water Traps, Reusable



717-0447-00	D-FEND water trap cartridge, Black. Replace at least once every two months or when occluded (91518)		10/Pack
-------------	--	--	---------



717-0448-00	D-FEND+ water trap cartridge, green Replace at least once every 24 hours or when occluded (91518)		10/Pack
-------------	--	--	---------


Water Traps, Single Use



378-0050-00	TruLink water trap (90518)		1 Each
-------------	----------------------------	--	--------

Multigas Supplies & Accessories

Nafion® Tubing



369014-003	Nafion tubing, single use Use it in-line with sampling line to reduce moisture and prevent clogging (90518)	61 cm/24 in	1 Each
------------	--	-------------	--------

Filter Cartridge



016-0428-00	Nafion tubing, single use Use it in-line with sampling line to reduce moisture and prevent clogging (90518)		1 Each
-------------	--	--	--------



050-1120-00	O ₂ sensor kit, micro-fuel cell (90518-B)		1 Each
-------------	--	--	--------

SUPPLIES & ACCESSORIES




PULSE OXIMETRY



For use with Spacelabs Command Modules (90496 & 91496-U).
Cables and connectors are color-coded blue.

Full line of oximetry products that support your needs. TruLink SpO₂ sensors are available in a wide range of sizes and shapes to fit all patients.











Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology			
Cables			
	700-0030-00 TruLink SpO ₂ cable	305 cm/10 ft.	1 Each
Reusable sensors			
	015-0660-00 TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
	015-0661-00 TruLink SpO ₂ Y multi-site sensor Includes: 1 Sensor 1 ear clip 10 adult finger/neonate feet/hand wraps 10 pediatric finger or infant toe adhesive wraps	Universal >30 kg/66 lbs <30 kg/66 lbs	1 Each

For use with Spacelabs Command Modules (90496 & 91496-U).
Cables and connectors are color-coded blue.

Full line of oximetry products that support your needs. TruLink SpO₂ sensors are available in a wide range of sizes and shapes to fit all patients.



Order Code	Description	Length/Size	Qty
Spacelabs TruLink Technology			
Disposable Sensors			
	015-0662-00 TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
	015-0663-00 TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
	015-0664-00 TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0665-00 TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0666-00 TruLink SpO ₂ cloth sensor	Infant 3-20 kg/7-44 lbs	24/Box
	015-0667-00 TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box
Dixtal Cables			
	700-0029-00 TruLink SpO ₂ cable compatible with Novametrix® sensors	30 cm/12 in.	1 Each
Dixtal Reusable Sensors			
	015-0321-01 Dixtal SpO ₂ universal sensor ear clip	Universal >30 kg/66 lbs or <30 kg/66 lbs	1 Each
Accessories			
	015-0330-00 SpO ₂ foam wrap	Universal	25/Pack

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
LNCS Masimo SET Technology			
LNCS Adapter Cable			
700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors (91496-M/91496-U)	305 cm/10 ft	1 Each
LNCS Reusable Sensors			
690-0247-00	Masimo LNCS DC-I, SpO ₂ digit sensor	Adult > 30 kg	1 Each
690-0248-00	Masimo LNCS DC-IP, SpO ₂ digit sensor	Pediatric 10-50 kg	1 Each
690-0249-00	Masimo LNCS TC-I, Tip-Clip SpO ₂ ear sensor	Adult > 30 kg	1 Each
690-0250-00	Masimo LNCS TF-I, Transflectance SpO ₂ forehead sensor	Adult > 30 kg	1 Each

For use with Ultraview SL™ Command Module (91496-M).
Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

LNCS Masimo SET Technology






LNCS Disposable Sensors



690-0251-00	Masimo LNCS Amtx, adhesive SpO ₂ sensor, 18 in.	Adult > 30 kg	20/Box
690-0252-00	Masimo LCNS Pdtx, adhesive SpO ₂ sensor, 18 in.	Pediatric 10 - 50 kg	20/Box
690-0258-00	Masimo LNCS Inf-L, adhesive SpO ₂ sensor, 3 ft.	Infant 3 kg > 20 kg	20/Box
690-0259-00	Masimo LNCS Neo-L, adhesive SpO ₂ sensor, 3 ft.	Neonatal < 3 kg or > 40 kg	20/Box






For use with Spacelabs Command Module (91496-M) .
Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
RD SET Masimo Technology			
RD SET Adapter Cable			
	010-1649-00 Masimo RD SET SpO2 Cable	12 ft.	1 Each
	010-2146-00 Masimo RD SET SpO2 Adapter Cable, LNCS to RD SET (Utilize existing LNOP sensors with new RD SET patient cable item no. 010-1649-00)	1.5 ft.	1 Each
RD SET Reusable Sensors			
	690-0267-00 Masimo RD SET DCI reusable SpO2 sensor	Adult > 30 kg	1 Each
	690-0270-00 Masimo RD SET DCI-P reusable SpO2 sensor	Pediatric 10-50 kg	1 Each
	690-0274-00 Masimo RD SET TC-I Tip-Clip reusable SpO2 sensor	Adult > 30 kg	20/Box

For use with Spacelabs Command Module (91496-M) .
Cables and connectors are color-coded gray.



Order Code	Description	Length/Size	Qty
RD SET MasimoTechnology			
RD SET Adhesive Disposable Sensors			
 RD SET Adt	690-0283-00 Masimo RD SET Adt, adhesive SpO2 sensor	Adult	20/Box
 RD SET Pdt	690-0284-00 Masimo RD SET Pdt, adhesive SpO2 sensor	Pediatric	20/Box
 RD SET Inf	690-0285-00 Masimo RD SET Inf, adhesive SpO2 sensor	Inant	20/Box
 RD SET Neo	690-0286-00 Masimo RD SET Neo, adhesive SpO2 sensor	Neo/Adult	20/Box
 RD SET NeoPt	690-0287-00 Masimo RD SET NeoPt, adhesive SpO2 sensor	Neonatal	20/Box

For use with Spacelabs Command Modules (90496,91496-N & 91496-U) and mCare Monitor 91220 and DM3 Monitor. Cables and connectors are color-coded purple.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Nellcor OxiMax Technology

Cable



700-0792-00	TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors (91496-N, 91220, 91330-N)		1 Each
-------------	---	--	--------



700-0030-00	TruLink SpO ₂ cable compatible with Nellcor OxiSensor II and OxiMax sensors (90496, 91496-U, 91330-S)		1 Each
-------------	--	--	--------

Reusable Sensors

Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien).



690-0003-01	Nellcor DS-100A, finger SpO ₂ sensor	Adult	1 Each
-------------	---	-------	--------



690-0004-01	Nellcor OXI-A/N, SpO ₂ sensor	Adult/Neonatal	1 Each
-------------	--	----------------	--------



690-0039-01	Nellcor OXI-P/I, SpO ₂ sensor	Pediatric/Infant	1 Each
-------------	--	------------------	--------



690-0215-01	Nellcor DS-Y, multi-site SpO ₂ sensor	Universal	1 Each
-------------	--	-----------	--------

SUPPLIES & ACCESSORIES TEMPERATURE



For use with Spacelabs Command Modules (90496 & 91496) when monitoring continuous temperature .

Full line of continuous core temperature monitoring products.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Continuous Temperature

Cables

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



700-0031-00	TruLink Dual Temperature Cable	30 cm/12 in.	1 Each
-------------	--------------------------------	--------------	--------



700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors	305 cm/10 ft	1 Each
-------------	--	--------------	--------

Reusable Probes/Sensors

Skin Surface



690-0028-01	YSI 700 series compatible with 6.5 mm DIA phone jack type connections	305 cm/10 ft	1 Each
-------------	---	--------------	--------

Esophageal / Rectal



690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------








118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each
-------------	--------------------------------------	--------------	--------

For use with Spacelabs Command Modules (90496 & 91496) when monitoring continuous temperature .

Full line of continuous core temperature monitoring products.



Order Code	Description	Length/Size	Qty
Single use Probes/Sensors - Requires YSI 400 series cable (700-4000-00)			
Skin			
	701-4000-00 Skin Temperature Sensor		20/Box
Esophageal/ Stethoscope			
	701-4009-00 9 French, YSI 400 series compatible		20/Box
	701-4012-00 12 French, YSI 400 series compatible		20/Box
	701-4018-00 18 French, YSI 400 series compatible		20/Box
	701-4024-00 24 French, YSI 400 series compatible		20/Box
Esophageal/ Rectal			
	701-4109-00 9 French, YSI 400 series compatible		20/Box
	701-4112-00 12 French, YSI 400 series compatible		20/Box
Foley catheter with temperature probe			
	701-4214-00 14 French, YSI 400 series compatible		20/Box
	701-4216-00 16 French, YSI 400 series compatible		20/Box
	701-4218-00 18 French, YSI 400 series compatible		20/Box
Tympanic temperature sensor			
	701-4300-00 Adult, YSI 400 series compatible		20/Box

SUPPLIES & ACCESSORIES

INVASIVE PRESSURE



For use with Spacelabs Command Modules (90496 & 91496) configured with the invasive pressure option. Option B supports two pressures; Option C supports four pressures.



(9X496-B) (9X496-C)
Command Modules

Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Integrated Adapter Cables (connects directly to transducer)

Edward Lifesciences



700-0293-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0292-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Abbott Laboratories/Hospira/ICU Medical



700-0295-01	Single pressure integrated cable	396 cm/13 ft	1 Each
-------------	----------------------------------	--------------	--------



700-0294-01	Dual pressure integrated cable	396 cm/13 ft	1 Each
-------------	--------------------------------	--------------	--------

Interface Cable & Adapters

Interface Cable



690-0021-00	Compatible with Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers. Requires adapter cable 700-0091-00 or 700-0028-00	366 cm/12 ft	1 Each
-------------	--	--------------	--------

Interface Adapter- Transducer interface cable required from transducer company



700-0791-00	Single pressure cable	122 cm/4 ft	1 Each
-------------	-----------------------	-------------	--------



700-0028-00	Dual pressure cable	122 cm/4 ft	1 Each
-------------	---------------------	-------------	--------

For use with Spacelabs Command Modules (90496 & 91496) configured with the invasive pressure option. Option B supports two pressures; Option C supports four pressures.



(9X496-B) (9X496-C)
Command Modules








Order Code	Description	Length/Size	Qty
Transducer Kits & Accessories			
Diaphragm Dome, Single Use			
690-0018-00	Diaphragm dome, single use Rotating luer lock fittings Compatible with Datex-Ohmeda TA101OND		12/Box
Pressure Transducer Kits, Single use			
690-0008-00	Transducer Kit, Patient Mount, DTX/Plus	152.5 cm/60 in	5/Box
690-0009-00	Transducer Kit, Pole Mount, DTX/Plus	152.5 cm/60 in	5/Box
690-0010-00	Transducer Kit, Pole Mount, DTX/Plus	182.9 cm/72 in	5/Box
690-0011-00	Pressure Transducer, Neonatal		10/Box
Pressure Transducer Accessories			
690-0012-00	Transducer Mounting Background		1 Each
690-0034-00	Transducer Holder, Reusable		1 Each

SUPPLIES & ACCESSORIES CARDIAC OUTPUT





For use with Spacelabs Command Modules (90496 & 91496) configured with the cardiac output option C.



Order Code	Description	Length/Size	Qty
Cardiac Output Cables & Probes			
Cables			
	700-0027-00 TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
	306655-003 Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
	306655-004 Cardiac Output cable compatible with Becton Dickinson/Datex-Ohmeda injectate temperature probe		1 Each
Injectate Probes			
	690-0023-00 In-line injectate temperature probe Compatible with Datex-Ohmeda SP4045		1 Each
	690-0024-00 Reference solution in-line injectate temperature probe. Compatible with Datex-Ohmeda SP5030		1 Each
	690-0033-00 In-line injectate temperature probe. Compatible with Datex-Ohmeda SP4042. Not compatible with Abbott or Baxter in-line injectate probes or delivery systems.		1 Each
	690-0043-01 In-line injectate sensor housing, single use, sterile. Compatible with Datex-Ohmeda SP5045		25/Box
Closed Loop Injectate Kits			
	690-0031-00 CritiKit SP4572 CO in-line room temperature injectate kit, closed-loop delivery system, single use, sterile		5/Box
	690-0032-00 CritiKit SP4500 CO in-line ice temperature injectate kit, closed-loop delivery system, single use, sterile		4/Box

SUPPLIES & ACCESSORIES RECORDING & PRINTING







High quality Z-fold paper for use with the Bedside Printer Module (90449) and Central Printer (90469).

Roll paper for use with the Standalone printer (91449) and Spacelabs Portable Monitors.



(91449)
Printer

Order Code	Description	Length/Size	Qty
Thermal Paper			
Z-fold Paper for use with 90449 and 90469			
 006-0196-00	TruLink Thermal Paper, Z-Fold Pads, 2-channel printing, blank paper, 230 Sheets/pad	27.4 m/90 ft	10 Pads/Box
 307568-001	TruLink Thermal Paper, Z-fold Pads, 2- or 4-channel printing, blank paper, 300 sheets/pad	45.7 m/150 ft	10 Pads/Box
Roll Paper for use with 91449, 91367, 91369, 91370 & 91390			
 305240-001	TruLink Thermal Paper, Rolls, single-channel printing, 40 mm/1.575 inch wide preprinted grid	61 m/200 ft	10 Rolls/Box
 307438-001	TruLink Thermal Paper, Rolls, single- or dual-channel printing, 50 mm/1.97 inch wide blank paper	30.1 m/100 ft	10 Rolls/Box

SUPPLIES & ACCESSORIES BISPECTRAL INDEX



For use with Spacelabs BISx™ Module (91482).






Order Code	Description	Qty
------------	-------------	-----

BISx Cables & Sensors

BISx Cables		
719-0004-01	BISx Patient Interface cable	1 Each
719-0006-20	BISx Pod	1 Each
719-0007-01	BISx host monitor cable	1 Each

BISx Sensors

Note: Sales of BISx sensors are limited to one box with the purchase of capital equipment. All future purchases should be made directly through Covidien.

719-0003-10	BISx Extend Sensor, single patient use	25/Box
		
719-0001-10	BISx Pediatric Sensor, single patient use	25/Box
		
719-0002-10	BISx Quatro Sensor, single patient use	25/Box
		

SUPPLIES & ACCESSORIES

ULTRAVIEW DM3 DUAL MODE VITAL SIGNS MONITOR



ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

TruLink Nylon Reusable BP Cuffs



Single Tube				
715-1121-10	Single tube, HP	Child	13 - 20 cm	1 Each
715-1131-10	Single tube, HP	Small Adult	18 - 26 cm	1 Each
715-1141-10	Single tube, HP	Adult	26 - 35 cm	1 Each
715-1151-10	Single tube, HP	Adult Long	29 - 38 cm	1 Each
715-1161-10	Single tube, HP	Large Adult	32 - 42 cm	1 Each
715-1171-10	Single tube, HP	Large Adult Long	35 - 44 cm	1 Each
715-1181-10	Single tube, HP	Thigh Cuff	42 - 50 cm	1 Each

TruLink Soft Disposable BP Cuffs



Single Tube				
715-1123-10	Single tube, HP	Child	13 - 20 cm	5/Bag
715-1133-10	Single tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1143-10	Single tube, HP	Adult	26 - 35 cm	5/Bag
715-1153-10	Single tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1163-10	Single tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1173-10	Single tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1183-10	Single tube, HP	Thigh Cuff	42 - 50 cm	5/Bag

TruLink Vinyl Disposable BP Cuffs



Single Tube				
715-1122-10	Single Tube, HP	Child	13 - 20 cm	5/Bag
715-1132-10	Single Tube, HP	Small Adult	18 - 26 cm	5/Bag
715-1142-10	Single Tube, HP	Adult	26 - 35 cm	5/Bag
715-1152-10	Single Tube, HP	Adult Long	29 - 38 cm	5/Bag
715-1162-10	Single Tube, HP	Large Adult	32 - 42 cm	5/Bag
715-1172-10	Single Tube, HP	Large Adult Long	35 - 44 cm	5/Bag
715-1182-10	Single Tube, HP	Thigh Cuff	42 - 50 cm	5/Bag



Single Tube Hose				
714-0018-00	Single tube	Adult	275 cm/9 ft	1 Each
714-0060-00	Single tube	Adult	366 cm/12 ft	1 Each

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

Spacelabs SpO₂ TruLink Technology - 91330 with SPO2 option S

Cables



700-0030-00	TruLink SpO ₂ cable (91330-S)		305 cm/10 ft.	1 Each
-------------	--	--	---------------	--------

Reusable sensors



015-0660-00	TruLink SpO ₂ finger sensor		Adult >30 kg/66 lbs	1 Each
-------------	--	--	------------------------	--------



015-0661-00	TruLink SpO ₂ Y multi-site sensor Includes: 1 Sensor, 1 ear clip, 10 adult finger/neonate feet/hand wraps and 10 pediatric finger or infant toe adhesive wraps		Universal >30 kg/66 lbs <30 kg/66 lbs	1 Each
-------------	--	--	---	--------

Disposable Sensors



015-0662-00	TruLink SpO ₂ foam sensor		Adult > 30 kg/66 lbs	24/Box
-------------	--------------------------------------	--	-------------------------	--------



015-0663-00	TruLink SpO ₂ vinyl sensor		Adult >30 kg/66 lbs	24/Box
-------------	---------------------------------------	--	------------------------	--------



015-0664-00	TruLink SpO ₂ foam sensor		Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	--------------------------------------	--	----------------------------------	--------



015-0665-00	TruLink SpO ₂ vinyl sensor		Pediatric 10-50 kg/22-110 lbs	24/Box
-------------	---------------------------------------	--	----------------------------------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

LNCS Masimo SET SpO₂ Technology - 91330 with SPO2 Option M

LNCS Adapter Cable



700-0906-01	TruLink SpO ₂ cable compatible with Masimo LNCS sensors		305 cm/10 ft	1 Each
-------------	--	--	--------------	--------

LNCS reusable sensors



690-0247-00	Masimo LNCS DC-I Digit Sensor		Adult > 30 kg	1 Each
-------------	-------------------------------	--	------------------	--------



690-0248-00	Masimo LNCS DC-IP Pediatric Digit Sensor		Pediatric 10-50 kg	1 Each
-------------	--	--	-----------------------	--------



690-0249-00	Masimo LNCS TC-I Tip-Clip Ear Sensor		Adult > 30 kg	1 Each
-------------	--------------------------------------	--	------------------	--------



690-0250-00	Masimo LNCS TF-I Transflectance Forehead Sensor		Adult > 30 kg	1 Each
-------------	---	--	------------------	--------

LNCS disposable sensors



690-0251-00	Masimo LNCS Amtx adhesive sensor, 18 in.		Adult > 30 kg	20/Box
-------------	--	--	------------------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

Masimo RD SET SpO₂ Technology - 91330 with SPO2 Option M

RD SET Adapter Cable

010-1649-00	Masimo RD SET SpO ₂ Cable		12 ft.	1 Each
010-2146-00	Masimo RD SET SpO ₂ Adapter Cable, RD SET to RD SET		1.5 ft.	1 Each

RD SET Reusable Sensors



690-0267-00	Masimo RD SET DCI reusable SpO ₂ sensor		Adult > 30 kg	1 Each
-------------	--	--	------------------	--------



690-0270-00	Masimo RD SET DCI-P reusable SpO ₂ sensor		Pediatric 10-50 kg	1 Each
-------------	--	--	-----------------------	--------



690-0274-00	Masimo RD SET TC-I Tip-Clip reusable SpO ₂ sensor		Adult > 30 kg	20/Box
-------------	--	--	------------------	--------

RD SET Adhesive Disposable Sensors



690-0283-00	Masimo RD SET Adt, Adult SpO ₂ adhesive Sensor		Adult	20/Box
-------------	---	--	-------	--------



690-0284-00	Masimo RD SET Pdt, Pediatric SpO ₂ adhesive Sensor		Pediatric	20/Box
-------------	---	--	-----------	--------



690-0285-00	Masimo RD SET Inf, Infant SpO ₂ adhesive Sensor		Infant	20/Box
-------------	--	--	--------	--------



690-0286-00	Masimo RD SET Neo, Neonatal/Adult SpO ₂ adhesive Sensor		Neo/Adult	20/Box
-------------	--	--	-----------	--------



690-0287-00	Masimo RD SET NeoPt, Neonatal SpO ₂ adhesive Sensor		Neonatal	20/Box
-------------	--	--	----------	--------

ULTRAVIEW DM3 VITAL SIGNS MONITOR SUPPLIES & ACCESSORIES



Ultraview DM3 supports measurement of SpO₂, pulse rate, NIBP and temperature, and can be used in any patient care environment where basic vital signs monitoring is needed.



91330

Order Code	Description	Size	Length/Size	Qty
------------	-------------	------	-------------	-----

Nellcor OxiMax SpO₂ Technology - 91330 with SPO2 Option N

Cable



700-0792-00	TruLink SpO ₂ cable compatible with Nellcor OxiMax sensors			1 Each
-------------	---	--	--	--------

Reusable Sensors

Note: Nellcor sensors are available for purchase through Spacelabs Healthcare only with the purchase of capital equipment. All future purchases should be made directly through Nellcor(Covidien).



690-0003-01	Nellcor OxiMax DS-100A finger sensor		Adult	1 Each
-------------	--------------------------------------	--	-------	--------



690-0004-01	Nellcor OXI-A/N sensor		Adult/Neonatal	1 Each
-------------	------------------------	--	----------------	--------



690-0039-01	Nellcor OXI-P/I sensor		Pediatric/Infant	1 Each
-------------	------------------------	--	------------------	--------



690-0215-01	Nellcor DS-Y multi-site sensor		Universal	1 Each
-------------	--------------------------------	--	-----------	--------

Temperature

Probes

019-0539-00	Replacement Temperature Probe, Oral/Axillary			1 Each
701-0122-00	Disposable Temperature Probe Covers			200/Box

SUPPLIES & ACCESSORIES

AMBULATORY BLOOD PRESSURE MONITOR



Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Reusable ABP Cuffs

UltraCheck Reusable Nylon Cuffs



016-1320-01	Self-applied, reusable	child	13 - 20 cm	1 Each
016-1826-01	Self-applied, reusable	small adult	18 - 26 cm	1 Each
016-2635-01	Self-applied, reusable	adult	26 - 35 cm	1 Each
016-3242-01	Self-applied, reusable	large adult	32 - 42 cm	1 Each

Spacelabs Reusable Cloth Cuffs



015-0118-04Q	Self-applied with quick connect rectus fitting	child	12 - 20 cm	1 Each
015-0067-04Q	Self-applied with quick connect rectus fitting	small adult	17 - 26 cm	1 Each
015-0068-05Q	Self-applied with quick connect rectus fitting	adult	24 - 32 cm	1 Each
016-0077-05Q	Self-applied with quick connect rectus fitting	large adult	32 - 42 cm	1 Each
016-0109-04Q	Self-applied with quick connect rectus fitting	XL adult	38 - 50 cm	1 Each

Disposable Cuff Covers

For use with Reusable ABP Cuffs

706-0182-00	Universal cuff cover	Standard	50/Box
706-0183-00	Universal cuff cover	Large	50/Box

Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

Disposable ABP Cuffs



Disposable Cuffs				
015-1826-00	Single tube, no fitting	Small Adult	18 - 26 cm	20/Box
015-2635-00	Single tube, no fitting	Adult	26 - 35 cm	20/Box
015-3242-00	Single tube, no fitting	Large Adult	32 - 42 cm	20/Box
015-4250-00	Single tube, no fitting	X-large Adult	42 - 50 cm	20/Box
712-1039-00	Quick Release Metal Fittings, female bayonet. For use with above single tube, disposable cuffs with no fitting.			3/Pack

ABP Accessories



Pouches - For use with APB Monitor 90207				
016-0263-00	Vinyl carrying pouch with Velcro closure			Each
016-0410-00	Cloth carrying pouch with Velcro closure			Each

Pouches - For use with APB Monitor 90217				
016-0340-00	Carrying pouch, cloth			Each

Pouches - For use with APB Monitor 90227				
706-0200-00	Disposable Pouch			50/Pack
015-0715-00	Reusable Pouch			Each
015-0716-00	Reusable Shoulder Strap (for use with Pouch 015-0715-00)			Each



Ambulatory blood pressure is rapidly becoming a leading diagnostic tool. It compensates for White Coat Syndrome and allows for a series of measurements rather than relying on one isolated reading.



(90207)



(90217)

Order Code	Description	Size	Range	Qty
------------	-------------	------	-------	-----

ABP Accessories

Support



015-0070-00	Cuff support harness			Each
-------------	----------------------	--	--	------



016-0080-00	Nylon web adjustable waist belt, d-ring adjustable fastener			Each
-------------	---	--	--	------



016-0262-00	Nylon web adjustable shoulder strap, Velcro fasteners For carrying the pouch and monitor over the shoulder			Each
-------------	---	--	--	------



344-0056-00	Garment clip with vinyl strap, 3.5"			Each
-------------	-------------------------------------	--	--	------

Patient Diaries and Paper




000-0027-02	ABP Patient Diary, English, 24 hour Monitoring			Each
-------------	--	--	--	------



006-0210-02	ABP Thermal paper, 57 mm x 32 ft.			6 Rolls/ Pack
-------------	-----------------------------------	--	--	---------------

SUPPLIES & ACCESSORIES

HOLTER



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Holter Hook-up Kits

Aria Hook-Up Kits			
23189-106-401	Aria hook-up kit with case, pouch and clear tape	Aria	20/Case
23189-106	4 electrode hook-up kit	Aria, evo	Each
23189-110-401	Hook-up Kit (no pouch)	Aria, evo	20/Case
50-0048	4 electrode hook-up kit	evo	
50-0049	4 electrode hook-up kit with disposable chest pouch	evo	
50-0190	4 electrode hook-up kit	evo	20/Case
50-0191	4 electrode hook-up kit with chest pouch	evo	20/Case

Lead Wires & Electrodes

Lead Wires			
25452-525	4-Lead wire set, snap, three 14" leads and one 7" white lead	Aria, evo	4LW/set
25452-535	4-Lead wire set, snap, three 14" leads and one 14" white lead	Aria, evo	4LW/set
25452-529	4-Lead wire set, snap, 40" leads, Aria and evo belt clip configuration	Aria, evo	4LW/set
25452-519	7-Lead wire set, snap, 18" (white/black/red/green/orange/blue/brown)	423, 423A, 459, 461, 461A, 464, 483, 485	7 LW/set
19-7509	6-Lead wire set, snap, 24 in.	Lifecard CF	6LW/Set
19-7510	7-lead wire set, 24 in.	Lifecard CF	7LW/Set
25450-187	Individual color lead wire, snap, white, 14 in.	Aria, evo	




For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Lead Wires & Electrodes

Electrodes



27565-9	ECG Electrode, clear tape	Holter	4/Pack
27565-9CASE	ECG Electrode, clear tape	Holter	150/Pack
717-0071-00	ECG Electrode, AMBU	Holter	5/Pack
27565-5	ECG Electrode, foam	Holter	5/Pack
27565-7	ECG Electrode, foam	Holter	7/Pack

Patient Cables, Pouches, Shoulder Straps & Belts

Patient Cables



700-0316-00	Patient cable, 3-Channel, 6-lead, 22 in.	Lifecard CF	Each
700-0317-00	Patient cable with Belt Clip, 3-Channel, 3-lead, 41 in.	Lifecard CF	Each
700-0318-00	Patient cable, 3-channel	Lifecard CF	Each
700-0319-00	Patient cable with Neck Pouch, 3-Channel, 3-lead, 26 in.	Lifecard CF	Each
700-0320-00	Patient cable, 2-Channel	Lifecard CF	Each
700-0321-00	Patient cable, 2-Channel, short	Lifecard CF	Each
700-0322-00	Patient cable, 3-Channel, 7-lead, 23 in.	Lifecard CF	Each
700-0323-00	Patient cable with Beth Clip, 2-Channel, 4-lead, long	Lifecard CF	Each
700-0324-00	Patient cable with Neck Pouch, 2-Channel, 4-lead, short	Lifecard CF	Each

For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Lifecard CF



evo



Aria

Order Code	Description	For Use With	Qty
Patient Cables, Pouches, Shoulder Straps & Belts			
Patient Cables			
700-0325-00	Lifecard 12 varios Yoke	Lifecard 12	Each
700-0326-00	Lifecard 12 varios Yoke, including 12-lead patient cable, IEC	Lifecard 12	Each
700-0327-00	Lifecard 12 varios Yoke, including 12-lead patient cable, AHA	Lifecard 12	Each
700-0328-00	Lifecard 12 patient cable, IEC	Lifecard 12	Each
700-0329-00	Lifecard 12 patient cable, AHA	Lifecard 12	Each
700-0330-00	Lifecard 12 patient cable, long, IEC	Lifecard 12	Each
700-0331-00	Lifecard 12 patient cable, long, AHA	Lifecard 12	Each
Pouches, Shoulder Straps & Belts			
27249-7-401	Aria disposable pouch	Aria	50/Pack
27380	Aria belt clip (replacement)	Aria	Each
50-0042	evo holder for waist belt	evo	Each
PK0001	Fabric waist belt	evo	Each
50-0050	Disposable chest pouch for evo recorder	evo	50/Pack
706-0189-00	Disposable pouch for evo	evo	100/Pack
46-0408	Fabric pouch set, 2 pouches with laces, cordlock & instructions	Lifecard CF	Each
19-7130	Laces for pouch set	Lifecard CF	Each
19-7131	cordlock for pouch set	Lifecard CF	Each
46-0092	Fabric pouch	Lifecard CF	Each



For use with Spacelab's Digital Holter Recorders
Lifecard CF, evo & Aria.



Order Code	Description	For Use With	Qty
------------	-------------	--------------	-----

Patient Cables, Pouches, Shoulder Straps & Belts

Pouches, Shoulder Straps & Belts



46-0512	Disposable pouch	Lifecard CF	100/Pack
27685-401	Disposable pouch	LC12	50/Pack
200-1508-001	Pouch & straps	485	Each

Recording

Recorder Memory

CFC90MB	Flashcard Memory, 90MByte	Lifecard CF	Each
CFC256MB	Flashcard LifeCard CF, 256MB	LC12	Each
46-0063	Plastic memory card case	CFCxxMB	Each
50-0220	LCF / evo pen	LCF / evo	Each
19-7118	Lifecard Compact USB Flashcard Reader	CFCxxMB	Each
27311	Aria download cable	Aria	Each
46-0413	LCF opener tool	LCF	Each

Patient Diaries

Patient Diaries

073-0267-00	Patient Diary, English CE0123	Aria, evo, Lifecard CF	Each
46-0459	Patient Diary, 7 day, English	Lifecard CF	Each
46-0573	Patient Diary, English, Italian, French, Spanish, Polish	Lifecard CF	Each
46-0435	Patient diary, 7 day, German	Lifecard CF	Each

SUPPLIES & ACCESSORIES GENERIC & SUNTECH





Order Code	Description	For Use With	Qty
Generic Hook-up Kits			
DC Hook-up			
13-9083	Razor with blade cover	Holter / Event	Each
392003-001	Skin prep pad		100/Box
Stress Test			
98-0002-50	Microphone pads, 5 PK	Suntech Tango	5/Pack
98-0006-00	K-Sound microphone	Suntech Tango	Each
98-0061-01	Cuff, Orbit-K, small adult	Suntech Tango	Each
Recorder Memory			
41-0080	CD12USB stress tets hook-up kit	CD12-S	Each

SUPPLIES & ACCESSORIES

12-LEAD ECG





CD12



CardioExpress

Order Code Description For Use With Qty

12-Lead ECG Patient Cables & Hook-up Sets

MS1R-107581-A0	ECG cable, snap connector, IEC	CardioExpress	Each
MS1R-107582-A0	ECG cable, snap connector, AHA	CardioExpress	Each
MS1R-107583-A0	ECG cable, grabber connector, IEC	CardioExpress	Each
MS1R-107584-A0	ECG cable, grabber connector, AHA	CardioExpress	Each
MS1R-106902-A0	ECG cable, 4MM, banana connector, IEC	CardioExpress	Each
MS1R-107048-A0	ECG cable, 4MM, banana connector, AHA	CardioExpress	Each
25-0027	10 Electrode patient cable, short, IEC	CD12	Each
25-0028	10 Electrode patient cable, short, AHA	CD12	Each
25-0004	10 Electrode patient cable, 4MM, IEC	CD12	Each
25-0005	10 Electrode patient cable, 4MM, AHA	CD12	Each
25-0020	10 Electrode patent cable	CD12	Each
41-0053EU	12-lead hook-up set, IEC	Generic	
41-0053US	12-lead hook-up set, AHA	Generic	

Miscellaneous 12-Lead ECG Accessories

Paper		
0315BOX	Nikotab Electrodes	500/Box
717-0072-00	Adult Foam Electrodes w/snap	25/Pack
25-0010	Limb Clamps, 4/set	4/Set
58100	Easiclip, Skintact, 10/pk	10/Pack
MS1R-107371-A0	Recording paper, Z-fold, 210 MM x 295 MM	6 Pcs/ Set
MS1R-32462	Recording paper, Z-fold, 210MM x 295MM	Each
MS1R-107560-A0	Recording paper, Z-fold, 110MM*140MM	
MS1R-107451-A0	Recording paper, Z-fold, 215MM x 280MM	
MS1R-32461	Recording paper, Roll, 210 MM x 30M	Each
MS1R-32463	Recording paper, Roll, 215 MM x 30M	Each
MS4R-111897	Paper roller, plastic, SL12	

SUPPLIES & ACCESSORIES

EVENT MONITORING



OnTrak

Initialization

Start New Test

Start	Interval	Tone
08	5 m	<input checked="" type="checkbox"/>
22	5 m	<input checked="" type="checkbox"/>

Display Results

Display Cuff Press.

Child Mode

Comfort Mode: 170

- Advanced Setup

SPACE LABS
HEALTHCARE



Order Code	Description	For Use With	Qty
Event Monitoring Supplies & Accessories			
Hook-Up Kits			
23-0255	Lanyard clip & assembly	CardioCall	Each
Lead Wires & Electrodes			
23-0059	Patient cable, 3 electrodes, 2	CardioCall	Each
23-0051	Patient cable, 2 electrodes, 1	CardioCall	Each
Programming & Download			
23-0245	Programing card & key	CCVS20	Each
23-0252	Programing card & key	CCST80	Each
23-0026	Download cable	CardioCall	Each

SUPPLIES & ACCESSORIES

ANESTHESIA



For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleaseFocus



BleaseSirius



ARKON

Order Code	Description	Qty
------------	-------------	-----

Anesthesia Supplies

Nomoline Gas Sampling

For use with ARKON Anesthesia Delivery System only



103-0234-00	Nomoline adapter with female Luer lock, multi-patient use	25/Box
166-7085-00	Sidestream sampling line with male Luer lock on both ends	2m/6.5 ft 25/Box
166-7121-00	Sidestream sampling line with male Luer lock on both ends	3m.10 ft 50/Box

CO₂ Canisters

For use with ARKON Anesthesia Delivery System only



015-0702-00	CO2 Canister w/Spherasorb, white to violet	6/Case
015-0701-00	CO2 Canister w/LoFloSorb, green to violet	6/Case

Patient Flow Sensors

For use with BleaseSirius and BleaseFocus Anesthesia Machines only



10110090	Adult Patient flow sensor	1 Each
54100063	Flow sensor caps	50/Box
376-0561-01	Pediatric Patient flow sensor	1 Each
10110219	1.5M Patient flow sensor tubing assembly	1 Each
10110218	2M Patient flow sensor tubing assembly	1 Each
10110122	2.5M Patient flow sensor tubing assembly	1 Each
10110102	3M Patient flow sensor tubing assembly	1 Each
10110103	6M Patient flow sensor tubing assembly	1 Each
10110104	1.8M Patient flow sensor tubing assembly kit, single ended	1 Each
10110106	1.8M Patient flow sensor & tubing assembly kit, single ended	1 Each
70300001	O2 probe with tee	1 Each

For use with Spacelabs ARKON Anesthesia Delivery System and BleaseSirius and BleaseSirius Anesthesia Machines.



BleaseFocus



BleaseSirius



ARKON

Order Code	Description	Qty
------------	-------------	-----

Anesthesia Supplies - not anesthesia machine specific

Reusable Breathing Circuits

M1050-4	Reusable patient "Y" circuit, adult	1 Each
ST7557	Magill breathing circuit w/face mask	1 Each
ST2111	2.0 L rebreathing bag, 22mm neck	1 Each
ST2670	Malleable face mask, size 4	1 Each
2100	Jackson-Rees pediatric circuit	1 Each
2150	Rendell Baker pediatric mask set, 4 sizes	4/Set
42100	Malleable face mask set, 7 sizes	7/Set

Disposable Breathing Circuits

54100035	Pediatric Y piece breathing circuit 2010-022	1 Each
54100036	Infant T piece mapleson f breathing circuit 2120	1 Each
54100031	1.8M Magill-Mapleson, open circuit	1 Each
54100032	Individual deluxe Bain coaxial circuit 2115	1 Each
54100033	1.6M Basic Y piece breathing circuit 200	1 Each
10400001	Basic Y piece breathing circuit w/limb and bag 2010	1 Each

Disposable Face Masks

54100040	Infant disposable face mask	1 Each
54100039	Small disposable face mask	1 Each
54100038	Medium disposable face mask	1 Each
54100037	Large disposable face mask	1 Each

Disposable Reservoir Bags

54100041	0.5 L anesthetic reservoir bag, 15mm female neck	1 Each
54100042	1.0 L anesthetic reservoir bag, 22mm female neck	1 Each
54100043	2.0 L anesthetic reservoir bag, 22mm female neck	1 Each
54100044	3.0 L anesthetic reservoir bag, 22mm female neck	1 Each

Connectors & Filters

54100026	22 mm connector - 30f scavenging 1971	1 Each
54100027	Filta-Guard 1944 standard bacterial/viral breathing filter,	1 Each

SUPPLIES & ACCESSORIES

ELANCE



The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.



Order Code Description Qty

UltraCheck Nylon Reusable Blood Pressure Cuffs



NIBP Cuffs					
714-4001-00	Single tube, MHP	Child		13 - 20 cm	1 Each
714-4002-00	Single tube, MHP	Small Adult		18 - 26 cm	1 Each
714-4003-00	Single tube, MHP	Adult		26 - 35 cm	1 Each
714-4004-00	Single tube, MHP	Adult Long		29 - 38 cm	1 Each
714-4005-00	Single tube, MHP	Large Adult		32 - 42 cm	1 Each
714-4006-00	Single tube, MHP	Large Adult Long		35 - 44 cm	1 Each
714-4007-00	Single tube, MHP	Thigh Cuff		42 - 50 cm	1 Each

Hose					
------	--	--	--	--	--



714-0018-00	NIBP hose, single tube	Adult		275 cm/9 ft	1 Each
-------------	------------------------	-------	--	-------------	--------

The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.









Order Code	Description		Qty
ECG Lead Wires & Cables			
ECG Integrated Cables and Lead Wires			
700-0008-22	5-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
700-0008-23	5-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
700-0008-24	3-lead Wire Set, Snap, Shielded, AAMI Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
700-0008-25	3-lead Wire Set, Snap, Shielded, IEC Color Coding, with Integrated Cable, Shielded	61cm/10in 305 cm/10 ft	1 Each
ECG Cables			
700-0008-06	Shielded ECG cable, 5-lead, AAMI	305 cm/10 ft	1 Each
700-0008-07	Shielded ECG cable, 5-lead, IEC	305 cm/10 ft	1 Each
700-0008-08	Shielded ECG cable, 3-lead, AAMI	305 cm/10 ft	1 Each
700-0008-09	Shielded ECG cable, 3-lead, IEC	305 cm/10 ft	1 Each
3-Lead ECG Lead Wire Sets, Snap			
700-0007-00	Snap shielded ECG lead wire set, 3-lead, AAMI	61cm/24in	1 Set
700-0007-02	Snap shielded ECG lead wire set, 3-lead, AAMI	102cm/40in	1 Set
700-0007-01	Snap shielded ECG lead wire set, 3-lead, IEC	61cm/24in	1 Set
700-0007-03	Snap shielded ECG lead wire set, 3-lead, IEC	102cm/40in	1 Set









The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.



Order Code	Description	Qty
ECG Lead Wires & Cables		
5-Lead ECG Lead Wire Sets, Snap		
	700-0007-08 Snap shielded ECG lead wire set, 5-lead, AAMI	61cm/24in 1 Set
	700-0007-10 Snap shielded ECG lead wire set, 5-lead, AAMI	102cm/40in 1 Set
	700-0007-09 Snap shielded ECG lead wire set, 5-lead, IEC	61cm/24in 1 Set
	700-0007-11 Snap shielded ECG lead wire set, 5-lead, IEC	102cm/40in 1 Set
3-Lead ECG Lead Wire Sets, Pinch		
	700-0006-00 Pinch shielded ECG lead wire set, 3-lead, AAMI	61cm/24in 1 Set
	700-0006-02 Pinch shielded ECG lead wire set, 3-lead, AAMI	102cm/40in 1 Set
	700-0006-01 Pinch shielded ECG lead wire set, 3-lead, IEC	61cm/24in 1 Set
	700-0006-03 Pinch shielded ECG lead wire set, 3-lead, IEC	102cm/40in 1 Set
5-Lead ECG Lead Wire Sets, Pinch		
	700-0006-08 Pinch shielded ECG lead wire set, 5-lead, AAMI	61cm/24in 1 Set
	700-0006-10 Pinch shielded ECG lead wire set, 5-lead, AAMI	102cm/40in 1 Set
	700-0006-09 Pinch shielded ECG lead wire set, 5-lead, IEC	61cm/24in 1 Set
	700-0006-11 Pinch shielded ECG lead wire set, 5-lead, IEC	102cm/40in 1 Set
Electrodes		
	369062-001 ECG electrode, transparent tape, adult, single use	Adult 100/Box
	690-0289-00 ECG electrode, transparent tape, adult, single use	Adult 30/Box
	690-0262-00 ECG electrode, transparent tape, general/stress	Adult 30/Box








The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O₂ saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO₂.



Order Code	Description		Qty
Spacelabs TruLink SpO₂ Technology			
SpO₂ Cable			
700-0287-00	SpO ₂ adapter cable	305 cm/10 ft	1 Each
Reusable sensors			
015-0660-00	TruLink SpO ₂ finger sensor	Adult >30 kg/66 lbs	1 Each
Disposable Sensors			
	015-0662-00 TruLink SpO ₂ foam sensor	Adult > 30 kg/66 lbs	24/Box
	015-0663-00 TruLink SpO ₂ vinyl sensor	Adult >30 kg/66 lbs	24/Box
	015-0664-00 TruLink SpO ₂ foam sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0665-00 TruLink SpO ₂ vinyl sensor	Pediatric 10-50 kg/22-110 lbs	24/Box
	015-0666-00 TruLink SpO ₂ finger sensor	Infant 3-20 kg/7-44 lbs	24/Box
	015-0667-00 TruLink SpO ₂ cloth sensor	Neonatal <3 kg/7 lbs	24/Box

The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.



Order Code	Description		Qty
LNCS Masimo SET SpO2 Technology			
LNCS Adapter Cable			
	690-0256-00 SpO2 Adapter Cable, Masimo, LNCS	305 cm/10 ft	1 Each
LNCS reusable sensors			
	690-0247-00 Masimo LNCS DC-I Digit Sensor	Adult > 30 kg	1 Each
	690-0248-00 Masimo LNCS DC-IP Pediatric Digit Sensor	Pediatric 10-50 kg	1 Each
	690-0249-00 Masimo LNCS TC-I Tip-Clip Ear Sensor	Adult > 30 kg	1 Each
	690-0250-00 Masimo LNCS TF-I Transflectance Forehead Sensor	Adult > 30 kg	1 Each
LNCS disposable sensors			
	690-0251-00 Masimo LNCS Amtx adhesive sensor, 18 in.	Adult > 30 kg	20/Box
	690-0252-00 Masimo LCNS Pmtx adhesive sensor, 18 in.	Pediatric 10 - 50 kg	20/Box

The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.



Order Code	Description	Qty
------------	-------------	-----

Nellcor OxiMax SpO2 Technology

Cable

700-0290-00	SpO2 Adapter Cable, Nellcor OxiMax	305cm/10ft	1 Each
690-0271-00	SpO2 Adapter Cable, Nellcor	305cm/10ft	1 Each



Reusable Sensors

690-0003-01	Nellcor OxiMax DS-100A finger sensor	Adult	1 Each
-------------	--------------------------------------	-------	--------

Temperature

Temperature Probe

690-0261-01	Skin/Surface Temp probe, reusable, YSI-400	305cm/10ft	1 Each
-------------	--	------------	--------

Invasive Blood Pressure

Interface Cable



690-0021-00	Compatible with Datex-Ohmeda TC-VTK or Becton Dickinson model 684085 pressure transducers.	366 cm/12 ft	1 Each
-------------	--	--------------	--------

Diaphragm Dome, Single Use



690-0018-00	Diaphragm dome, single use. Rotating luer lock fittings Compatible with Datex-Ohmeda TA101OND		12/Box
-------------	---	--	--------

The élance family of integrated monitors offers you high-performance monitoring of essential parameters including ECG (3- and 5-lead), respiration, noninvasive blood pressure, peripheral O2 saturation and two temperatures. The “i” models add two channels of invasive blood pressure, and the “c” models add microstream end tidal CO2.



Order Code	Description	Qty
------------	-------------	-----

Invasive Blood Pressure

Pressure Transducer Kits, Single use

690-0008-00	Transducer Kit, Patient Mount, DTX/Plus	152.5 cm/60 in	5/Box
690-0009-00	Transducer Kit, Pole Mount, DTX/Plus	152.5 cm/60 in	5/Box
690-0010-00	Transducer Kit, Pole Mount, DTX/Plus	182.9 cm/72 in	5/Box
690-0020-00	Pressure Monitoring Kit, Single Use		1 Each

Pressure Transducer Accessories

690-0012-00	Transducer Mounting Background		1 Each
690-0034-00	Transducer Holder, Reusable		1 Each

Capnography

Sampling Lines

704-0167-00	FilterLine H Set	Adult/Pediatric	25/PK
704-0168-00	FilterLine Set	Adult/Pediatric	25/PK
704-0170-00	Smart CapnoLine	Pediatric	25/PK
704-0171-00	Smart CapnoLine Plus		25/PK

Recording & Printing

Paper



307568-001	TruLink Thermal Paper, Z-fold Pads, 2- or 4-channel printing, blank paper, 300 sheets/pad	45.7 m/150 ft	10 Pads/Box
------------	---	---------------	-------------