

INFU-SURG® CLARITY

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

The latex free Infu-Surg® Clarity Pressure Infusion Bag is an innovative, user-friendly device for intravenous and intra-arterial infusion treatment. Designed with a clear panel, the patient label and barcode on the solution bag can be easily scanned. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand

INFU-SURG® CLARITY PRESSURE INFUSION BAG

- Transparent panel for clear visibility of solution bag
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mmHg)
- Color-coded gauge provides accurate pressure monitoring (0-300 mmHg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Suitable for use with external inflation sources
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested

CODE	DESCRIPTION	UOM
4050	Ethox Clarity Pressure infusion bag 500ml Single Patient use	Box of 5
4050/25	Ethox Clarity Pressure infusion bag 500ml Single Patient use	Box of 25
4100	Ethox Clarity Pressure infusion bag 1000ml Single Patient use	Box of 5
4100/25	Ethox Clarity Pressure infusion bag 1000ml Single Patient use	Box of 25
4030	Ethox Infusurg Pressure infusion bag 3000ml Single Patient use	Box of 5
4030/10	Ethox Infusurg Pressure infusion bag 3000ml Single Patient use	Box of 25