

Infrared Forehead Thermometer



Key Features

- Non contact forehead or object measurement
- Accuracy: +/- 0.2°C
- Measuring range: 32.0°C - 43.0°C
- Response time: approximately 1 second
- Beeps when reading is complete
- Fever alarm
- Stores 10 readings
- Shock resistant
- Automatic shut-off
- Battery: 2 x AAA included
- Battery life: 6000 readings (approx. 1 year)
- Low battery indicator

PRODUCT CODE	DESCRIPTION	UOM
DIA-THE-1001	ECOMED non contact infrared forehead thermometer with 2 x AAA batteries	each

AUSTRALIA

T: 1800 67 5432 | F: 1800 67 5433
E: enquiries@ecomед.com.au | www.ecomed.com.au

NEW ZEALAND

T: 0800 300 421 | F: 0800 300 422
E: sales@ecomед.co.nz | www.ecomed.co.nz