

MS FLEX

MS FLEX is a minor surgery and advanced exploration luminaire. The lamp is manufactured with a LED light source and it is mounted with a highly precise flexible arm. Its arm allows the luminaire to be adjusted in any position.

The quality of the materials used in its manufacture makes it especially suitable for intensive use in surgery and examinations.

Exceptional quality of LED light

The powerful LED light source offers exceptional light quality and a lifespan of more than 50,000 hours.

The high performance of the colour Ra 97% depicts the exceptional quality of white light for an excellent distinction of fabrics and its reproducing all visible colours, naturally and precisely.

Ergonomic design

The circular shape of the head and its handle in the lower part facilitate its exact positioning, fulfilling the ergonomic requirements for any use.

Designed with clean and soft shapes which facilitate hygiene and cleanliness in medical environments.

Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

MS FLEX

| Technical Details | (Code M183500FL) |
|----------------------------|------------------|
| Light source | LED 8W |
| Lux a 50 centimeters | 25.000 lux |
| Adjustable light intensity | No |
| Colour temperature | 4.500 °K |
| Colour rendering | > 97% |
| Lamp lifetime | 50.000 hours |
| Arm length | 70 cm |
| Colour | White |





Wall support
Code B0101800



AH table clamp
Code B010190



Rail Plus bracket
Code 82026



Extension wall arm
Code B010150



Trolley stand
8,8 kg Code 09651

