

# seca 201 & seca 203

## Ergonomic Measuring Tape

new



- Retractable measuring tape for precise body measurements.
- WHR calculator for quick determination of body fat distribution (only 203).
- WHR scale of metal for long life (only 203).

## seca 201 & seca 203: Two high-quality medical products.

In the seca 201 and 203 you've found exactly the right products if you're in search of a measuring tape that is clearly a cut above the rest. Or if you're looking for a measuring tape with millimeter markings that can accurately measure to five millimeters.



### seca 201:

The seca 201 is a true medical product. It's very easy to pull out the tape, which just as easily clicks back into place. Handling is extraordinarily convenient. Made of durable materials, the scale features a tape with a top measurement of 205 centimeters and its own low weight of only 55 grams.


### seca 203:

The seca 203 also has an extra instrument to calculate the Waist-to-Hip Ratio (WHR). This ratio is an indirect method of determining the amount of fatty tissue in the abdomen and consequently, the

distribution of fat reserves. The high-quality scale is made of metal so that you'll be able to depend on this valuable additional function for a long time.



seca 201:  
High-quality materials for precise functionality.



Millimeter graduation for accurate measurements.

seca 203:  
Display to determine body fat distribution.

## seca 201 & seca 203

### Technical Data

#### seca 201:

- Measurement range: 205 cm
- Division: 1 mm
- Diameter: 70 mm
- Height: 22 mm
- Weight of device: about 55 g

#### seca 203:

- Measurement range: 205 cm
- Division: 1mm
- Diameter: 70 mm
- Height: 28 mm
- Weight of device: about 75 g
- Function: WHR calculator



☎ 1800 67 5432  
 📠 1800 67 5433  
 ✉ sales@ecomед.com.au  
 🌐 www.ecomed.com.au