

The scale's high-quality printing is especially tear-resistant. Measurements are precisely displayed and easy to read even after years of use.



For use with infants and toddlers, the measuring rod has head and foot pieces that slide.


seca 207

Baby measuring rod with large callipers

For precise measuring when baby is lying down. The measuring rod can be firmly screwed to the changing table (wall mounting also possible) or used as a mobile rod. Large head and foot positioners make it easy to determine the length and the callipers fold for ease of storage. Made of aluminium for durability.

- Measuring range: 0 – 99 cm
- Graduation: 1 mm

A scale on each side permits measuring of head circumference and facial symmetry.



seca 212

Measuring tape for head circumference of babies and toddlers

Precise and easy-to-read measurements made possible by clear red markings.

This measuring tape of non-stretch Teflon synthetic material measures both the head circumference (front side) and the symmetry of the face (reverse side). A practical dispenser (containing 15 tapes) ensures that they are stored safely.

- Measuring range: 5 – 59 cm
- Graduation: 1 mm

ECOMED

Serving the HealthCare Industry

T 1800 67 5432

F 1800 67 5433

E enquiries@ecomед.com.au

W ecomед.com.au