

# seca 218

Disposable measuring tape  
in practical wall dispenser

new



- Disposable measuring tape of high-quality paper to prevent cross infection.
- Uncomplicated removal of measuring tape for immediate use.
- At hand at all times from exclusive wall dispenser.
- Extra space to write down name, height, weight and date.

## seca 218:

### The measuring tape for daily use.

Our experts are constantly striving to make your routine work as easy as possible. That's why seca now offers the environmentally friendly and high-quality disposable measuring tape seca 218 in a practical wall dispenser for improved hygienic conditions. Now you have the desired measuring tape at hand when you need it. The small, space-saving wall dispenser is easy to install and to refill.



#### Always at hand in a wall dispenser.

The seca 218 is available when needed. One hundred measuring tapes fit into a small wall dispenser that releases them one at a time. The dispenser can be mounted easily anywhere and refilled quickly and simply.




Simple refilling of up to 100 measuring tapes.

#### Optimum handling.

Thanks to the carefully chosen weight of the paper, the measuring tape is strong enough to do the job and light enough to tear easily from the wall dispenser. It's an environmentally friendly solution that can be recycled simply.

#### Hygienic one-time use.

The disposable measuring tape seca 218 fulfills demanding medical hygiene standards. One-time use of the stable measuring tapes to prevent cross infection.



Extra space to write down name, height, weight and date.

## seca 218

### Technical Data

- Measurement range: 0–100 cm / 0–36"
- Graduation: 1 mm / 1/8"
- Weight of device: 4 g / 0.14 oz
- Dimensions of wall dispenser (WxHxD): 50 x 73 x 36 mm / 1.97 x 2.87 x 1.42"
- Quantity: 500 / 1.000



1800 67 5432  
 1800 67 5433  
[sales@ecomед.com.au](mailto:sales@ecomед.com.au)  
[www.ecomed.com.au](http://www.ecomed.com.au)