



Low weight and practical handle make transport effortless.

seca 334

Mobile electronic baby scales

The weighing tray has a more generous and comfortable design than many other models, but the easy-to-transport seca 334 weighs only 3.5 kg and requires surprisingly little space. With the large, clear LCD figures and the practical carrying handle – which can also be used to hang the scales – this battery-operated scale for babies and infants makes daily weighing very easy. The breast-milk-intake function shows the exact amount of milk the baby has consumed.

- Capacity: 20 kg
- Graduation: 5 g < 10 kg > 10 g
- Functions: Automatic switch-off, TARE, Auto-HOLD, BMIF, kg/lbs switch-over
- Optional: Measuring rod seca 232, head positioner seca 419, carrying case seca 428



Accessories



seca 232

Measuring rod for baby scales 334

Just as easy to fit as it is to read: this baby measuring rod with integrated head and foot positioners makes the seca 334 even more efficient because now weighing and measuring can be done in one step.



The optional measuring rod permits weighing and measuring in one step.

- Measuring range: 35 – 80 cm
- Graduation: 1 mm

Accessories



seca 414

Carrying case for seca baby scales

The seca 334 baby scales can be transported safely and compactly in this practical carrying case of water-repellent material.

- With shoulder strap and tailor-made pads

The baby scales can be easily stored and transported.



Serving The HealthCare Industry

1800 67 5432

1800 67 5433

enquiries@ecomед.com.au

www.ecomed.com.au