



The weighing tray can be removed quickly and easily at the touch of a button.



seca 354

Electronic baby scales double as flat scales for children

With its dual function, this baby scale is a small sensation: at the touch of a button, the tray can be detached quickly and easily from the base, leaving you with robust flat scales for children weighing up to 20 kg. As the two elements are securely connected, it is absolutely safe to lift and transport the scales. The softly rounded shape of the tray gives the baby a feeling of safety and security. The breast-milk-intake function shows the exact amount of milk the baby has consumed.

- Capacity: 20 kg
- Graduation: 10 g < 10 kg > 20 g
- Functions: TARE, Auto-HOLD, BMIF, kg/lbs switch-over
- Optional: Carrying case seca 413

Accessories

seca 413

Carrying case for seca 354 baby scales

With this case, seca 354 baby scales can be transported comfortably and safely. The exact fit of the seca 413 provides optimum protection and the large opening facilitates the packing and unpacking of the scales. An adjustable shoulder strap makes transport easier and the water-repellent nylon surface keeps everything dry.



- Weight: 500 g

ECOMED

Serving the HealthCare Industry
T 1800 67 5432
F 1800 67 5433
E enquiries@ecomed.com.au
W ecomed.com.au