



The unique V-shape automatically places the baby in the right measuring position.



seca 416

Infantometer for measuring babies and infants

The completely novel shape of the infantometer makes it quick and easy to measure babies and infants up to two years old. The board is generously designed and the raised and softly rounded sides hold the baby easily and safely in the right position. The footpiece, fitted into two guide rails, slides smoothly along the scale and, as it is lockable, the baby can be removed and the result noted later. All parts are tough and durable.



- Measuring range: 33 – 100 cm
- Graduation: 1 mm

The head and foot positioners lock into place so the indication can be read after measuring.

ECOMED

Serving the HealthCare Industry

T 1800 67 5432

F 1800 67 5433

E enquiries@ecomед.com.au

W ecomед.com.au