

Petite Elite Portable Suction Pump



Key Features

- Suitable for use anywhere in or out of the hospital
- High performance >25Lpm at 500+ mmHg
- Lightweight and portable with optional carry bag
- Single x 1000ml standard jar assembly, autoclaveable
- Robust moulded case construction
- Encapsulated bacterial/hydrophobic
- Filter disk
- Luminous touch control panel
- Maximum vacuum override button
- Battery level/charge indicator
- 115v - 240V AC
- Battery duration up to 1hr 20m
- Battery 12v dc, 4.3 Mah, rechargeable
- 4.59kg

ITEM NO.	DESCRIPTION
PSP001	Petite Elite portable suction pump with reusable autoclavable canister
PSP004	Petite Elite portable suction pump with disposable VacSax canister
PPSP050	Petite Elite carry bag
PSP500	1000ml standard jar assembly
PSP600	Filter disks

AUSTRALIA

T: 1800 67 5432 | F: 1800 67 5433
E: enquiries@ecomed.com.au | www.ecomed.com.au

NEW ZEALAND

T: 0800 300 421 | F: 0800 300 422
E: sales@ecomed.co.nz | www.ecomed.co.nz