



Serving The HealthCare Industry

Digital Thermometer



DIA-THE-1003



DIA-THE-1002

Key Features

- Oral, rectal or under-arm measurement
- Large easy to read screen
- Accuracy: +/- 0.1°C
- Plastic storage case included
- Beeper function
- Memory: last measuring reading
- Automatic power off
- Battery: one 1.5V(LR/SR41) included
- Battery life: 100 hours under continuous operation

PRODUCT CODE	DESCRIPTION	UOM
DIA-THE-1002	ECOMED digital medical thermometer with 1.5V battery and plastic storage case	each
DIA-THE-1003	Disposable covers, digital thermometer	box of 100

AUSTRALIA

T: 1800 67 5432 | F: 1800 67 5433
E: enquiries@ecomed.com.au | www.ecomed.com.au

NEW ZEALAND

T: 0800 300 421 | F: 0800 300 422
E: sales@ecomed.co.nz | www.ecomed.co.nz