

Rapid Digital Thermometer



DIA-THE-1003



DIA-THE-1029

TGA ARTG Number 354821

Key Features

- Measuring range: 32.0°C - 42.9°C
- Measuring accuracy: $\pm 0.1^\circ\text{C}$
- Response: 10 seconds
- Water-proof
- Flexible tip
- Beeper function
- Memory: Last measuring reading
- Automatic power off: within 10 Mins on after non-use
- Battery: One 1.5V (LR/SR41) button battery
- Fever alarm for temp greater than 37.8 °C

PRODUCT CODE	DESCRIPTION	UOM
DIA-THE-1029	Rapid Digital Thermometer	Each
DIA-THE-1003	Disposable covers, digital thermometer	Box of 100