



Reusable NIBP Cuffs



Specifications:

- Non Invasive Blood Pressure cuffs for patient monitors
- Reusable for multiple patient use
- Sized to meet AHA,AAMI / ANSI/SP10 and EN1060 requirements
- Fitting different brands of monitors if proper connector attached
- Latex free
- CE marked

Warranty:


6 months warranty after delivery unless special requirements








Colorful PU Leather BP Cuffs



- ◆ Reusable for multiple patient use.
- ◆ PU leather jacket with different colors to clearly signify the patient size applied.
- ◆ Individual TPU bladder with 50cm single/dual tubes.
- ◆ Convenient and easy to clean.
- ◆ Transparent PE bag package

Cuff without Connector	Air Tube	Patient size	Color	Size Range	Part Number	Package
 Connectors should be ordered separately if necessary	Single-tube 	Adult		25-35cm	LF1103S	1pc/bag
		Pediatric		18-26cm	LF1102S	1pc/bag
		Infant		10-19cm	LF1101S	1pc/bag
		Neonate		6-11cm	LF1100S	1pc/bag
		Large Adult		33-47cm	LF1104S	1pc/bag
		Thigh		46-66cm	LF1105S	1pc/bag
		Adult long		25-35/45cm	LF1107S	1pc/bag
	Double-tube 	L- adult long		33-47/55cm	LF1108S	1pc/bag
		Adult		25-35cm	LF1103D	1pc/bag
		Pediatric		18-26cm	LF1102D	1pc/bag
		Infant		10-19cm	LF1101D	1pc/bag
		Neonate		6-11cm	LF1100D	1pc/bag
		Large Adult		33-47cm	LF1104D	1pc/bag
		Thigh		46-66cm	LF1105D	1pc/bag
		Adult long		25-35/45cm	LF1107D	1pc/bag
		L- adult long		33-47/55cm	LF1108D	1pc/bag

Part Number	Patient type	Color	Size Range	Package	Cuff with Connector
LF1103S-12	Adult		25-35cm	1pc/bag	Connector : LC-B12 
LF1102S-12	Pediatric		18-26cm	1pc/bag	
LF1101S-12	Infant		10-19cm	1pc/bag	
LF1100S-12	Neonate		6-11cm	1pc/bag	
LF1104S-12	Large Adult		33-47cm	1pc/bag	
LF1105S-12	Thigh		46-66cm	1pc/bag	
LF1107S-12	Adult long		25-35/45cm	1pc/bag	
LF1108S-12	Large adult long		33-47/55cm	1pc/bag	
LF1103S-03	Adult		25-35cm	1pc/bag	Connector : LC-B03 
LF1102S-03	Pediatric		18-26cm	1pc/bag	
LF1101S-03	Infant		10-19cm	1pc/bag	
LF1100S-03	Neonate		6-11cm	1pc/bag	
LF1104S-03	Large Adult		33-47cm	1pc/bag	
LF1105S-03	Thigh		46-66cm	1pc/bag	
LF1107S-03	Adult long		25-35/45cm	1pc/bag	
LF1108S-03	Large adult long		33-47/55cm	1pc/bag	
LF1103S-18	Adult		25-35cm	1pc/bag	Single-tube LC-B18 
LF1102S-18	Pediatric		18-26cm	1pc/bag	
LF1101S-18	Infant		10-19cm	1pc/bag	
LF1100S-18	Neonate		6-11cm	1pc/bag	
LF1104S-18	Large Adult		33-47cm	1pc/bag	
LF1105S-18	Thigh		46-66cm	1pc/bag	
LF1107S-18	Adult long		25-35/45cm	1pc/bag	
LF1108S-18	Large adult long		33-47/55cm	1pc/bag	
LF1103D-17/18	Adult		25-35cm	1pc/bag	Double-tube LC-B17+18 
LF1102D-17/18	Pediatric		18-26cm	1pc/bag	
LF1101D-17/18	Infant		10-19cm	1pc/bag	
LF1100D-17/18	Neonate		6-11cm	1pc/bag	
LF1104D-17/18	Large Adult		33-47cm	1pc/bag	
LF1105D-17/18	Thigh		46-66cm	1pc/bag	
LF1107D-17/18	Adult long		25-35/45cm	1pc/bag	
LF1108D-17/18	Large adult long		33-47/55cm	1pc/bag	

Part Number	Patient type	Color	Size Range	Package	Cuff with Connector
LF1103D-08	Adult		25-35cm	1pc/bag	Double-tube LC-B08+08 
LF1102D-08	Pediatric		18-26cm	1pc/bag	
LF1101D-08	Infant		10-19cm	1pc/bag	
LF1100D-08	Neonate		6-11cm	1pc/bag	
LF1104D-08	Large Adult		33-47cm	1pc/bag	
LF1105D-08	Thigh		46-66cm	1pc/bag	
LF1107D-08	Adult long		25-35/45cm	1pc/bag	
LF1108D-08	Large adult long		33-47/55cm	1pc/bag	
LF1103D-48	Adult		25-35cm	1pc/bag	Double-tube LC-B48 
LF1102D-48	Pediatric		18-26cm	1pc/bag	
LF1101D-48	Infant		10-19cm	1pc/bag	
LF1100D-48	Neonate		6-11cm	1pc/bag	
LF1104D-48	Large Adult		33-47cm	1pc/bag	
LF1105D-48	Thigh		46-66cm	1pc/bag	
LF1107D-48	Adult long		25-35/45cm	1pc/bag	
LF1108D-48	Large adult long		33-47/55cm	1pc/bag	

Connectors:

Part Number	LC-B12-P	LC-B17	LC-B17	LC-B18	LC-B08
Image					
Mate fitting	Fitting LC-B15	Fitting LC-B15	Fitting LC-B18	Fitting LC-B17	Fitting LC-B09

Part Number	LC-B03	LC-B2T-P	LC- B2T	LC-B05-P	LC-B32
Image					
Mate fitting	Fitting LC-B06	Fitting air tubes	Fitting air tubes	Fitting LC-B20	Fitting LC-B31

Part Number	LC-B54	LC-B55	LC-B48	LC-B..	LC-B..
Image					
Mate fitting	Fitting LC-B55	Fitting LC-B54	Fitting LC-B47		