

M70C FINGERTIP PULSE OXIMETER



OLED screen with 6 display modes



4 direction display



Resistance to surrounding light and man-made light



Auto shut off without signal after 8s



5 levels adjustable brightness



Bluetooth (Optional)



SPECIFICATIONS:

Weight: about 21g (exclude batteries)

Dimensions: 60mm(length) × 35mm(width) × 35mm (height)

Power Requirements: Battery: AAx2

Parameters: SpO2, PR (Pulse Rate)

SpO2

Measuring Range: 0~100%

Resolution: 1%

Accuracy: At 70%~100%, ±2%; At 0~69%, unspecified

PR

Measuring Range: 25 bpm ~250 bpm

Resolution: 1 bpm

Accuracy: ±1% or ± 1 bpm, whichever is greater

