

M70D

FINGERTIP PULSE OXIMETER



LCD screen with
backlight



Resistance to surrounding
light and man-made light



Auto shut off without
signal after 8s



Streamline design

SPECIFICATIONS:

Weight: about 21g (exclude batteries)

Dimensions: 60mm(length) × 37mm(width) × 31mm (height)

Power Requirements: Battery: AAAx2

Parameters: SpO2, PR(Pulse Rate)

SpO2

Measuring Range: 0~100%

Resolution: 1%

Accuracy: At 70%~100%, ±2%; At 0~69%, unspecified

PR

Measuring Range: 25 bpm ~250 bpm

Resolution: 1 bpm

Accuracy: ±1% or ± 1 bpm, whichever is greater