

INTRODUCING THE RIESTER USB-C HIGH PERFORMANCE BATTERY

USB-C Handle

- ▶ High capacity, 10 hour operating time
- ▶ Battery status indicator visible on handle (bicolor LED indicates fully charged, charging, or battery needs to be charged)
- ▶ Wide range charger (100 - 240 V)
- ▶ Universal USB type C charging technology
- ▶ Fully charged in 3.5 hours



ri-derma[®] OTOSCOPES



ri-scope[®] L DERMATOSCOPES



ri-scope[®] OPHTHALMOSCOPES



Type C, USB-C

Rechargeable Battery Handle with rheotronic®

The Riester Type C, USB-C rechargeable battery handle with rheotronic® is quickly and easily charged by convenient USB-C charging cable in 3.5 hours. This solid metal handle is made from solid black metal and equipped with a practical bayonet attachment mechanism for quickly interchanging instrument heads.

- › Made from solid metal
- › Practical bayonet attachment method for quickly replacing instrument heads
- › Handle diameter: 28 mm
- › Easy to charge via USB-C cable provided
- › Usable for all **ri-scope®** and **ri-derma®** 2.5 V and 3.5 V instrument heads
- › Charger and rechargeable battery included
- › Battery status LEDs
- › Variety of country specific plugs
(We offer EU, UK, US, Australia as standard, however on request - 10 further plugs are available)
- › Rechargeable battery: 18650 Li-ion Battery, 3,7 V 2600 mAh 9,62 Wh
- › Charging time: Approx. 3.5 h
- › Working time: Approx. 10 h with LED illumination, 3.5 h with Xenon illumination

Status LEDs:

- › Red: Battery needs to be recharged
- › Green flashing: Charging
- › Green constant: Battery fully charged
- › Cable length: 1m

rheotronic®

The patented rheotronic® technology for LEDs and light bulbs allows for continuous light intensity adjustment directly on the handle.

- › Turn on by turning the knob to +
- › Shut down by turning the knob to -
- › Continuous adjustment by turning toward + or -
- › Automatic safety shutdown after 120 seconds
- › Variety of country specific plugs
(We offer EU, UK, US, Australia as standard, however on request - 10 further plugs are available by supplier).

| | |
|--------------------------|--|
| Art.No. 10704.003 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, US |
| Art.No. 10704.002 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, UK |
| Art.No. 10704.001 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, AU |
| Art.No. 10704 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, EU |

