

Convenience and Security Options



Wallmount 134201
(TAT-5000 not included)



Quick Release Security System:
124277 (Core) 124278 (Oral Equiv)
Includes TAT-5000, Wallmount
134201. Latex-free coiled cable



Quick Release Security System:
124286 (Core) 124287 (Oral Equiv)
Includes TAT-5000, Wallmount
134201, 8' Latex free nylon
covered steel cable



Quick Release Security System:
124292 (Core) 124293 (Oral Equiv)
Includes TAT-5000, Wallmount
134306, Latex free coiled cable



Quick Release Security System:
124288 (Core) 124289 (Oral Equiv)
Includes TAT-5000, Wallmount
134305, Latex free coiled cable



Quick Release Security System:
124290 (Core) 124291 (Oral Equiv)
Includes TAT-5000, Wallmount
134305, 8' Latex free nylon
covered steel cable



Dual Security System
5 or 20 Temps before Lockout
5: 124280 (Core) - 124281 (Oral Equiv)
20: 124392 (Core) - 124393 (Oral Equiv)
Includes: TAT-5000, Requires
Keyless Locking Wallmount 134308



Dual Security System
5 or 20 Temps before Lockout
5: 124280 (Core) - 124281 (Oral Equiv)
20: 124392 (Core) - 124393 (Oral Equiv)
Includes: TAT-5000, Requires
Keyless Locking Wallmount 134307



Quick Release Security System:
124294 (Core) 124295 (Oral Equiv)
Includes TAT-5000, Wallmount
134306, 8' Latex free nylon



8' Security Cable
Part #124311 - Latex Free Coiled Cable



Nylon Covered Steel Cable
Latex Free
6' Cable - Part #134302
8' Cable - Part #134030

- **Quick Release System:** A discrete method for easily removal and reattachment of the cable when required.
- **Dual Security System:** A system requiring return of the instrument to the holder after 5 or 20 temperature measurements.
- **Suggestions for Selection:**
 - If instruments are to be secured at the nurses' station, consider the keyless locking wallmount or Dual Security System.
 - If instruments are to be mounted one per bed, or in a dedicated area (e.g. triage), consider Quick Release System.
 - If instruments are to be attached to an existing rolling stand, or a vital signs monitor, consider the 6' vinyl covered steel cable 134302.