

# CryOmega®

MULTI-USE DISPOSABLE  
CRYOSURGICAL DEVICE



## PERFECTLY PRECISE CRYOSURGERY



COMMON  
WART BEFORE  
TREATMENT



IMMEDIATELY  
AFTER  
TREATMENT



TWO WEEKS  
AFTER  
TREATMENT



EASIER CARTRIDGE  
ACTIVATION SYSTEM

PINPOINT DELIVERY  
OF CRYOGEN



- CryOmega® dispenses liquid nitrous oxide (-89°C) the coldest portable cryogen and it is cleared for 16 common lesions including verrucae, skin tags, and seborrheic keratoses.
- CryOmega® is an ideal choice in portable cryosurgery for physician offices looking for convenience and great first-time outcomes.
- Portability and easy storage are “must-haves” when it comes to cryosurgery in your office.

Clinically **effective** direct application of liquid **nitrous oxide** (-89°C) in a **convenient** and **economical** device.

CryOmega® is supported by an extensive library of how-to videos which can be found on our website

Item No.	Product Code	CryOmega® Product Description
160-2001	CryOmega® Single Pack INTL	Contains: 1 Standard Pen with 16g of N <sub>2</sub> O cryogen 1 Instructions for Use
160-2004	CryOmega® Twin Pack INTL	Contains: 2 Standard Pens with 16g of N <sub>2</sub> O cryogen 1 Instructions for Use

 **CryoConcepts**  
THE CRYO EXPERTS

**ECOMED**  
Serving The HealthCare Industry