

Histofreezer[®] **STAT**

ADVANCED
PORTABLE CRYOSURGERY

**Designed with One
Goal in Mind—Better
First-time Results for
Your Patients.**



CONES



BUDS

NEW!

BENEFITS:

- ✓ **Better Design**
- ✓ **Better Outcomes**
- ✓ **Better for the Environment**



“ **There is so much to like about
the new Histofreezer[®] STAT!!** ”

~ Dr. David White, Pediatrician, Nashville, TN

Histofreezer[®] STAT

ADVANCED
PORTABLE CRYOSURGERY

Breakthrough Features:



Uses CONES
OR BUDS



Effective
Freezing



High
Efficacy



Immediately
Cold




Single-Handed
Design




Higher
Safety



50 or 80 Count Mixed Kit




Histofreezer[®] STAT
closed cones deliver
precise freezing of
the lesion.



Histofreezer[®] STAT buds are ideal for
hard-to-treat locations and can be used
on or off the device.

HISTOFREEZER[®] STAT CLOSED CONE DESIGN FREEZING TEMPERATURE OVER TIME



— Histofreezer[®] FLEX

Closed cones deliver and maintain a critical
freezing temperature (-50C).

Item No.	Kit Product Code	Histofreezer [®] STAT Product Description
ALL KITS CONTAIN: One Release Clip, One Cone Selection Guide, One Practice Pad, and One Instructions for Use		
1500-1008	50SM Mixed Kit	Kit Contains: 50 Treatments 1 Canister of cryogen 30 Buds 40 Cones (3, 5, 7, 9, and 12mm) 4 Cone Adapters
1500-1009	80SM Mixed Kit	Kit Contains: 80 Treatments 1 Canister of cryogen 50 Buds 40 Cones (3, 5, 7, 9, and 12mm) 4 Cone Adapters

 **CryoConcepts**
THE CRYO EXPERTS

ECOMED
Serving The HealthCare Industry