


PE-1204

Bluetooth ECG




The tablet not included


• Product Features

- Simultaneous 12 leads ECG acquisition and display
- Lightweight and ultra-compact design
- Sampling box supports wired and wireless ECG acquisition
- Compatible with Windows and Android operating system
- Automatic ECG measurements and interpretation
- Data access with HIS/ EMR/ PACS based on DICOM/ FTP standard protocol
- Multi- format reports export, such as PDF, JPEG, BMP, EM-XML, SCP*, DICOM, FDA-XML*
- The comprehensive patient case management supports quick retrieve, query or combined query, delete, edit, upload and export
- Patient information can be retrieved from the worklist*

- Software Interface




12 Leads Mode




Lead Off Alarm

- Flexible and Convenient Data Transmission Solution



- Multiple Report Formats with Automatic ECG Measurements and Interpretation



3CHx4+3R

6CHx2+1R

Measurement Parameter Report

3CHx4+1R

Rhythm Lead