

freezpen®

cryosurgery
made simple



Freezpen is a powerful, precise cryosurgery instrument for the removal of benign skin lesions.

Simply powerful

High power (50-55 bar) microjet achieves in-tissue temperature of $-89^{\circ}\text{C}/-128^{\circ}\text{F}$, for complete tissue necrosis in seconds.

Simply precise

Ultrafine sharp jet for total precision and no damage to healthy tissue.
Range of applicators.

Simple procedure

Assess the lesion, choose the right applicator, point, freeze, wipe, freeze again. Done!

Simply painless

No pain, no cuts or bleeding, no anesthesia - even for children.
Slight sting signals when treatment is complete.

Simply perfect

Perfect results with no scarring. Healthy tissue is unaffected; frozen tissue disappears in a few weeks.

Simply practical

Easy handling and transport. Handy disposable cartridges; store at room temperature.

It's simply good business!

Each treatment takes just seconds. Return on investment in 10 weeks at just six treatments/week.



Freezpen can remove most common skin pathologies in just a few seconds.



Pigmented Spots



Hemangioma



Warts



Fibroma



Acrochordon




Skin Tags



Pigmented Spots



Hemangioma



Keratoses



Solar Lentigo



Plantar Wart



Age Spots

Before treatment



After treatment





Freezen 16g instrument set - 3 applicators:
Blue, white and green applicators.
User manual and tweezers for filter placement.



Box with 12 N₂O cartridges of 16g
each carton with 12 boxes.



Freezen 35g instrument set - 3 applicators:
Green, orange long and white long applicators.
User manual and tweezers for filter placement.



Box with 6 N₂O cartridges of 35g
each carton with 12 boxes.



Applicators

- 1-3mm applicator |
- 2-6mm applicator ||
- 4-8mm applicator |||

power

Long applicators - only for the Freezen 35g

- 1-4mm long applicator |
- 2-6mm long applicator ||

power