


Disposable BP cuffs Series NO.202

Neonate BP cuffs for single patient use.
Soft fiber with non-woven fabrics.
5 sizes to fit different patient range.
20cm single/ dual tubes with connector.
Packed in disposable paper-plastic bag.
Out of pack warranty, cuff can be replaced if found to be faulty when used for the first time.



Size Number	Part Number	Limb circumference	Cuffs Size (length*width)	Unit Price	Tube& Connector
Size #1	LF1411S-05	3-6cm	11*3.2cm		
Size #2	LF1412S-05	4-8cm	14*3.6cm		
Size #3	LF1413S-05	6-11cm	20*4.8cm		
Size #4	LF1414S-05	7-13cm	24.8*5.6cm		
Size #5	LF1415S-05	8-15cm	27*6.3cm		
Size #1	LF1411S-12	3-6cm	11*3.2cm		
Size #2	LF1412S-12	4-8cm	14*3.6cm		
Size #3	LF1413S-12	6-11cm	20*4.8cm		
Size #4	LF1414S-12	7-13cm	24.8*5.6cm		
Size #5	LF1415S-12	8-15cm	27*6.3cm		
Size #1	LF1411S-55	3-6cm	11*3.2cm		
Size #2	LF1412S-55	4-8cm	14*3.6cm		
Size #3	LF1413S-55	6-11cm	20*4.8cm		
Size #4	LF1414S-55	7-13cm	24.8*5.6cm		
Size #5	LF1415S-55	8-15cm	27*6.3cm		
Size #1	LF1411D-05	3-6cm	11*3.2cm		
Size #2	LF1412D-05	4-8cm	14*3.6cm		
Size #3	LF1413D-05	6-11cm	20*4.8cm		
Size #4	LF1414D-05	7-13cm	24.8*5.6cm		
Size #5	LF1415D-05	8-15cm	27*6.3cm		