


# Measuring tape for head circumference of babies and toddlers


## seca 212

- + Hygienic, disposable measuring tape of non-stretch Teflon material
- + Double function: head circumference with the front side of the tape and facial symmetry with the rear side
- + Easy-to-read, colour-coded precise results
- + Practical dispenser (content: 15 tapes)



seca 212

Measuring tape for head circumference of babies and toddlers



Measuring tape for head circumference of babies and toddlers

The measuring tape is made of non-stretch Teflon synthetic material and measures both the head circumference (front side) and the symmetry of the face (reverse side). A practical dispenser (containing 15 tapes) ensures that they are stored safely.

Measuring range	3 – 59 cm / 1 – 23 in
Graduation	1 mm / 1/16 in
Dimensions (WxHxD)	25 × 0.2 × 663 mm / 1 × 0.01 × 26.1 in
Net weight	0.1 kg / 0.2 lbs
Medical Device Regulation (MDR) certified	