

AUDIOMETER



STEP INTO A WIRELESS WORLD

AUDIOLYSER® ADL CONNECT

HEARING SCREENING AND MONITORING

CLASS 4 TONE AUDIOMETER

- ▶ Wired and Bluetooth connectivity
- ▶ High-reliability answer button
- ▶ Time-saving: up to 4 exams in parallel
- ▶ Multifunctional docking station

Model: FF1166


CE
0459






Holmco
PD81
Headphones

SECURITY AND INTEROPERABILITY IN MIND

Audiolyser® ADL Connect operates **without a local database**, enhancing data security. It offers **extensive interoperability** with GUEST (old and new systems), iodata and GDT formats, making it easy to integrate into existing systems. **User two-factor authentication** and **integration with electronic medical records**  will enhance the device starting in Q4 2025. This solution is designed to meet the growing demands for cybersecurity and interoperability in the medical environment.


MULTILINGUAL VOICE ASSISTANT



The audiometer is equipped with a Multilingual Voice Assistant that allows you to **communicate test instructions** directly to the patient **in the language of your choice**. Tests can be performed **automatically, saving the operator time**.

More than 25 languages included: Arabic, German, English, Chinese, Spanish, French, Italian, Dutch, Norwegian, Polish, Portuguese, Thai, among others.

AUDIOWIN® SOFTWARE



THE +

- ▶ Ergonomic interface with simplified navigation.
- ▶ Integrated start-up wizard and self-training module.
- ▶ Automatic software updates.
- ▶ Customisable PDF reports
- ▶ Backward compatibility with previous systems (Audiolyser® ADL20).
- ▶ Automated daily check with logging for greater traceability. (Q4 2025)
- ▶ Ambient noise monitoring during examination (Q4 2025)

Sigycop
report
example



▶ Numerous index calculations:

Asymmetric Hearing Loss, Fletcher (weighted), Fletcher, EWI, HSE, MPB2002, MERLUZZI 1979, Mean Hearing Loss, SNCF criteria, SIGYCOP: evaluation of auditory thresholds under the letter O, used to determine an individual's medical suitability for service in the French armed forces or fire and rescue services (SDIS).

Software supplied with the device (1 license per work station)



83% of French people
bothered by noise are at
their place of work

(Ifop survey – JNA 2022 Noise and health)

Hearing loss is a type of deafness recognised as an **occupational disease**. It is included in table n°42 of occupational diseases of Decree n°63-405 of 10/04/1963, updated 25/09/2003. It is diagnosed by an audiometry showing, on the best ear, **a deficit of at least 35dB**, performed after cessation of noise exposure for at least 3 consecutive days. (INRS).

YOUR DISTRIBUTOR



1 FREE BOX
OF EARPAD
COVERS*

ACCESSORIES

(included with your device)

- 1 Box of biocompatible, single-use earpad covers (preserves earpad life time, is more hygienic for your earpad, and does not muffle the sound)
- 1 Travel bag

2
YEARS
WARRANTY
ADL CONNECT



CE
0459