

## EASY-TO-USE CARTRIDGES

One of the unique features of Freezpen is its handy N<sub>2</sub>O cartridges. Cartridges can be safely stored at room temperature until needed, and loading a new cartridge takes just seconds.



## PACKAGING

**HS code**  
28112930

### 8g N<sub>2</sub>O Cartridges for Freezpen 8g

#### Box

Contains 21 cartridges

Dimensions:  
11 x 17 x 13 cm

Weight: 785g

Product reference:  
CL-N2O-8G-NBOX

#### Carton

Contains 16 boxes or 336 cartridges

Dimensions:  
34 x 44 x 26.5 cm

Weight: 12,900g

Product reference:  
CL-N2O-8G-NCARTON

### 16g N<sub>2</sub>O Cartridges for Freezpen 16g

Contains 12 cartridges

Dimensions:  
11 x 12.5 x 14.5 cm

Weight: 848g

Product reference:  
CL-N2O-16G-NBOX

Contains 12 boxes or 144 cartridges

Dimensions:  
32 x 26.5 x 30 cm

Weight: 10,400g

Product reference:  
CL-N2O-16G-NCARTON

### 35g N<sub>2</sub>O Cartridges for Freezpen 35g

Contains 12 cartridges

Dimensions:  
7.5 x 14.3 x 17.5 cm

Weight: 925g

Product reference:  
CL-N2O-35G-BOX

Contains 12 boxes or 72 cartridges

Dimensions:  
44 x 30.5 x 20.5 cm

Weight: 11,500g

Product reference:  
CL-N2O-35G-CARTON

## PRODUCT INFORMATION

### Storage

Once punctured, a cartridge maintains full pressure for several weeks. Store at room temperature. Cartridges have a shelf life of 36 months.

### Cartridge pressure

Depending on the surrounding temperature, the cartridge pressure is 49-53 bars / 710-770 psi.

### N<sub>2</sub>O liquid temperature

The temperature of the dispensing liquid is -89°C.

### Disposal

Used cartridges can be disposed of as metal waste, in the same way as metal cans.