

eVAC

Hand Operated Suction Pump

The eVAC is a simple and convenient hand operated *emergency* suction pump, offering an ergonomic design and ease of operation. With its robust, lightweight construction you are assured the eVAC will perform reliably, generating high vacuum when needed most.

The eVAC comes complete with two catheters suitable for adult and child pharyngeal suctioning and a single use canister that releases effortlessly from the pump for hygienic disposal.

- High Vacuum
- Overflow protection
- Lightweight
- Simple to operate
- Maintenance free
- Easy to store
- Compact design
- Disposable jars and catheters



In emergencies insist on **eVAC**

eVAC Hand Operated Suction Pump

TECHNICAL SPECIFICATIONS FOR eVAC HAND OPERATED SUCTION PUMP	
Model Catalogue Number	eVAC SUC 82015
Vacuum	0 to -73 kPa (0 to -550 mmHg)
Flow Rate	Maximum at least 20 litre/min
Disposable container volume	300ml
External diameter of cap	17mm
External diameter of catheter connection	13mm and 6mm
Operating Environmental Temperature	-20 °C to +50 °C
Storage Environmental Temperature	-40 °C to +60 °C
Unit Dimensions	198 x 169 x 74mm
Shipping Dimensions	230 x 180 x 85mm
Unit Weight	247g
Shipping Weight	359 g

SPARES & ACCESSORIES	
SUC 82015001	Disposable collection jar
SUC 82015002	Disposable suction catheter kit (includes Adult and Child)

All specifications are subject to change without notice. Please check the exact specification when ordering.

The Clements range of products are renowned worldwide for quality, reliability and durability.



ECOMED
Serving The HealthCare Industry

t: 1800 675 432
f: 1800 675 433
e: enquiries@ecomед.com.au
w: www.ecomed.com.au

