

Breathe *Eazy* Spacer

ANTI-STATIC CHAMBER SPACER FOR METERED DOSE INHALER WITH VISION VALVE

- ☆ For the efficient delivery of asthma and COPD medication to the lungs.
- ☆ Easy to use with all standard 'puffers' & metered dose inhalers
- ☆ Conveniently store puffer inside spacer
- ☆ Compact, portable, durable
- ☆ Latex Free
- ☆ Easy to Clean
- ☆ ***Baby, Child, and Adult Mask available separately.***



Breathe Eazy Spacer
Breathe Eazy
Adult Mask



Breathe Eazy
Child Mask






Breathe Eazy
Baby Mask

Breathe *Eazy* Spacer

ANTI-STATIC CHAMBER SPACER FOR METERED DOSE INHALER WITH VISION VALVE



Part Number	Description	Included in Package
ANR-SPA-0006	Breathe Eazy Spacer with Mouthpiece. Age 6+ years	
ANR-SPA-0012	Breathe Eazy Baby Mask. Age 0-3 years	
ANR-SPA-0009	Breathe Eazy Child Mask. Age 3-6 years	
ANR-SPA-0014	Breathe Eazy Adult Mask. Age 6+ years	