

Medical Regulators Model 801

The Model 801 Single Stage Preset Regulator possesses the sensitivity and flow capacity of much larger regulators in a compact lightweight design. Preset at 400kPa, the Model 801 is engineered to provide equipment downstream with constant pressure.

- Features include - Wt. 1.3kg
- ▶ Proven seat design
 - ▶ Fully chromed
 - ▶ Large, easy to read gauges
 - ▶ Extensive monel filter
 - ▶ Complies to AS 2896 & AS 4267

801Z258 Medical Nit/Ox with Yoke

801Z259 Medical Air with Yoke

801Z260 Medical Oxy with Yoke

Medical Yokes

Pin Index Yoke, constructed of a high quality brass and chrome plating, complete with extensive filter and neoprene/metal seal.

Y01870 Oxy Yoke

Y01910 Nitrous Oxide Yoke

Y01940 CO2 Yoke

Y01950 Air Yoke



Medical Flowmeters

Harris Flowmeters feature an integrated one-piece design with well defined calibrations which accurately measures flow.

- Features include - Wt. 0.4kg
- ▶ Precision engineered
 - ▶ Pressure compensated design accurately measures flow
 - ▶ Easy to read
 - ▶ Scratch resistant
 - ▶ Virtually unbreakable transparent cover tube
 - ▶ Colour coded handwheel

86615LOX Oxygen flowmeter 0-15LPM

86615LAIR Air flowmeter 0-15LPM

86630LOX Oxygen flowmeter 0-30LPM

86630LAIR Air flowmeter 0-30LPM



801Z259



801Z260

