

empire® mercury sphygmomanometer


The mercury sphygmomanometer featuring the **latest Riester measuring technology**. Safe and accurate.

- **Maximum error tolerance of +/- 3 mmHg**
- **99.99% pure mercury**
- **Precision glass column**
Internal diameter 4.2 mm +/- 0.2 mm
 - Easily legible scale up to 300 mmHg
 - Mercury lock
 - Special seal against mercury contamination
- **empire® cotton bandage cuff adults**
- **empire® cotton hook cuff adults**
- **empire® calibrated nylon velcro adults**
- **Universal clamp incl. screws**
- **Wall support incl. screws (for conversion to wall-mounted model)**
- **Latex bulb with chromium plated release valve**
- **Precision air release valve, wear-free and with fine regulation**
- **Microfilter for long life**
- **Cuff with two-tube latex bladder**
 - Coil tubing extendible up to 3 m
 - Cleaning device for glass column
 - Spacious basket for tubes and cuffs
 - ABS casing

Art. No.: 1130
Art. No.: 1131
Art. No.: 1132
Art. No.: 10384
Art. No.: 54663

empire® desk design

SPHYGMOMANOMETERS



empire® mercury sphygmomanometer

The mercury sphygmomanometer featuring the **latest Riester measuring technology**. Safe and accurate.

- **Maximum error tolerance of +/- 3 mmHg**
- **99.99% pure mercury**
- **Precision glass column**
Internal diameter 4.2 mm +/- 0.2 mm

- Easily legible tilted scale up to 300 mmHg
- Coil tubing extendible up to 3 m
- Mercury lock
- Special seal against mercury contamination

- empire® cotton bandage cuff adults
- empire® cotton hook cuff adults
- empire® calibrated nylon velcro adults

- **Latex bulb with chromium plated release valve**
- **Precision air release valve, wear-free and with fine regulation**
- **Microfilter for long life**
- **Cuff with two-tube latex bladder**

- Resistant plastic housing
- Cleaning device for glass column
- ABS casing

Art. No.: 1090

Art. No.: 1091

Art. No.: 1092



Rudolf Riester GmbH | P.O. Box 35 | Bruckstr. 31 | D-72417 Jungingen | Germany
Tel.: (+49) +7477-9270-0 | Fax: (+49) +7477-9270-70 | info@riester.de | www.riester.de




empire® mercury sphygmomanometer

The mercury sphygmomanometer featuring the **latest Riester measuring technology**. Safe and accurate.

- **Maximum error tolerance of +/- 3 mm Hg**
- **99.99% pure mercury**
- **Precision glass column**
Internal diameter 4.2 mm +/- 0.2 mm
- Easily legible tilted scale up to 300 mmHg
- Mercury lock
- Special seal against mercury contamination
- Resistant plastic housing
- Coil tubing extendible up to 3 m
- **empire® cotton bandage cuff adults**
- **empire® cotton hook cuff adults**
- **empire® calibrated nylon velcro adults**
- Basket small: L250 x B90 x H200 mm
- Basket large: L230 x B230 x H250 mm
- **Latex bulb with chromium plated release valve**
- **Precision air release valve, wear-free and with fine regulation**
- **Microfilter for long life**
- **Cuff with two-tube latex bladder**
- Height adjustable from 0.75 to 1.20 m
- Diameter 53 cm
- Cleaning device for glass column
- Easy to move due to low friction castors

Art. No.: 1110
Art. No.: 1111
Art. No.: 1112
Art. No.: 10392
Art. No.: 10393



empire® mercury sphygmomanometer

The mercury sphygmomanometer featuring the **latest Riester measuring technology**. Safe and accurate.

- **Maximum error tolerance of +/- 3 mmHg**
- **99.99% pure mercury**
- **Precision glass column**
Internal diameter 4.2 mm +/- 0.2 mm
- Easily legible scale up to 300 mmHg
- Mercury lock
- Special seal against mercury contamination
- Coil tubing extendible up to 3 m
- Inclusive of screws and dowels for all mounting
- **empire® cotton bandage cuff adults**
- **empire® cotton hook cuff adults**
- **empire® calibrated nylon velcro adults**
- Universal clamp incl. screws
- Rail support (for conversion to the anesthesia model)
- **Latex bulb with chromium plated release valve**
- **Precision air release valve, wear-free and with fine regulation**
- **Microfilter for long life**
- **Cuff with two-tube latex bladder**
- Cleaning device for glass column
- Rotating up to 45°
- Spacious basket for tubes and cuffs
- ABS casing

Art. No.: 1100
Art. No.: 1101
Art. No.: 1102
Art. No.: 10384
Art. No.: 54665



 1800 67 5432

 1800 67 5433

 sales@ecomед.com.au

 www.ecomed.com.au