

graeve, cusco, speculight

 **Riester**

The familiar way

ECOMED

Serving The HealthCare Industry



CUSCO



speculight



graeve

The gynaecologist can select from two vaginal specula models and **different sizes**. Regardless which quality instrument is chosen, all specula excel **through their particular reliability and durability**. The **speculight provides for the right illumination** – and is the optimal addition to all vaginal specula.

graeve vaginal speculum, external polished finish
cusco vaginal speculum, external matt finish

- No reflections thanks to matt inner surface
- Made of stainless steel type 18/8
- Can be sterilized up to 134°C in an autoclave
- Each speculum is supplied in a cardboard box

speculight

- Battery handle C, chrome-plated
- For 2 batteries type C
- 85 cm cable with clip and 2.5 V standard illumination (vacuum)
- Light intensity can be regulated at will using the rheostat
- Supplied boxed

graeve Fig. 1, 75 x 20 mm
graeve Fig. 2, 95 x 35 mm
graeve Fig. 3, 115 x 35 mm

No. 5201
No. 5202
No. 5203

cusco virgins, 75 x 17 mm
cusco Fig. 1, 75 x 32 mm
cusco Fig. 2, 85 x 35 mm
cusco Fig. 3, 100 x 35 mm


No. 5230
No. 5231
No. 5232
No. 5233

speculight, lighting attachment
for vaginal specula

No. 5240

Distributed by



 1800 67 5432

 1800 67 5433

 sales@ecomед.com.au

 www.ecomed.com.au