

# metpak

PRESSURE INFUSION INSTRUMENTS



metpak

Pressure infusion devices are an essential item in any hospital! All the more so when it comes to the high-quality products from **Riester**, which excel by their particularly high degree of **reliability and durability**. Pressure cuffs in combination with a manometer for producing pressure provide **for optimal pressure infusions**. For solutions or blood bags from 500 ml, 1000 ml, 3000 ml or 5000 ml.

- **Sturdy cotton cuff in blue, with suspension loop for infusion stands, washable up to 60°C**
- **Precise reading of the level in the bag thanks to a transparent net on the cuff**
- **Chrome-plated manometer with easy-to-read 49 mm diameter aluminium scale up to 300 mm Hg**
- **Specially hardened copper-beryllium membrane**
- **No zero point setting required**


- **Latex-free pad, resistant to pressure**
- **Latex-free inflation bulb**

metpak, 500 ml  
metpak, 1000 ml  
metpak, 3000 ml  
metpak, 5000 ml

Art. No.: 5270  
Art. No.: 5275  
Art. No.: 5270-536  
Art. No.: 5270-537

Distributed by



 1800 67 5432

 1800 67 5433

 [sales@ecomед.com.au](mailto:sales@ecomед.com.au)

 [www.ecomed.com.au](http://www.ecomed.com.au)