



## CAPNOGRAPHY SUPPLIES & ACCESSORIES

### Capno Pod (92516), Capnography Module (92517) & Multigas Module (92518)

The Nomoline multi-patient use adapter for sidestream monitoring offers a cost effective solution for both intubated and non-intubated patients. The sophisticated technology removes water from the sample line and reduces potential for water intrusion and cross contamination. Utilizing the Nomoline multi-patient use adapter in conjunction with single use CO<sub>2</sub> cannulas, we are able to offer a capnography solution at a much lower cost.

Mainstream monitoring is also available utilizing disposable airway adapters with a reusable sensor that is extremely lightweight at <25g!\*

#### SIDESTREAM MONITORING SUPPLIES

- Nomoline multi-patient use adapter

##### Non-intubated Patients

- CO<sub>2</sub> Nasal Cannulas
  - Adult, Infant and Pediatric
- CO<sub>2</sub>/O<sub>2</sub> Divided Nasal Cannulas
  - Adult, Infant and Pediatric
- CO<sub>2</sub>/O<sub>2</sub> Oral/Nasal Cannulas
  - Adult and Pediatric

##### Intubated Patients

- Sidestream Sampling Lines
  - 6.5 ft and 10 ft
- Sidestream T-Adapter

#### MAINSTREAM MONITORING SUPPLIES

##### Intubated Patients Only

- Mainstream Sensor
- Mainstream Airway Adapters
  - Adult / Pediatric and Infant

#### HIGHLIGHTS

- Utilize same Nomoline supplies across all three capnography monitoring devices
- Nomoline multi-patient use adapter with CO<sub>2</sub> Cannulas helps to reduce costs while still protecting your investment
- Mainstream Sensor is extremely small and lightweight\*

\* Available with the 92517 module only.

## CAPNOGRAPHY SUPPLIES & ACCESSORIES

PART NUMBER	DESCRIPTION	TYPE	QTY
<b>SIDESTREAM MONITORING ADAPTER</b>			
103-0234-00	Sidestream Nomoline Adapter with Female Luer Lock	Multi-Patient Use	25/Box
<b>SIDESTREAM MONITORING SUPPLIES - NON-INTUBATED PATIENTS</b>			
704-0176-00	CO <sub>2</sub> Nasal Cannula, Adult, 7 ft.	Single Use	25/Box
704-0177-00	CO <sub>2</sub> Nasal Cannula, Pediatric, 7 ft.	Single Use	25/Box
704-0178-00	CO <sub>2</sub> Nasal Cannula, Infant, 7 ft.	Single Use	25/Box
704-0181-00	CO <sub>2</sub> /O <sub>2</sub> Divided Nasal Cannula, Adult, 7 ft.	Single Use	25/Box
704-0180-00	CO <sub>2</sub> /O <sub>2</sub> Divided Nasal Cannula, Pediatric, 7 ft.	Single Use	25/Box
704-0179-00	CO <sub>2</sub> /O <sub>2</sub> Divided Nasal Cannula, Infant, 7 ft.	Single Use	25/Box
704-0183-00	CO <sub>2</sub> /O <sub>2</sub> Oral/Nasal Cannula, Adult, 7 ft.	Single Use	25/Box
704-0182-00	CO <sub>2</sub> /O <sub>2</sub> Oral/Nasal Cannula, Pediatric, 7 ft.	Single Use	25/Box
<b>SIDESTREAM MONITORING SUPPLIES - INTUBATED PATIENTS</b>			
166-7085-00	Sidestream Sampling Line with Male Luer Lock on Both Ends, 6.5 ft.	Single Use	25/Box
166-7121-00	Sidestream Sampling Line with Male Luer Lock on Both Ends, 10 ft.	Single Use	50/Box
704-0175-00	Sidestream T-Adapter Requires a Sampling Line for Use	Single Use	25/Box
<b>MAINSTREAM MONITORING SUPPLIES - INTUBATED PATIENTS ONLY</b>			
010-1980-00	Mainstream Sensor	Reusable	1 Each
704-0173-00	Mainstream Airway Adapter, Adult/Pediatric	Single Use	25/Box
704-0174-00	Mainstream Airway Adapter, Infant	Single Use	25/Box



Multi-Patient use  
Nomoline Adapter



Nomoline Adapter with  
Sampling Line and T-Adapter



Mainstream CO<sub>2</sub> Sensor



Airway Adapters (Infant, Adult/Pediatric)

Contact your Spacelabs Healthcare representative to discuss our broad supplies and accessories portfolio.

Spacelabs Healthcare reserves the right to make changes to the products at any time without notice. All brands and product names are trademarks of their respective owners. Product is CE Marked.

Product may not be available for sale in all locations. Check with your local Spacelabs representative.

35301 SE Center Street, Snoqualmie, WA 98065 | T +1 425 396 3300 | +1 800 522 7025

© 2013 Spacelabs Healthcare Specifications subject to change without notice. S4111-Ju14

Distributed by



1800 67 5432

1800 67 5433

sales@ecomед.com.au

www.ecomed.com.au



**SPACE LABS**  
HEALTHCARE

www.spacelabshealthcare.com