

Its intuitive touchscreen interface makes it easy to access patient information. Its LCD widescreen display adds clarity and style. Its elegant design optimizes space without sacrificing features. Most important, élance offers robust and reliable patient-side monitoring wherever it is needed.



Distributed by



1800 67 5432  
 1800 67 5433  
 sales@ecomед.com.au  
 www.ecomed.com.au



35301 SE Center Street Snoqualmie, Washington 98065 | T +1 425 396 3300 | +1 800 522 7025 • 1 Harford Court, John Tate Road, Hertford, SG13 7NW, UK | T +44 (0) 1992 507700 | F +44 (0) 1992 501213  
 www.space-labs-healthcare.com © 2012 Space Labs Healthcare Specifications subject to change without notice. M1600-Ju13






# Small Monitor. Big Heart.

We know a hospital stay is stressful for both patients and their families. We understand how important it is to provide the best possible experience along with the highest level of care. That's why we developed the élançe monitor.

Slim and lightweight, this patient-side monitor has been meticulously engineered to meet the needs of caregivers and their patients. Its intuitive design provides the highest level of performance and reliability, supporting the best of care.

	élançe 5	élançe 5i	élançe 5c	élançe 7	élançe 7i	élançe 7c
Patient Type						
Pediatric	▼	▼	▼	▼	▼	▼
Adult	▼	▼	▼	▼	▼	▼
Display						
10.2" Widescreen LCD	▼	▼	▼			
12.1" Widescreen LCD				▼	▼	▼
Touchscreen Controls	▼	▼	▼	▼	▼	▼
800 x 480 resolution	▼	▼	▼			
1280 x 800 resolution				▼	▼	▼
Waveforms	up to 4	up to 5	up to 5	up to 6	up to 6	up to 6
Parameters						
ECG - 3 & 5 Lead	▼	▼	▼	▼	▼	▼
Resp	▼	▼	▼	▼	▼	▼
NIBP	▼	▼	▼	▼	▼	▼
SpO <sub>2</sub>						
Spacelabs Std	Opt	Opt	Opt	Opt	Opt	Opt
Masimo	Opt	Opt	Opt	Opt	Opt	Opt
Nellcor	Opt	Opt	Opt	Opt	Opt	Opt
Temp - 2 ch	▼	▼	▼	▼	▼	▼
IBP - 2 ch		▼	▼		▼	▼
EtCO <sub>2</sub> - Microstream			▼			▼
Features						
ST Segment Analysis	Opt	Opt	Opt	Opt	Opt	Opt
Enhanced Arrhythmia Analysis	Opt	Opt	Opt	Opt	Opt	Opt
120 Hrs. Trends	▼	▼	▼	▼	▼	▼
Ext. 2 Ch Recorder	Opt	Opt	Opt	Opt	Opt	Opt
Networking	Opt	Opt	Opt	Opt	Opt	Opt
Power						
Internal Power Supply	▼	▼	▼	▼	▼	▼
Li-ion, Internal Battery	▼	▼	▼	▼	▼	▼
5-Hour Battery Time	▼	▼	▼	▼	▼	▼



## REFINED RELIABILITY

élançe provides high-performance patient monitoring with essential parameters throughout the care continuum. Advanced design features, such as its fanless technology, enhance reliability.

More than a bedside monitor, élançe is a patient-side monitor. Lightweight, its built-in handle and simple, quick-release mounting options make it highly mobile.



## ROBUST PERFORMANCE

élançe offers a robust feature set to provide the best possible patient care. All essential parameters such as ECG, Respiration, NIBP, SpO<sub>2</sub>, and Temp are monitored with precision, and are clearly displayed on the ultra-bright LCD widescreen.


Whether you choose a 10.2" or 12.1" wide flat screen, 5 or 6 waveforms, the sophisticated system architecture and algorithms ensure the highest levels of performance and reliability.



## EXCEPTIONAL USER EXPERIENCE

The élançe monitor has an ultra-bright and vivid widescreen display, a touch screen user interface and a simple, single-level menu structure. All of this makes it easier and more intuitive to use. Its contemporary slim design is available in either natural white or elite black.

Microsoft® Windows® CE technology makes it easy to network élançe monitors directly to a Windows based central station, bringing together information for up to 32 patients seamlessly.



## PERFORMANCE WITH STYLE