

UltraCheck® & SoftCheck®

Your requirements for cuffs vary by care area and patient need. Statcorp Medical cuffs may be ordered with connectors that make them compatible with most brands of automated blood pressure monitors and wall or portable manometers.

Durable, leak-proof, comfortable, accurate, and economical

Three styles are available:

UltraCheck Reusable

- PVC / DEHP free
- Anti-microbial protection built into the fabric
- Fluid barrier for isolation areas

Features

- Meets AHA/ AAMI Standards
- Color coded
- 8 sizes for accurate measurements
- Genuine Velcro® closures
- Latex free
- Fabric fold at elbow is more comfortable than a weld
- Statcorp adapters allow for compatibility with monitors worldwide





Fittings



DM

HP



ML

FL



MQ

MQML

MQMF (MQ + MQML)Manuf

Manual bulbs



Statcorp Medical

[illegible]

*--To order cuffs for a specific monitor, add the appropriate suffix to the product number.
For Critikon Dinamap®, add DM, for Datascope®, Spacelabs®, and Philips®, add HP, for GE Healthcare®, add MQ, and for Datex-Ohmeda®, add MQMF.

Statcorp, Inc.
an OSI Optoelectronics company
14476 Duval Place West, Suite 303
Jacksonville, FL 32218
800-992-0014, 904-786-5113
www.statcorpmedical.com

Product is  Marked

 ISO 13485
 Certificate Number FM 40458

ECOMED
Serving The HealthCare Industry

 1800 67 5432
 1800 67 5433
 enquiries@ecomed.com.au
 www.ecomed.com.au

Statcorp Medical reserves the right to make changes to this sheet and the product at any time without notice.
Velcro® is a registered trademark of VELCRO USA, Inc. All other product names mentioned herein may be trademarks and/or registered trademarks of their respective owners.

BPCSHT Rev 02