



ACCESSORIES + Advance

Canisters

- 3833-132 - 1 Litre 'V' Mount Advance Canister
- 3833-134 - 1 Litre Parallel Mount Advance Canister
- 3833-090 - 2 Litre 'V' Mount Advance Canister
- 3833-098 - 2 Litre Parallel Mount Advance Canister
- 3833-041 - 3 Litre 'V' Mount Advance Canister
- 3833-042 - 3 Litre Parallel Mount Advance Canister
- 3833-055 - 2 Litre Advance Cascade 'V' Mount Canister
- 3833-058 - 2 Litre Advance Cascade Parallel Mount Canister

Vacuum Tubing

- 4160-007 - 1m Yellow Vacuum Tubing Set
- 4160-002 - 2m Yellow Vacuum Tubing Set

Tissue Collection/Sample

- 9410-001 - Tissue Collection Device Pack of 12
- 9970-006 - Critical Measuring Vessel Pack of 12
- 9970-007 - Sample Sock Pack of 50

Tubing Adaptors (for different size patient tubing)

- 9910-040 - Tubing Adaptor pack of 30

Disposal

- 2940-035 - Disposal Bin (Rigid Plastic) up to 7 Liners Pack of 10
- 2940-036 - Disposal Box for up to 6 x 1L or 2 L Liners Pack of 10
- 9940-340 - Disposal Box for up to 2 x 1L or 2 L Liners Pack of 10
- 9310-341 - Disposal Box for up to 3 x 3L Liners Pack of 10

Pipeline Protection

- 9510-100 - Pipeline Protector (Filter) Bacterial Pack of 10
- 3786-024 - Viral Filter Pack of 10



All liners manufactured by VacSax are made from High Density PE and do not contain any PVC and produce only carbon dioxide and water upon 100% combustion helping to save the environment.

Every effort has been made to ensure the accuracy of information contained within this literature.

However freedom from errors and omissions cannot be guaranteed.

VacSax policy is one of continued product improvement.

The company reserves the right to change the information herein without prior notice.



Certificate No. 542
ISO 9001:2008 & ISO 13485:2003



All VacSax products are CE marked



Class 1 Measuring Applies only to CMV

ECOMED

Serving The HealthCare Industry

t: 1800 675 432

f: 1800 675 433

e: enquiries@ecomед.com.au

w: www.ecomed.com.au



RANGE + Advance

Disposable Suction Liners

The Advance Range from VacSax is a closed disposable suction system, the ideal choice for all areas within the Hospital. Advance is available in 1 Litre, 2 Litre and 3 Litre versions.

Careful consideration was given to the design ensuring users have the most easy to use and safe product. Designed with just two Ports, patient and vacuum, both of different dimensions ensuring cross connection is eliminated.

Each Liner with the exception of the cascade version has a hydrophobic bacterial shut off valve located in the cap for complete overflow protection.

The unique patented design of the bag is folded for ease of storage and begins to inflate immediately the vacuum is attached. After use the integral carry handles on the liner allows for the safe removal from the canister and the bag will stand upright whilst full of fluid.

Within the Advance range is the GelSax Advance PreGelled Liner. This liner has exactly the right amount of PreGel which will quickly solidify the contents of the liner.

PreGel has been formulated not to expand when wet, giving the user confidence in the approximate volume of the Liner.



www.vacsax.com



Advance

VacSax Advance

Part Number

- 9910-340 - 1 Litre Advance Liner Box of 25
 - VAL - 201 - 2 Litre Advance Liner Box Of 30
 - VAL - 301 - 3 Litre Advance Liner Box of 20
-
- VAL - 203 - 2 Litre Advance Cascade Liner Box of 30 including 30 cascade tubing sets
 - VAL - 303 - 3 Litre Advance Cascade Liner Box of 20 including 20 cascade tubing sets
-
- 9910-350 - 4 Canister Advance Cascade Stand including Suction (others available on request)

GELSAX Advance

GelSax Advance Pre-Gelled Liner

Liners with just the right amount of PreGel

- 9810-340 - 1 Litre GelSax Advance Liner Box of 25
 - GAL - 211 - 2 Litre GelSax Advance Liner Box of 30
 - GAL - 311 - 3 Litre GelSax Advance Liner Box of 20
-
- GAL - 213 - 2 Litre GelSax Advance Cascade Liner Box of 30 including 30 cascade tubing sets
 - GAL - 313 - 3 Litre GelSax Advance Cascade Liner Box of 20 including 20 cascade tubing sets
-
- 9910-350 - 4 Canister Advance Cascade Stand including Suction (others available on request)

