



the difference is clear

Precision Syringes

Overview

Vitalograph precision syringes provide quick and easy accuracy checks with reliability and fast verification.

Our precision syringes are available in 1-L and 3-L volumes and can be used with other manufacturers' equipment. Adaptors may be required.

3-L Precision Syringe

Spirometry guidelines jointly published by the American Thoracic Society and the European Respiratory Society advocate the use of a 3-L precision syringe for daily accuracy checks.

Vitalograph has designed a NEW, extremely accurate, robust and lightweight 3-Litre precision syringe to give you the ultimate confidence in your spirometry measurements.

Ordering Information: 36020 3-L Precision Syringe.



1-L Precision Syringe

Our convenient 1-L version with built in valves allows for fast accuracy checks and for linearity checks on volume displacement spirometers.

Ordering Information: 20408 1-L Precision Syringe

