

the difference is clear

Vitalograph myPEF Peak Flow Meter

Mechanical Peak Flow Meter with secure online data capture

Web enabled myPEF™ is supplied with a unique ID card allowing users access to a secure website where they can enter their PEF readings.

- ❖ Eliminates the need for paper records, supporting data accuracy
- ❖ Easy and convenient to use, encouraging greater compliance
- ❖ Data trending and customised colour zones can provide pre-symptomatic indications of potential asthma exacerbations
- ❖ Supports more effective patient coaching and asthma management
- ❖ Easy to read graphical data display and report generation featuring plotting, trending and time stamping functions
- ❖ Clinicians can be granted access to view patients' PEF results online for reliable, cost effective remote monitoring
- ❖ Suitable for use with any web enabled device. Non web enabled version also available
- ❖ Attractively priced
- ❖ Ergonomically designed for maximum comfort and ease of use

Vitalograph is a world leading provider of outstanding quality cardio-respiratory diagnostic devices, clinical trial services and medical equipment servicing. With a pioneering heritage of excellence spanning half a century Vitalograph continues to make valuable contributions to effective medical care and enhanced quality of life.





Technical Specification:

Product: myPEF™ Peak Flow Meter

Model Number: 4300

Test Types: Single breath

Parameters Measured: PEF

Device Dimensions:

138mm x 56mm x 23mm (unpackaged)

142mm x 81mm x 26mm (individually packaged)

Nett Weight:

79g (unpackaged)

96g (individually packaged)

Gross Weight and Size (box of 25):

Weight: 2.54 kg

Packed Size: 320mm x 170mm x 150mm

Environmental Data:

Storage Humidity: 10% to 95%

Storage Temperature: 0°C to 50°C

Recommended Operating Temperature Range:

17 - 37°C; Design Limits 10-40°C

Ambient Humidity: 0% to 99%

Performance Standards:

EN/ISO 23747:2009; ATS/ERS:2005; AS/NZS4237

Medical Safety Standard: Medical Devices Directive 93/42/EEC 1993 (as amended)

Designed & Manufactured Under:

ISO 13485:2003, FDA 21CFR820, CMDR & JPAL

Ordering Info:

- 43950 4300 myPEF web enabled Peak Flow Meter
- 43900 4300 myPEF non web enabled Peak Flow Meter
- 29805 myPEF Mouthpiece Adaptors (96)
- 20242 2024 SafeTway Mouthpieces (200)
- 20980 2024 Mini SafeTway Mouthpieces (50)
- 20303 2030 Noseclip Disposable (200)

Distributed by



☎ 1800 67 5432

📄 1800 67 5433

✉ sales@ecomед.com.au

🌐 www.ecomed.com.au

Vitalograph®

Your respiratory partner

UK Sales

Vitalograph Ltd.
Maids Moreton Buckingham
MK18 1SW England
Tel: 01280 827110
Fax: 01280 823302
e-mail: sales@vitalograph.co.uk
www.vitalograph.co.uk

Export Sales

Vitalograph Ltd.
Maids Moreton Buckingham
MK18 1SW England
Tel: +44 1280 827120
Fax: +44 1280 823302
e-mail: sales@vitalograph.co.uk
www.vitalograph.eu

Vitalograph GmbH

Rellinger Straße 64a
20257 Hamburg
Germany
Tel: (040) 547391-0
Fax: (040) 547391-40
e-mail: info@vitalograph.de
www.vitalograph.de

Vitalograph Inc.

13310 West 99th Street
Lenexa
Kansas 66215 USA
Tel: (913) 888 4221
Fax: (913) 888 4259
e-mail: vitcs@vitalograph.com
www.vitalograph.com

Vitalograph (Ireland) Ltd.

Gort Road Business Park
Ennis
Co Clare Ireland
Tel: (065) 6864100
Fax: (065) 6829289
e-mail: sales@vitalograph.ie
www.vitalograph.ie

www.vitalograph.co.uk

Vitalograph® and myPEF™ are trademarks or registered trademarks of Vitalograph Ltd.



PRINT REF: 19264_1